

---

**Probe for EPMA v. 9.4.3**

**Advanced Topics**  
*Extreme Edition*



***By Daniel T. Kremser, Ph.D.  
for Probe Software, Inc.  
© Copyright 1994-2013***

# Contents

<b>ACKNOWLEDGEMENTS .....</b>	<b>1</b>
<b>DISCLAIMER.....</b>	<b>1</b>
<b>INTRODUCTION.....</b>	<b>1</b>
<b>ELEMENT SETUPS .....</b>	<b>3</b>
<b>SAMPLE SETUPS.....</b>	<b>15</b>
<b>FILE SETUPS .....</b>	<b>31</b>
<b>WAVESCAN.....</b>	<b>41</b>
QUICK WAVESCAN ACQUISITION.....	41
CALIBRATED MULTI-ELEMENT WAVESCAN .....	45
<b>POLYGON GRIDGING – USING SURFER OPTION .....</b>	<b>57</b>
<b>STAGE BIT MAPS AND PICTURE SNAP! FEATURE .....</b>	<b>82</b>
<b>MODAL ANALYSIS .....</b>	<b>100</b>
<b>DEADTIME CALCULATIONS .....</b>	<b>118</b>
<b>CALCULATION OPTIONS .....</b>	<b>126</b>
<b>LINEAR CALIBRATION - CURVE METHOD.....</b>	<b>142</b>
<b>TIME DEPENDENT INTENSITY (TDI) CORRECTIONS.....</b>	<b>154</b>
<b>ADVANCED INTERFERENCE CORRECTIONS .....</b>	<b>172</b>
<b>LIGHT ELEMENT ANALYSIS - EMPIRICAL APFS.....</b>	<b>184</b>
<b>REFERENCES.....</b>	<b>192</b>

(This page intentionally left blank)

# Acknowledgements

Many people have aided the author both directly and indirectly in the preparation of this Advanced Topics manual. First and foremost, John Donovan and Paul Carpenter, they have freely given their time, expertise, and even data files for incorporation into this manual. Their technical discussions as well as their friendship over the years have made this manual better.

## Disclaimer

IN NO EVENT SHALL DAN KREMSEK, ADVANCED MICROBEAM OR PROBE SOFTWARE BE LIABLE TO ANY PARTY FOR DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, INCLUDING LOST PROFITS, ARISING OUT OF THE USE OF THIS DOCUMENTATION, EVEN IF DAN KREMSEK, ADVANCED MICROBEAM OR PROBE SOFTWARE HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

DAN KREMSEK, ADVANCED MICROBEAM, AND PROBE SOFTWARE SPECIFICALLY DISCLAIMS ANY WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THIS DOCUMENTATION PROVIDED HEREUNDER IS ON AN "AS IS" BASIS, AND DAN KREMSEK, ADVANCED MICROBEAM, AND PROBE SOFTWARE HAVE NO OBLIGATIONS TO PROVIDE MAINTENANCE, SUPPORT, UPDATES, ENHANCEMENTS, OR MODIFICATIONS.

## Introduction

PROBE FOR EPMA is a very versatile and powerful 32-bit acquisition, automation and analysis package for WDS and EDS electron microprobe analysis running under Microsoft Windows operating system.

One of the strengths of PROBE FOR EPMA is the wide variety of options and features for many different tasks that are available to the probe operator. The aim of this manual then is to document some of the more advanced features usually skipped over in an introductory text. And as always, the path taken to cover a feature may not be the only avenue to approach the subject.

This manual was originally produced on the Washington University (Earth and Planetary Sciences) JEOL 733 Superprobe equipped with three wavelength dispersive spectrometers and modified using PROBE FOR EPMA in demo mode.

(This page intentionally left blank)

# Element Setups

Individual element analytical configurations for a specific element, x-ray, spectrometer, and reflecting crystal may be saved to the SETUP.MDB database for use in creating new sample setups within a probe run, for use in future runs or for documentation and performance evaluation purposes. The example below will illustrate how to create element setups from within a typical eight-element olivine routine and store them in a new SETUP.MDB database.

Open a new PROBE FOR EPMA run in the usual manner. From the **Acquire!** window, create a new unknown sample from the **New Sample** dialog box, then click the **Elements/Cations** button. Next, enter the elements of interest into the **Analyzed and Specified Elements** window in the usual manner. Below is the completed **Analyzed and Specified Elements** window after the entry of all seven elements plus oxygen.

**Analyzed and Specified Elements**

**Selected Samples**

Un 1 \* setup for olivines

OK Cancel

Save Element Setup

Save Sample Setup

Add/Remove Standards

Reload Standard Assignments

Remove TDI Correction


1 2 3 4 5 6

Click Element Row to Edit Element/Cations Parameters (click empty row to add)

Channel	Element	X-Ray	Analyzed	Cations	Oxygens	Off-Peak
1	si	ka	Yes	1	2	Linear
2	ti	ka	Yes	1	2	Linear
3	al	ka	Yes	2	3	Linear
4	fe	ka	Yes	1	1	Linear
5	mn	ka	Yes	1	1	Linear
6	mg	ka	Yes	1	1	Linear
7	ca	ka	Yes	1	1	Linear
8	o		No	1	0	

Go through the calibration process; find new peak positions and standardize to acquire intensity data on each standard. Normally one should save the element setup of an element that is assigned as the standard for that element. This is done because in that case the x-ray intensity data, P/B data, PHA parameters, and other information will also be saved in the SETUP.MDB database. This information is very useful for documentation and evaluation purposes.

After completing the calibration and standardization process, open the **Analyze!** window. Choose the element setup to be stored and highlight the standard (iron in fayalite, in this example).



Click the **Elements/Cations** button.

This opens the **Analyzed and Specified Elements** window.

**Analyzed and Specified Elements**

**Selected Samples**

St 203 Set 1 Fayalite

**OK** **Cancel**

**Save Element Setup**

**Save Sample Setup**

**Add/Remove Standards**

Reload Standard Assignments

Remove TDI Correction

1 2 3 4 5 6

**Click Element Row to Edit Element/Cations Parameters (click empty row to add)**

Channel	Element	X-Ray	Analyzed	Cations	Oxygens	Off-Peak
1	fe	ka	Yes	1	1	Linear
2	si		No	1	2	
3	ti		No	1	2	
4	al		No	2	3	
5	mn		No	1	1	
6	mg		No	1	1	
7	ca		No	1	1	
8	o		No	1	0	

Click the **Save Element Setup** button.



The **Element Setup Database** opens.

Element Setup Database

Current Sample:  
St 203 Set 1 Fayalite

fe ka Spectro 3 LIF (134.626)  
si (specified)  
ti (specified)  
al (specified)  
mn (specified)  
mg (specified)  
ca (specified)  
o (specified)

Double click Analyzed Element List to see Element Setups

<< Add to Sample

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

Close

Element Setup Data From SETUP.MDB Database

◀◀ Data Cursor ▶▶

Delete from Database

Enter Search Element >>

Total Records = 0

Elem/Order/CatOxd

Spec/Crystal/2d

User Name

Sample Name

Date - Time

Probe Data File

On/Hi/Lo Pos

BgdType/Offset

S-Hi/Lo/Exponen

Base/Win/Gain/Bias

KeV/TO/DT/DIFF

Standard Intensity Data

Std/PeakToBgd

On/Hi/Lo sec

On/Hi/Lo cps

Beam/Abs current

Stage Position

SETUP.MDB SETUP2.MDB (MAN) SETUP3.MDB (Interf.)

Import Export

Wavescan and Peaking Parameters

Wavescan Hi/Lo/Points/Time

Peakscan Hi/Lo/Points/Time

Start/Stop/PB/Count/Attempts

Highlight the specific element to save and click the **Add To Database>>** button. In this standard, the iron intensity is the only element to be archived, select *fe ka Spectro 3 LIF (134.626)*.

Record number 1 has been stored as illustrated below. Note that the *Standard Intensity Data* and *Wavescan and Peaking Parameters* are stored as well.

**Element Setup Database**

**Current Sample:**  
St 203 Set 1 Fayalite

fe ka Spectro 3 LIF (134.626)  
 si (specified)  
 ti (specified)  
 al (specified)  
 mn (specified)  
 mg (specified)  
 ca (specified)  
 o (specified)

Double click Analyzed Element List to see Element Setups

<< Add to Sample

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

Close

**Element Setup Data From SETUP.MDB Database**

◀◀ Data Cursor ▶▶

Delete from Database

Enter Search Element >>

Total Records = 1

**Elem/Order/CatOxd** fe ka 1 1 1

**Spec/Crystal/2d** 3 LIF 4.0267 0.000058

User Name Dan Kremser

Sample Name Fayalite

Date - Time 10/15/2005 11:35:22 AM

Probe Data File C:\Documents and Settings\Dan\Desktop\July

On/Hi/Lo Pos 134.626 137.471 131.78

BgdType/Offset OFF LIN 2.845 -2.846

S-Hi/Lo/Exponen 1 1 1

Base/Win/Gain/Bias 1 10 200 1700

KeV/TO/DT/DIFF 15 40 2.66E-06 0

**Standard Intensity Data**

Std/PeakToBgd 203 213.5893

On/Hi/Lo sec 10 5 5

On/Hi/Lo cps 2392.2 11.6 10.8

Beam/Abs current 31.061 0

Stage Position 6.394809 27.67546 11.0132 1

SETUP.MDB  SETUP2.MDB (MAN)  SETUP3.MDB (Interf.)

Import Export

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time 138.9923 130.4557 100 4

Peakscan Hi/Lo/Points/Time 137.9252 131.5228 40 4

Start/Stop/PB/Count/Attempts 8.173325E-02 5.448884E-03 3 10 30


Click the **Close** button. The **Analyzed and Specified Elements** window reappears. Click the **OK** button. The **Analyze!** window returns.

The other element setups from this calibrated and standardized run or other probe runs may be entered into the database in a similar manner for future use.

The **Save Setups** button in the **Analyze!** window will allow the user to save all element setups in the selected sample (highlighted in the *Sample List*) to the element database (SETUP.MDB).

To recall an element setup from the SETUP.MDB database for a new sample setup follow the procedure outlined below. Open a new PROBE FOR EPMA run. This process will also be applicable if the user simply wants to add an element to an existing sample setup. This example will illustrate recalling elements from the database for the setup of a new pyroxene run.

From the **Acquire!** window, click the **New Sample** button.



This opens the **New Sample** window. Edit the *New Sample Name* text field.

**New Sample**

**New Sample Type**

Standard

**Unknown**

Wavescan

**Add/Remove Standards**

**OK** **Cancel**

**Load Element Setups**

Load Sample Setup

**Load File Setup**

Load Multiple Setup

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

Load Wavescan From Another Probe Run

**New Sample Name**

setup for pyroxene analyses

**New Sample Description** **Add <cr>**

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.

Click the **Load Element Setups** button in the **New Sample** window.

This opens the **Element Setup Database**.

**Element Setup Database**

**Current Sample:**  
Un 1 setup for pyroxene analyses

Double click Analyzed Element List to see Element Setups

<< Add to Sample

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

Close

**Element Setup Data From SETUP.MDB Database**

◀◀ Data Cursor ▶▶

Enter Search Element >>

Delete from Database

Total Records = 21

Elem/Order/CatOxd	fe	ka		1	1
Spec/Crystal/2d	3	LIF	4.0267	0.000058	
User Name	Dan Kremser				
Sample Name	Fayalite				
Date - Time	3/16/1999 2:55:43 PM				
Probe Data File	E:\Probe Operators\Manual Files\olivine.mdb				
On/Hi/Lo Pos	134.5533	137.3983	131.7073		
BgdType/Offset	OFF	LIN	2.845	-2.846	
S-Hi/Lo/Exponen	0	0	1		
Base/Win/Gain/Bias	1	10	200	1700	
KeV/TO/DT/DIFF	15	40	0.000001	0	


**Standard Intensity Data**

Std/PeakToBgd	203	172.0335			
On/Hi/Lo sec	20	4	4		
On/Hi/Lo cps	3340.317	19.5	19.33333		
Beam/Abs current	3051.167	0			
Stage Position	7.343216	28.27236	10.95967	1	

SETUP.MDB  
 SETUP2.MDB (MAN)  
 SETUP3.MDB (Interf.)  
  

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time	138.8216	130.285	100	6	
Peakscan Hi/Lo/Points/Time	137.7545	131.3521	40	8	
Start/Stop/PB/Count/Attempts	0.1089777	5.448884E-03	3	10	30

Scroll through the list of elements (records) and find the desired element and setup using the data cursor. Use the left, right arrows (top, center) to move through the database. To see all the setups for a particular element, enter the element symbol into the *Search Element* text field and use the arrow keys as before. To see all element setups again, simply clear the *Search Element* text field. To view the most recent addition(s) to the SETUP database, click the  button on the data cursor.

Here, the user browses through the records and selects the appropriate silicon (si) entry as the first element setup to load. The output list order of elements will follow this list.

**Element Setup Database**

**Current Sample:**  
Un 1 setup for pyroxene analyses

Double click Analyzed Element List to see Element Setups

**<< Add to Sample**

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

**Close**

**Element Setup Data From SETUP.MDB Database**

◀◀ **Data Cursor** ▶▶

Delete from Database

Enter Search Element >>  Total Records = 21

<b>Elem/Order/CatOxd</b>	si	ka		1	2
<b>Spec/Crystal/2d</b>	1	PET	8.75	0.000144	
<b>User Name</b>	Dan Kremser				
<b>Sample Name</b>	Orthopyroxene				
<b>Date - Time</b>	3/16/1999 4:50:31 PM				
<b>Probe Data File</b>	E:\Probe Operators\Manual				
<b>On/Hi/Lo Pos</b>	228.2325	231.0645	225.3995		
<b>BgdType/Offset</b>	OFF	LIN	2.8320000000	-2.833	
<b>S-Hi/Lo/Exponen</b>	0	0	1		
<b>Base/Win/Gain/Bias</b>	0.5	10	50	1700	
<b>KeV/TO/DT/DIFF</b>	15	40	0.000001	0	
<b>Standard Intensity Data</b>					
<b>Std/PeakToBgd</b>	206	394.676			
<b>On/Hi/Lo sec</b>	20	4	4		
<b>On/Hi/Lo cps</b>	3288.967	6.083333	10.58333		
<b>Beam/Abs current</b>	3053.2	0			
<b>Stage Position</b>	23.88083	19.99947	10.99982	1	

SETUP.MDB  
 SETUP2.MDB (MAN)  
 SETUP3.MDB (Interf.)  
  

**Wavescan and Peaking Parameters**

<b>Wavescan Hi/Lo/Points/Time</b>	232.4814	223.9836	100	6	
<b>Peakscan Hi/Lo/Points/Time</b>	231.4192	225.0458	40	8	
<b>Start/Stop/PB/Count/Attempts</b>	0.1084835	5.424177E-03	3	10	30

Click the << **Add to Sample** button to add the element setup to the current sample.

The **Element Properties** window for silicon appears.

**Element Properties**

Enter Element Properties For: **si ka**

Element	X-Ray Line	Bragg Order	Cations / Oxygens	
si	ka	1	1	2

Leave the X-ray Line Blank to Indicate an Un-Analyzed Element (Specified, by Difference or Stoichiometry)

Charge: 4.000

Disable Acq  Disable Quant

---

**Parameters (note that Background Type can differ for Standards and Unknowns)**

**Background Type**

Off Peak  MAN  Multi-Point

**Off-Peak Entry**

Absolute Position  Relative Offset

Check All Interfering Elements

Hi Off-Peak Interferences  
Low Off-Peak Interferences

---

Spectrometer	Crystal	On-Peak	High Off-Peak	Low Off-Peak
1	PET	228.232	2.71600	-2.9490

BaseLine	Window	Gain	Bias	Deadtime (us)
.50	10.00	50.00	1700.	1.00

Use Differential PHA Mode

Slit Size	Slit Position	Detector Mode

**Integrated Intensity Scan**

Use Integrated Intensities  Use Inverted Intensity Steps

Initial Step Size	Minimum Step Size	Specified APF
.141625	.028325	1.00000

---

**Off Peak Correction Type**

Linear  Average  High Only  Low Only

Exponential	1.0000	Position1	Position2	Position3
		.000000	.000000	.000000

Slope (Hi)  Slope (Lo)  Polynomial

1.0000	1.0000	50.716	Coeff1	Coeff2	Coeff3
			.000000	.000000	.000000

---

Multi-Point

Acquire High	Iterate High	Acquire Low	Iterate Low
4	2	4	2

Fit Type: Linear

High Multi-Point Positions		Low Multi-Point Positions	
1	3.855011	1	-4.070999

Edit if required, then click the **OK** button to accept these values.

The silicon record is then listed in the text field under the previously defined sample name.

**Element Setup Database**

**Current Sample:**  
 Un 1 setup for pyroxene analyses  
 si ka Spectro 1 PET (228.232)

Double click Analyzed Element List to see Element Setups

**<< Add to Sample**

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

**Close**

**Element Setup Data From SETUP.MDB Database**

◀◀ **Data Cursor** ▶▶

Delete from Database

Enter Search Element >>

Total Records = 21

Elem/Order/CatOxd	fe	ka		1	1
Spec/Crystal/2d	3	LIF	4.0267	0.000058	

User Name: Dan Kremser

Sample Name: Fayalite

Date - Time: 3/16/1999 2:55:43 PM

Probe Data File: E:\Probe Operators\Manual Files\olivine.mdb

On/Hi/Lo Pos	134.5533	137.3983	131.7073	
BgdType/Offset	OFF	LIN	2.845	-2.846
S-Hi/Lo/Exponen	0	0	1	
Base/Win/Gain/Bias	1	10	200	1700
KeV/TO/DT/DIFF	15	40	0.000001	0

**Standard Intensity Data**

Std/PeakToBgd	203	172.0335		
On/Hi/Lo sec	20	4	4	
On/Hi/Lo cps	3340.317	19.5	19.33333	
Beam/Abs current	3051.167	0		
Stage Position	7.343216	28.27236	10.95967	1

SETUP.MDB  
  SETUP2.MDB (MAN)  
  SETUP3.MDB (Interf.)

Import   Export

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time	138.8216	130.285	100	6	
Peakscan Hi/Lo/Points/Time	137.7545	131.3521	40	8	
Start/Stop/PB/Count/Attempts	0.1089777	5.448884E-03	3	10	30

Continue browsing the element setup database and add all required element setups desired to the sample.



A typical pyroxene element setup list is shown below.

**Element Setup Database**

**Current Sample:**  
Un 1 setup for pyroxene analyses

Double click Analyzed Element List to see Element Setups

**<< Add to Sample**

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

**Close**

**Element Setup Data From SETUP.MDB Database**

◀◀ **Data Cursor** ▶▶

Delete from Database

Enter Search Element >>

Total Records = 21

Elem/Order/CatOxd	k	ka		2	1
Spec/Crystal/2d	1	PET	8.75	0.000144	

User Name: Dan Kremser

Sample Name: Microcline

Date - Time: 3/16/1999 4:43:13 PM

Probe Data File: E:\Probe Operators\Manual

On/Hi/Lo Pos: 119.7471 123.6551 115.8381

BgdType/Offset: OFF LIN 3.908 -3.9090000000

S-Hi/Lo/Exponen: 0 0 1

Base/Win/Gain/Bias: 0.5 10 50 1700

KeV/TO/DT/DIFF: 15 40 0.000001 0

**Standard Intensity Data**

Std/PeakToBgd	116	147.9095
On/Hi/Lo sec	20	4 4
On/Hi/Lo cps	2859.583	16.16667 22.5
Beam/Abs current	3045.2	0
Stage Position	4.122887 19.90012	11.12511 1

SETUP.MDB  SETUP2.MDB (MAN)  SETUP3.MDB (Interf.)

Import Export

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time	125.6102	113.8839	100	6
Peakscan Hi/Lo/Points/Time	124.1444	115.3497	40	8
Start/Stop/PB/Count/Attempts	0.1496983	7.484914E-03	3	10 30

Click the **Close** button when finished. The **New Sample** window reappears.

Click the **OK** button, returning to the **Acquire!** window. Don't forget to add oxygen as a specified element to the list for stoichiometry and other calculations. Note specified elements cannot be saved to the SETUP.MDB database.

## Sample Setups

Normally, PROBE FOR EPMA uses the sample setup from the last unknown (or standard if there are no unknown samples) to create the next new sample setup. Sample setups on the other hand are designed to allow the user to easily recall a previous sample setup within a current run. This allows the user to create and re-use multiple setups comprised of different groups of elements **within** a single run. In the example below, sample setups for pyroxene and olivine will be created, each with a different set of elements and conditions, that may be recalled at anytime during the current probe run.

The saving of a sample setup actually saves only a pointer to the sample selected. All of this sample's acquisition and calculation options, elements/cations, standard assignments, etc will be utilized when a new sample is created based on this sample setup. However, because counting time and the associated unknown count factor is treated by PROBE FOR EPMA as data as opposed to setup information, it is necessary (if the user wants this information to be carried over) to acquire at least one data point with the sample setup prior to saving it as a sample setup.

A new PROBE FOR EPMA run is opened in the usual manner. Ten elements and appropriate standards for pyroxenes are loaded from the SETUP.MDB database and the STANDARD.MDB database, respectively. Each element is then calibrated and standardized. Count times, acquisition and calculation options are adjusted to optimize the analyses and output requirements. And finally, a pyroxene standard is run once to update count time information.

To save the just calibrated pyroxene sample as a sample setup, start by clicking the **Add To Setup** button from the **Analyze!** window.

The screenshot shows the 'Analyze!' software interface. The 'Sample List (multi-select) (double-click to see intensity data)' is visible, with 'Un 2 Pyroxene Sample Setup' selected. The 'Standard Assignments' tab is active, showing analysis parameters for 'Un 2 Pyroxene #164 check'.

Analysis Parameters:

43.837	Total Oxygen	99.445	Total Weight %
43.837	Calculated Oxygen	12.363	Z - Bar
.000	Excess Oxygen	21.726	Atomic Weight

Results in Oxide Weight Percent:

Copy	SiO2	TiO2	Al2O3	Cr2O3	FeO	MnO	MgO	CaO	Na2O	K2O
Average:	49.723	.560	7.714	.914	4.641	.118	17.630	17.293	.846	.006
Std Dev:	.000	.000	.000	.000	.000	.000	.000	.000	.000	.000
ZAF Corr:	1.2835	1.1938	1.4370	1.1798	1.1910	1.2113	1.4777	1.0865	1.9293	1.0956
Std Err:	.000	.000	.000	.000	.000	.000	.000	.000	.000	.000
%Rel SD:	.00	.00	.00	.00	.00	.00	.00	.00	.00	.00
Minimum:	49.723	.560	7.714	.914	4.641	.118	17.630	17.293	.846	.006
Maximum:	49.723	.560	7.714	.914	4.641	.118	17.630	17.293	.846	.006

The bottom table shows the current analysis line:

Copy	SiO2	TiO2	Al2O3	Cr2O3	FeO	MnO	MgO	CaO	Na2O	K2O
10 G	49.723	.560	7.714	.914	4.641	.118	17.630	17.293	.846	.006

To review the saved sample setup, click the **Elements/Cations** button.

The **Analyzed and Specified Elements** dialog box appears.

**Analyzed and Specified Elements**

**Selected Samples**

Un 2 Pyroxene Sample Setup

**OK** **Cancel**

**Save Element Setup**

**Save Sample Setup**

**Add/Remove Standards**

Reload Standard Assignments

Remove TDI Correction


1 2 3 4 5 6

**Click Element Row to Edit Element/Cations Parameters (click empty row to add)**

Channel	Element	X-Ray	Analyzed	Cations	Oxygens	Off-Peak
1	si	ka	Yes	1	2	Linear
2	ti	ka	Yes	1	2	Linear
3	al	ka	Yes	2	3	Linear
4	cr	ka	Yes	2	3	Linear
5	fe	ka	Yes	1	1	Linear
6	mn	ka	Yes	1	1	Linear
7	mg	ka	Yes	1	1	Linear
8	ca	ka	Yes	1	1	Linear
9	na	ka	Yes	2	1	Linear
10	k	ka	Yes	2	1	Linear
11	o		No	1	0	

Click the **Save Sample Setup** button.


The **Sample Setup Selection** window opens.



Edit the *Sample Setup Description* text box as desired.

Click the << **Add To Sample Setups** button.


The **Sample Setup Selection** window appears as below. Note: the first number (12) represents the sample's row number and can be seen listed using the **Run | List Sample Rows, Names, Conditions** menu.



Click the **OK** button, returning to the **Analyzed and Specified Elements** window.

Click the **OK** button of the **Analyzed and Specified Elements** window returning to the **Analyze!** window.

Return to the **Acquire!** window to create a new sample.



Click the **New Sample** button.

Edit the *New Sample Name* text field. Here, the user will establish an olivine sample setup.

**New Sample**

**New Sample Type**

Standard

Unknown

Wavescan

Add/Remove Standards

OK Cancel

Load Element Setups

Load Sample Setup

Load File Setup

Load Multiple Setup

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

Load Wavescan From Another Probe Run

**New Sample Name**

setup for olivine analyses

**New Sample Description** Add <cr>

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.

81 Albite

116 Microcline

203 Fayalite

204 Forsterite

206 Orthopyroxene

207 Kyanite

210 Wollastonite

Several paths may be taken from here to load new elements for the olivine sample. To enter an entirely new list of elements and parameters it might be easier to click the **OK** button and follow the **Element/Cations** button of the **Acquire!** window to the **Element Properties** dialog.

If (as in this example) only minor changes to a sample are required then from the **New Sample** window, click the **Load Element Setups** button.



The **Element Setup Database** opens.

**Element Setup Database**

**Current Sample:**  
 Un 2 Pyroxene Sample Setup

si ka Spectro 1 PET (228.132)  
 ti ka Spectro 3 LIF (191.114)  
 al ka Spectro 2 TAP (89.5182)  
 cr ka Spectro 3 LIF (159.211)  
 fe ka Spectro 3 LIF (134.634)  
 mn ka Spectro 3 LIF (146.149)  
 mg ka Spectro 2 TAP (106.846)  
 ca ka Spectro 1 PET (107.534)  
 na ka Spectro 2 TAP (129.282)  
 k ka Spectro 1 PET (119.786)  
 o (specified)

Double click Analyzed Element List to see Element Setups

**<< Add to Sample**

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

**Close**

**Element Setup Data From SETUP.MDB Database**

Enter Search Element >>  Total Records = 21

Elem/Order/CatOxd	fe	ka		1	1
Spec/Crystal/2d	3	LIF	4.0267	0.000058	

User Name: Dan Kremser

Sample Name: Fayalite

Date - Time: 3/16/1999 2:55:43 PM

Probe Data File: E:\Probe Operators\Manual Files\olivine.mdb

On/Hi/Lo Pos: 134.5533 137.3983 131.7073

BgdType/Offset: OFF LIN 2.845 -2.846

S-Hi/Lo/Exponen: 0 0 1

Base/Win/Gain/Bias: 1 10 200 1700

KeV/TO/DT/DIFF: 15 40 0.000001 0

**Standard Intensity Data**

Std/PeakToBgd	203	172.0335		
On/Hi/Lo sec	20	4	4	
On/Hi/Lo cps	3340.317	19.5	19.33333	
Beam/Abs current	3051.167	0		
Stage Position	7.343216	28.27236	10.95967	1

SETUP.MDB
  SETUP2.MDB (MAN)
  SETUP3.MDB (Interf.)

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time	138.8216	130.285	100	6
Peakscan Hi/Lo/Points/Time	137.7545	131.3521	40	8
Start/Stop/PB/Count/Attempts	0.1089777	5.448884E-03	3	10 30

Edit the previous pyroxene list, in this example chromium, sodium, potassium are eliminated from the list by highlighting each element and clicking the **Delete from Sample** button. If additional elements are required, recall them at this time (nickel is added in this example).

After editing, the window appears as below.

**Element Setup Database**

**Current Sample:**  
Un 2 Pyroxene Sample Setup

si ka Spectro 1 PET (228.132)  
 ti ka Spectro 3 LIF (191.114)  
 al ka Spectro 2 TAP (89.5182)  
 fe ka Spectro 3 LIF (134.634)  
 mn ka Spectro 3 LIF (146.149)  
 mg ka Spectro 2 TAP (106.846)  
 ca ka Spectro 1 PET (107.534)  
 ni ka Spectro 3 LIF (115.190)  
 o (specified)

Double click Analyzed Element List to see Element Setups

**<< Add to Sample**

Delete from Sample

Delete All Elements From Sample

Add To Database >>

Load Standard Intensity From Database To Current Run

**Close**

---

**Element Setup Data From SETUP.MDB Database**

Enter Search Element >>  Total Records = 21

Elem/Order/CatOxd	ni	ka		1	1
Spec/Crystal/2d	3	LIF	4.0267	0.000058	

User Name: Dan Kremser

Sample Name: Ni-Olivine

Date - Time: 3/16/1999 3:36:47 PM

Probe Data File: E:\Probe Operators\Manual Files\olivine.mdb

On/Hi/Lo Pos	115.1904	118.2654	112.1164
BgdType/Offset	OFF	LIN	3.075
S-Hi/Lo/Exponen	0	0	1
Base/Win/Gain/Bias	1	10	200
KeV/TO/DT/DIFF	15	40	0.000001

**Standard Intensity Data**

Std/PeakToBgd	222	126.0353
On/Hi/Lo sec	10	2
On/Hi/Lo cps	3213.9	27
Beam/Abs current	3056.467	0
Stage Position	3.937658	27.15654

SETUP.MDB  
  SETUP2.MDB (MAN)  
  SETUP3.MDB (Interf.)  
  


---

**Wavescan and Peaking Parameters**

Wavescan Hi/Lo/Points/Time	119.8027	110.5782	100	6
Peakscan Hi/Lo/Points/Time	118.6496	111.7313	40	8
Start/Stop/PB/Count/Attempts	0.1177592	5.887959E-03	3	10

Click the **Close** button of the **Element Setup Database**, returning to the **New Sample** window.

Click the **OK** button in the **New Sample** window, returning to the **Acquire!** window.



If different standard choices are required they should be added from the STANDARD.MDB database at this point. Use the **Standard | Add/Remove Standards To/From Run** menu in the main PROBE FOR EPMA log window.

Recalibrate and standardize all new elements and adjust count times, acquisition, and calculation options to optimize for olivine analysis. Run an olivine standard to obtain data and check the calibration.

From the **Analyze!** window, click the **Add To Setup** button.

The screenshot shows the 'Analyze!' software window. The 'Sample List (multi-select) (double-click to see intensity data)' section is active, with 'Un 4 Olivine Sample Setup' selected. The 'Add To Setup' button is highlighted. Below this, the 'Standard Assignments' tab is selected, showing parameters for 'Un 4 Olivine #117 check' (TO = 40, KeV = 15, Beam = 30, Size = 10, Mode = Analog Spot). A table displays 'Results in Oxide Weight Percent' for elements SiO2, TiO2, Al2O3, FeO, MnO, MgO, CaO, NiO, and O. The 'Average' row shows values: SiO2: 40.686, TiO2: .011, Al2O3: .017, FeO: 9.584, MnO: .097, MgO: 49.227, CaO: .088, NiO: .379, O: .000, Total: 100.087. Below the table are buttons for 'Delete Selected Line(s)', 'Undelete Selected Line(s)', and 'Analyze Selected Line(s)'. A second table shows the '14 G' sample with the same oxide weight percentages as the average. The 'Cancel' and 'Next' buttons are at the bottom right.

Copy	SiO2	TiO2	Al2O3	FeO	MnO	MgO	CaO	NiO	O	Total
Average:	40.686	.011	.017	9.584	.097	49.227	.088	.379	.000	100.087
Std Dev:	.000	.000	.000	.000	.000	.000	.000	.000	.000	.000
ZAF Corr:	1.3962	1.1682	1.6831	1.1851	1.2042	1.4293	1.0939	1.1850		
Std Err:	.000	.000	.000	.000	.000	.000	.000	.000	.000	.000
%Rel SD:	.00	.00	.00	.00	.00	.00	.00	.00	.00	.00
Minimum:	40.686	.011	.017	9.584	.097	49.227	.088	.379	.000	100.087
Maximum:	40.686	.011	.017	9.584	.097	49.227	.088	.379	.000	100.087


  

Copy	SiO2	TiO2	Al2O3	FeO	MnO	MgO	CaO	NiO	O	Total
14 G	40.686	.011	.017	9.584	.097	49.227	.088	.379	.000	100.087

Click the **Elements/Cations** button.

The **Analyzed and Specified Elements** window opens, click the **Save Sample Setup** button.

The **Sample Setup Selection** window appears. The Olivine Sample Setup has been added to *Available Sample Setups* list box. If the user wishes to modify the *Sample Setup Description* either change the label using the **Name/Description** button in the **Analyze!** window before clicking the **Add To Sample Setups** button or from the **Sample Setup Selection** window, delete the just loaded selection (highlight the appropriate selection) and then edit the *Sample Setup Description* text field and click the << **Add To Sample Setups** button, storing the Olivine setup along with the previously stored Pyroxene setup.



Click the **OK** button of the **Sample Setup Selection** window returning to the **Analyzed and Specified Elements** window.

Click the **OK** button to go back to the **Analyze!** window.

Any number of sample setups can be created as described above.

The user now has two calibrated sample setups available to analyze any pyroxene or olivine in the samples supplied for microprobe analysis. The olivine setup (last) is currently active however to recall any other sample setup, follow the steps outlined below.

Bring forward the **Acquire!** window. Move to the next unknown analysis spot, in this example the user wishes to analyze several pyroxene grains.

The screenshot shows the 'Acquire!' software window. At the top, there are columns labeled 1, 2, 3, X, Y, Z, W. Below these are several rows of numerical data. A 'Current Sample' field contains 'Un 4 Olivine Sample Setup'. A prominent yellow button reads 'Start Standard or Unknown Acquisition'. Below this is a 'Start Wavescan' button. The bottom half of the window is a grid of buttons for various acquisition options. On the right side, there is a 'Progress/Beam Deflection' section with numerical values for X and Y coordinates, and a table of instrument parameters.

1	2	3	X	Y	Z	W
228.128	89.5225	115.280	3.44996	27.5378	10.9485	1.00000
<b>1-PET</b>	<b>2-TAP</b>	<b>3-LIF</b>	<b>Faraday</b>			
4.00	4.00	4.00	1.00			
91.	95.	181.	3000.01			

Current Sample: Un 4 Olivine Sample Setup

Normal Acquisition Unknown

Data Rows: 1    Good Data Rows: 1

**Start Standard or Unknown Acquisition**

Start Wavescan

New Sample	PHA	Move	Stage
Elements/Cations	Peak/Scan Options	Acquisition Options	Imaging
Analytical Conditions	Count Times	Special Options	Peaking Options
Combined Conditions	Locate	Rate Meter	Start Peaking

Progress/Beam Deflection

	3.44996	27.5378
um	.000000	.000000
px	0	0

Magnification: 400  
 Beam Mode: Analog Spot  
 Kilovolts: 15  
 Beam Current: 30  
 Beam Size: 10  
 Acquire String:

Click the **New Sample** button.

The **New Sample** window opens. Enter the appropriate text into the *New Sample Name* and *New Sample Description* fields.

**New Sample**

**New Sample Type**

Standard

Unknown

Wavescan

**Add/Remove Standards**

**OK** **Cancel**

**Load Element Setups**

**Load Sample Setup**

**Load File Setup**

**Load Multiple Setup**

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

**Load Wavescan From Another Probe Run**

**New Sample Name**

Pyroxene S2-123 C1-1

**New Sample Description** **Add <cr>**

Stillwater

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.

81 Albite

116 Microcline

203 Fayalite

204 Forsterite

206 Orthopyroxene


207 Kyanite

210 Wollastonite

Click the **Load Sample Setup** button.

This opens the **Sample Setup Selection** window.

Select the Pyroxene setup, highlighting it allows the operator to view the element list.




Click the **OK** button of the **Sample Setup Selection** window to load the sample setup.

The program returns to the **New Sample** window. Click the **OK** button.



The **Acquire!** window reappears.




Double check your spot selection and focus and click the **Start Standard or Unknown Acquisition** button to initiate data acquisition.

The availability of multiple sample setups during the course of automated unknown analysis gives the user tremendous flexibility. Upon activation of the **Use Digitized Sample Setups** button in the **Automate!** window, each unknown analysis may be based on a different sample setup that was specified when the unknown sample position was digitized. See the User's Guide and Reference documentation for more details.

# File Setups

To load any sample setup from a previous probe run file, the file setup option is provided. These file setups are old Probe database files that contain old sample setups and may or may not contain standardization count intensity data.

The example below will illustrate how to use the file setup option to easily import two different (an olivine and a pyroxene) sample setups into the current new probe run. Open a new PROBE FOR EPMA run and click the **New Sample** button from the **Acquire!** window.



The **New Sample** window appears. Edit the *New Sample Name* text box.

**New Sample**

**New Sample Type**

Standard

**Unknown**

Wavescan

**Add/Remove Standards**

**OK** **Cancel**

**Load Element Setups**

Load Sample Setup

**Load File Setup**

Load Multiple Setup

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

Load Wavescan From Another Probe Run

**New Sample Name**


olivine setup

**New Sample Description** **Add <cr>**

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.


Click the **Load File Setup** button.

The **Load File Setup** window opens and will list all available PROBE FOR EPMA files that can be loaded. The initial available Probe Run Files directory pointer is the location specified when opening a new probe database file earlier. Move to another directory location if necessary. The last file listed in the available Probe Run Files along with the last entry in the *Samples List* will be shown by default.




In this example, the file STILLWATER OL-PX.MDB was newly created and the user will load both an olivine and a pyroxene sample setup into this file.

Scroll through the available Probe Run Files list and highlight the file to load from. The last sample setup will be displayed in the *Samples List* and *Element List* text field. Next, select the sample setup that you wish to load into the new probe run. All of the run parameters and options for that sample setup will be loaded. The only parameters not loaded are the nominal beam current and the volatile element assignments since they are unknown sample specific.




Click the **OK** button to load in this sample setup of interest.

The **NewFileSetupLoadStandardData** window appears next, asking whether the user wants the previous standard intensity data to be loaded as well.



Selecting **Yes** would load the old standard intensity data from the file setup into this new run. Depending on the stability of your instrument, it may or may not be necessary to re-standardize some or all of the standards. In this case, the user chooses to load the standard intensity data, selecting the **Yes** button.

The **NewFileSetupLoadStandardData** window appears.



Click the **OK** button.

The **New Sample** dialog box reappears.

**New Sample**

**New Sample Type**

Standard

Unknown

Wavescan

**Add/Remove Standards**

**OK** **Cancel**

**Load Element Setups**

**Load Sample Setup**

**Load File Setup**

**Load Multiple Setup**

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

**Load Wavescan From Another Probe Run**

**New Sample Name**

olivine setup

**New Sample Description** **Add <cr>**

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.

Click the **OK** button to complete the loading of the olivine sample setup from the old probe run.

The program now returns to the fully active **Acquire!** window.

The screenshot shows the 'Acquire!' software window with the following data and controls:

1	2	3	X	Y	Z	W
228.128	89.5225	115.280	3.44996	27.5378	10.9485	1.00000

1-PET	2-TAP	3-LIF	Faraday
.00	.00	.00	.00
.	.	.	.000000

Current Sample: Un 1 \* olivine setup

Normal Acquisition Unknown

Data Rows: 0      Good Data Rows: 0

Buttons: **New Sample**, PHA, Move, Stage, Elements/Cations, Peak/Scan Options, Acquisition Options, Imaging, Analytical Conditions, Count Times, Special Options, Peaking Options, Combined Conditions, Locate, Rate Meter, Start Peaking

Start Standard or Unknown Acquisition (highlighted)

Start Wavescan

Progress/Beam Deflection

	3.44996	27.5378
um	.000000	.000000
px	0	0

Magnification: 400

Beam Mode: Analog Spot

Kilovolts: 15

Beam Current: 30


Beam Size: 10

Acquire String:



Normally the user would check the calibration by running a secondary standard or two to verify the composition, repeaking and/or collecting standard intensities as required.

The user then opens the **Analyze!** window to save this olivine setup as a sample setup in this current probe run.




Click the **Add To Setup** button.

Click the **Elements/Cations** button, opening the **Analyzed and Specified Elements** window.

Click the **Save Sample Setups** button.


The **Sample Setup Selection** window opens. Edit the *Sample Setup Description* text box and click the << **Add To Sample Setups** button, resulting in the following window.




Click the **OK** button, returning to the **Analyzed and Specified Elements** window. Click this **OK** button to return to the **Analyze!** window.

If another, previously created sample setup is needed for this current probe run, open the **New Sample** window and follow the instructions of the past eight pages.

Remember to save each sample setup in the **Sample Setup Selection** window as described above.



Returning to the **Acquire!** window, the user can now employ either sample setup for probe work.




# Wavescans

## Quick Wavescan Acquisition

This feature is useful if an EDS detector is not available or WDS resolution over the entire spectrometer range is required. The program will move each spectrometer currently assigned to it's upper limit and then continuously scan each spectrometer to it's lower travel limit while acquiring simultaneous count data. The count time used for the **Quick Wavescan Acquisition** is specified in the **Count Times** dialog box, opened from the **Acquire!** window. The current sample setup specifies which spectrometer and reflecting crystal to use. The program uses the spectrometer calibration of the first acquired element (order = 1) in the sample.

From an open PROBE FOR EPMA run, containing an unknown sample and the appropriate unknown under the crosshairs, click the **Special Options** button from the **Acquire!** window.



This opens the **Volatile Calibration and Quick Wavescan Samples** window. Note the default acquisition option is *Normal Acquisition*.

**Volatile Calibration and Quick Wavescan Samples**

**Special Sample Acquisition Options**

**Normal Acquisition**  
Acquire a normal standard, unknown or wavescan sample.

**Self Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample that uses a calibration curve based on itself. The TDI assignment will be made automatically. This method will append a new analysis line for each acquisition on an existing unknown sample.  
**Consider using "synchronous" spectrometer acquisition for TDI to minimize beam exposure before starting count integration.**

Note: time dependent intensity (TDI) elements must be the first element acquired on each spectrometer. The acquisition order may be modified in the Acquisition Options dialog.

**Time Dependent Intensity (TDI) Count Time Intervals**   
 Acquire TDI Data on Standard Samples

**Assigned Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample for use as a TDI calibration curve for another sample(s). This assignment is made subsequently in the Standard Assignments dialog in the Analyze! window. This method will create a new time dependent intensity calibration sample for each acquisition.

**Time Dependent Intensity (TDI) Sample Name**

**Time Dependent Intensity (TDI) Count Time Intervals**   
**Stage X Increment (um)**   
**Stage Y Increment (um)**

**Quick Wavescan Acquisition**  
Quick wavescan samples are created based on the first element on each spectrometer used in the last unknown sample (or wavescan sample if present) or a selected sample setup. This method will create a new quick wavescan sample for each acquisition.

**Quick Wavescan Sample Name**

**Quick Scan Speed % (.001 - 100)**   
 Use ROM Based Spectrometer Scan

**OK**  
**Cancel**  
**Setup**

Select the *Quick Wavescan Acquisition* dialog button. Enter a *Quick Wavescan Sample Name* and *Quick Scan Speed* into the text fields. The smaller the scan speed percentage the slower the spectrometer will travel per second and of course each instrument would require different settings.

**Volatile Calibration and Quick Wavescan Samples**

**Special Sample Acquisition Options**

**Normal Acquisition**  
Acquire a normal standard, unknown or wavescan sample.

**Self Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample that uses a calibration curve based on itself. The TDI assignment will be made automatically. This method will append a new analysis line for each acquisition on an existing unknown sample.

Note: time dependent intensity (TDI) elements must be the first element acquired on each spectrometer. The acquisition order may be modified in the Acquisition Options dialog.

**Time Dependent Intensity (TDI) Count Time Intervals**

Acquire TDI Data on Standard Samples

**Assigned Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample for use as a TDI calibration curve for another sample(s). This assignment is made subsequently in the Standard Assignments dialog in the Analyze! window. This method will create a new time dependent intensity calibration sample for each acquisition.

**Time Dependent Intensity (TDI) Sample Name**

**Time Dependent Intensity (TDI) Count Time Intervals**

**Stage X Increment (um)**

**Stage Y Increment (um)**

**Quick Wavescan Acquisition**  
Quick wavescan samples are created based on the first element on each spectrometer used in the last unknown sample (or wavescan sample if present) or a selected sample setup. This method will create a new quick wavescan sample for each acquisition.

**Quick Wavescan Sample Name**

**Quick Scan Speed % (.001 - 100)**

Use ROM Based Spectrometer Scan

**Consider using "synchronous" spectrometer acquisition for TDI to minimize beam exposure before starting count integration.**


**OK**

**Cancel**

**Setups**

Click the **OK** button to return to the **Acquire!** window.

To initiate the quick wavescan acquisition, click the **Start Wavescan** button in the **Acquire!** window. A new wavescan sample is automatically started using the sample name just supplied. The spectrometers move to their respective upper limits and proceed with the wavescan. The **Wavescan Acquisition** window opens and real time data display is viewable. A completed three-spectrometer **Wavescan Acquisition** window appears below.



The size of each graph maybe expanded (not shown) by clicking on the relevant wavescan. Upon completion of the quick wavescan, the data may be exported via the **Export Data** button to an ASCII file or examined in more detail along with KLM marker overlay capabilities from the **Plot!** window. Printing of the quick wavescan is possible by selecting the **Print** option under the **Graph Data** window (see next section for a specific example).

## Calibrated Multi-Element Wavescans

Another unique feature of PROBE FOR EPMA is the ability to acquire calibrated multi-element wavescans. This provides an easy and rapid method to scan **all** elements in a sample for off-peak interferences. The example below will illustrate calibrated wavescans on a ten-element pyroxene sample and the adjustment of off-peak background positions.

Open a new PROBE FOR EPMA run in the usual manner. Confirm motor and crystal positions as well as setting the beam current to the appropriate value. Click the **New Sample** button and create a sample using the elements of interest. Next, re-peak the elements using either manual or automatic peaking on the appropriate standards. This calibrates the spectrometer motors. And finally, move to the sample to perform the calibrated wavescan.

From the **Acquire!** window, click the **New Sample** button to create a wavescan sample.

The screenshot shows the 'Acquire!' software interface. At the top, there are columns labeled 1, 2, 3, X, Y, Z, and W. Below these, a row of numerical values is displayed: 228.140, 89.5139, 191.106, 3.44996, 27.5378, 10.9485, 1.00000. Below this, there are columns for 1-PET, 2-TAP, 3-LIF, and Faraday, with values .00, .00, .00, and .00 respectively. A row below shows a value of .000000. The 'Current Sample' field contains 'Un 1 \* template for wavescan'. A yellow button labeled 'Start Standard or Unknown Acquisition' is prominent. Below it is a 'Start Wavescan' button. The interface includes a grid of buttons for 'New Sample', 'PHA', 'Move', 'Stage', 'Elements/Cations', 'Peak/Scan Options', 'Acquisition Options', 'Imaging', 'Analytical Conditions', 'Count Times', 'Special Options', 'Peaking Options', and 'Combined Conditions', 'Locate', 'Rate Meter', and 'Start Peaking'. On the right side, there is a 'Progress/Beam Deflection' section with a red crosshair and numerical values for X (3.44996) and Y (27.5378) coordinates, and um and px units. Below this, a table lists acquisition parameters: Magnification (400), Beam Mode (Analog Spot), Kilovolts (15), Beam Current (40), and Beam Size (10). The 'Acquire String' field is empty.

1	2	3	X	Y	Z	W
228.140	89.5139	191.106	3.44996	27.5378	10.9485	1.00000
1-PET	2-TAP	3-LIF	Faraday			
.00	.00	.00	.00			
			.000000			

um	3.44996	27.5378
px	0	0

Magnification	400
Beam Mode	Analog Spot
Kilovolts	15
Beam Current	40
Beam Size	10
Acquire String	



The **New Sample** window opens. Select the *Wavescan* check button as the *New Sample Type*. Edit the *New Sample Name* and *New Sample Description* text fields.

**New Sample**

**New Sample Type**

Standard

Unknown

**Wavescan**

**Add/Remove Standards**

**OK** **Cancel**

**Load Element Setups**

**Load Sample Setup**

**Load File Setup**

**Load Multiple Setup**

Note that a new standard sample element setup is based by default on the last unknown sample in the run. To change the analyzed elements in a run, either click one of the Load Setup buttons above or first create a new unknown sample and then make any necessary changes to the element setup.

**Load Wavescan From Another Probe Run**

**New Sample Name**


Wavescan on Hedenbergite #66

**New Sample Description** **Add <cr>**

To add standards to the standard list below, cancel this dialog, then click the Standard | Add Standards to Run menu item from the main menu.

Click the **OK** button.

The program returns to the **Acquire!** window.



To modify the wavescan range and/or number of data points to be collected, click on the **Peak/Scan Options** button in the **Acquire!** window.

Select the *Wave Scan Limits* check button under *Display:* and click on the appropriate element row to edit the parameters. The stage may also be moved (incremented) during the acquisition using the *Stage Step During Peakscan/Wavescan or Peaking (X and Y Axis)* check box and *Increment Size (microns)* text field.

**Peak and Scan**

Click Element Row to Edit Peak and Scan Parameters

Channel	Element	Spectro	Crystal	Wave-Hi	Wave-Lo	Wave-Pts	Hi-Off	Lo-Off
1	si ka	1	PET	232.479	223.985	100	4.34789	-4.1470
2	ti ka	3	LIF	194.704	187.536	100	3.58954	-3.5777
3	al ka	2	TAP	99.8831	80.5209	100	10.3649	-8.9972
4	cr ka	3	LIF	163.062	155.208	100	3.85172	-4.0029
5	fe ka	3	LIF	138.824	130.282	100	4.18993	-4.3521
6	mn ka	3	LIF	150.248	142.052	100	4.09908	-4.0970
7	mg ka	2	TAP	116.179	98.4268	100	9.33347	-8.4188
8	ca ka	1	PET	113.696	101.318	100	6.16216	-6.2153
9	na ka	2	TAP	137.608	121.450	100	8.32600	-7.8316
10	k ka	1	PET	130.000	110.000	100	10.2140	-9.7860

Display :

On/Off Peaks  
 **Wave Scan Limits**  
 Peak Scan Limits  
 Peaking Parameters

Spectrometer "Offsets" are the difference between the theoretical or predicted peak position (from x-ray tables) and the actual or measured peak position. If the spectrometer peak position has not had a peak center procedure performed then the "Offsets" value will usually be close to zero. The calculation is "Offset = Predicted - Actual"

Use ROM Based Spectrometer Scanning

Increment Stage During Peakscan/Wavescan or Peaking (X and Y Axis)

Use Increment During Scanning  
 Use Increment During Peaking

Increment Size (microns)  
 X  Y

Increment Interval (seconds)

OK Cancel

The user wishes to adjust the spectrometer start and stop values for k ka, click on the row of that element.

This opens the **Peak and Scan Properties** window. Adjust the appropriate values.

**Peak and Scan Properties**

Enter Peak and Scan Properties For: k ka 1 PET

	Spectrometer	Angstrom
On Peak	119.786	3.74243
Hi Off Peak	123.694	3.86454
Low Off Peak	115.877	3.62029
Wavescan Hi Limit	130.000	4.06157
WaveScan Low Limit	110.000	3.43666
Wavescan Points	100	.202020 per step
Peakscan Hi Limit	124.145	3.87863
Peakscan Low Limit	115.349	3.60379
Peakscan Points	40	.225538 per step
Peaking/ROM Start Size	.112287	.003508
Peaking Stop Size	.007486	.000234
Minimum P/B	3.00	
Minimum Peak Counts (cps)	10.0	
Maximum Peaking Attempts	30	

**Enter Positions in**

Spectrometer Units  
 Angstrom Units

**Display Positions in**

Absolute Position  
 Relative Offset

The on-peak element in a wavescan sample is simply there for spectrometer calibration purposes. Specifically the spectrometer unit to angstrom conversion. However, the high and low scan limits may be arbitrarily set to any value in the total spectrometer range.

**Set Wavescan Full Range**  
**Set Wavescan Normal**  
**Move To On Peak**  
**Set On Peak To Default**  
**Set On Peak To Current**

Click the **OK** button of the **Peak and Scan Properties** window when done editing.

Then click the **OK** of **Peak and Scan** window to close.

Wavescan count times for each element are adjusted via the **Count Times** button in the **Acquire!** window.

**Count Times**

Click Element Row to Edit Count Times

Channel	Element	Spectro	Crystal	On-Peak	Hi-Peak	Lo-Peak	MaxCount	Factor	Wave	Peak	Quick
1	si ka	1	PET	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
2	ti ka	3	LIF	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
3	al ka	2	TAP	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
4	cr ka	3	LIF	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
5	fe ka	3	LIF	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
6	mn ka	3	LIF	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
7	mg ka	2	TAP	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
8	ca ka	1	PET	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
9	na ka	2	TAP	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50
10	k ka	1	PET	20.00	4.00	4.00	10000000	1.00	4.00	8.00	.50

Beam Averages:

Nominal Beam (nA):

Change the Nominal Beam to modify the normalization constant used for the x-ray intensity display. For example, enter 1 (nA) for cps/nA intensity display.

Return To On-Peak Time:

Crystal Flip Time:

Set Column (TKCS) Time:

2083 secs  
34 min

Calculated Spectrometer Motion and Acquisition Time

0 secs

1	2	3
PET	TAP	LIF
K	Na	Mn
Ca	Mg	Fe
Si	Al	Ti

OK


Cancel

Measure Nominal Beam

Click on the appropriate element row to edit the wavescan time. Edit the **Count Time Properties** dialog box and then click it's **OK** button to close.

Click the **OK** button to close the **Count Times** window.

Click the **Start Wavescan** button in the **Acquire!** window to initiate the calibrated multi-element wavescan. The **Wavescan Acquisition** window opens. The program will automatically start acquiring the wavescan ranges selected. If more than one element is assigned to a given spectrometer, the program will automatically go to the next element's wavescan range after the previous wavescan element range is completed. The order of acquisition is defined in the **Acquisition Options** window. Below illustrates a completed wavescan acquisition.




As the wavescan is acquiring data, the wavescan graph may be viewed in greater detail by clicking on the graph to toggle/expand the display size.

The *Position* (spectrometer units), *Angstroms*, and *Counts* in any channel may be read by placing the cursor on the graph. Selecting the *Load Xray* check box and clicking the graph, loads the NIST x-ray database.


After all wavescans have been acquired on the sample, the user would typically inspect off-peak interferences and background locations by using the **Plot!** window. Note that if more than 100 points were acquired in a wavescan be sure to highlight all of the “continued” samples associated with the wavescan.

Select an *X-Axis* parameter (normally a specific spectrometer) and a *Y-Axis* parameter (normally the associated wavescan counts). The number (X) after the element in each *List* designates the spectrometer employed to collect the data. Finally, click the *Line* check button under *Graph Type*.



Click the **Output** button to graph the wavescan.


The **Plot Graph Data** window opens displaying the plotted components. The currently selected off-peak positions for background measurements are also indicated (green).



To evaluate potential interferences select a *KLM Markers* option (*Analyzed Elements* check button, for instance) to view the KLM markers or use the **Load Xray Database** button.




Click and drag the mouse to **Zoom** in on any portion of the graph. The **Plot Graph Data** window open below illustrates this powerful feature and the identification of the small x-ray peaks (satellite lines) to the high-energy side of the main silicon x-ray peak.




The off-peak positions for background determinations for quantitative samples are adjusted with the **Low** and **High** buttons (located lower right of **Plot Graph Data** window).

In the next screen capture, detailing *K (1) Wavescan Counts (PET cps) versus K (1) Spectrometer*, the user will note that the high background position falls on top of the signal from Fe  $K\alpha_{1,2}$  second order. This could give an anonymously high background reading on this sample, if iron is present. Therefore, the user chooses to move the high background slightly to a lower spectrometer position.




Click the **High** button, move the mouse cursor (note it appears as a cross on the graph) to an appropriate new position and click the graph.

The graph will update the new off-peak position (purple is the old position, green the new position).



Click the **OK** button of the **Graph Data** window. The **GetPeakSave** window appears.

The new parameters will take effect on the next new sample that is started




Click the **OK** button to accept and store the new high background position for K K $\alpha$  x-rays.


## Polygon Gridding – Using Surfer Option

Another useful feature of PROBE FOR EPMA is the ability to perform automated polygon gridded analyses of unknowns. After acquiring the digitized data set, PROBE FOR EPMA can create a script file (if the SURFER.BAS file option is selected in the **Plot!** window) for use with SURFER for Windows to automatically generate contour, surface and \*.GRD concentration files of your data. These \*.GRD files can be imported into Probe Image for viewing in false color. The images will be quantitatively registered during the import process so that color represents elemental or oxide concentration.


In this example an unknown and complexly exsolved pyroxene (see image below) will be gridded and digitized, then run quantitatively. Move to the unknown grain location.



Click the **Digitize** button of the **Automate!** window.




The **Digitize Sample Positions** window opens.



Select the *Unknown* check button from the *Sample Type* choices.


Enter a new sample name in the *Unknown or Wavescan Position Samples* text field, and click the **Add New Unknown To Position List** button.

A digitized polygon area grid will now be setup on the unknown grain.



Click the **Polygon Grid** button at the bottom of the **Digitize Sample Positions** window.

The **Polygon Grid Parameters** window opens.



The user will outline the perimeter of the grain to be gridded. An easy way to accomplish this is to image the grain with backscattered electrons, at any magnification, and trace around the grain boundary. Start in one corner and on a recognizable feature, click the **Add Polygon Boundary Coordinate** button and then move linearly toward another feature or edge, clicking the **Add Polygon Boundary Coordinate** button to outline this portion of the grain. Continue to trace line segments around the grain, clicking the **Add Polygon Boundary Coordinate** button to enclose



another portion of the grain. Eventually, returning to the starting point, completing the enclosure.

In this example, twenty line segments were used to enclose the grain of interest. Each end point is listed in the *Polygon Grid Stage Coordinate Boundary Points* text box. If a mistake is made or you simply wish to remove the previous boundary point, click the **Remove Last Polygon Position** button.

**Polygon Grid Parameters**

**Polygon Grid Stage Coordinate Boundary Points**

12	27.0759	31.5605	11.1451
13	27.2444	31.5994	11.1171
14	27.0605	31.8654	11.1419
15	27.0567	31.0791	11.1420
16	27.3036	31.8187	11.1434
17	27.1575	32.0145	11.1472
18	27.4033	31.8310	11.1011
19	27.1011	31.2082	11.1176
20	27.1516	31.9513	11.1189

**Add Polygon Boundary Coordinate**

**Plot Polygon Boundary**      **Remove Last Polygon Position**

**OK**      **Cancel**

**Test**      **Load**

**Polygon Grid Z Interpolation Options**

Use Averaged Z Position of Polygon Boundary Coordinates  
 Use 3D Plane Fit Z Position of Polygon Boundary Coordinates


**Grid Stage Step Sizes**

X Grid Step Size In Microns: .000000  
Y Grid Step Size In Microns: .000000  
X Number Of Points In Grid:   
Y Number Of Points In Grid:   
**Calculate Number Of Points In Polygon >**  
 Stagger the Grid Points

**Plot Coordinates and Boundary**

Click the **Plot Polygon Boundary** button to inspect the perimeter just drawn.

To start over and re-draw the perimeter outline again, click the **Close** button on the **Polygon Plot** window, click the **Cancel** button of the **Polygon Grid Parameters** window, and the click the **Polygon Grid** button of the **Digitize Sample Positions** window.



When satisfied with the outline of the grid, click the **Close** button of the **Polygon Plot** window.

Enter *Grid Stage Step Sizes* (in microns) for both X and Y.

**Polygon Grid Parameters**

**Polygon Grid Stage Coordinate Boundary Points**

12	27.0759	31.5605	11.1451
13	27.2444	31.5994	11.1171
14	27.0605	31.8654	11.1419
15	27.0567	31.0791	11.1420
16	27.3036	31.8187	11.1434
17	27.1575	32.0145	11.1472
18	27.4033	31.8310	11.1011
19	27.1011	31.2082	11.1176
20	27.1516	31.9513	11.1189

**Add Polygon Boundary Coordinate**

**Plot Polygon Boundary**      **Remove Last Polygon Position**

**Polygon Grid Z Interpolation Options**

Use Averaged Z Position of Polygon Boundary Coordinates  
 Use 3D Plane Fit Z Position of Polygon Boundary Coordinates

**Grid Stage Step Sizes**

X Grid Step Size In Microns: 20  
Y Grid Step Size In Microns: 20  
X Number Of Points In Grid: 37  
Y Number Of Points In Grid: 47

**Calculate Number Of Points In Polygon >**


Stagger the Grid Points

**Plot Coordinates and Boundary**


**OK**      **Cancel**      **Test**      **Load**

Click the **Calculate Number of Points in Polygon>** button to determine how many data points will be digitized. Readjust the *X and Y Grid Step Sizes* if necessary. Select a method of Z determination from the two option buttons under *Polygon Grid Z Interpolation Options*.

The **DigitizeSaveGrid** window appears with the number of points in an ideal rectangular grid.




Click the **Yes** button to calculate the total number of points.



When the appropriate gridding parameters have been set, click the **OK** button, closing the **Polygon Grid Parameters** window.

The **DigitizeSaveGrid** window re-appears, click the **Yes** button.



The program automatically digitizes each of the number of points in the polygon and returns to the **Automate!** window.

**Automate!**

Position List (multi-select) (double-click to see data)

Standards  
 Unknowns  
 Wavescans  
 All Samples

Un 1 Fid 0 SW-123 Pyroxene

Move Stage

Digitize

Plot

Fiducials

Replicates

Conditions

Sample Setups

File Setups

Multiple Setups

Automation Actions

Confirm Standard Positions  
 Confirm Unknown Positions  
 Confirm Wavescan Positions

Peak Spectrometers **Peaking**

Acquire Standard Samples  
 Acquire Unknown Samples  
 Acquire Wavescan Samples  
 Acquire Standard Samples (again)

Automation Options

Peak on Assigned Standards  
 Use "Quick" Standards  
 Use Filament Standby Afterwards  
 Use Confirm During Acquisition  
 Use Beam Deflection For Position  
 Suppress ROM Based Backlash

Confirm All Positions In Sample  
 Combine Multiple Sample Setups

Use ROM Auto Focus  
 New Sample  Every Point  
 Digitized  Interval

Standard Points To Acquire

Automate Confirm Delay (sec)

Standard X Increment (um)

Re-Standard Y Increment (um)

Re-Standard Interval (hrs)

Use Last Unknown Sample  
 Use Digitized Conditions  
 Use Digitized Sample Setups  
 Use Digitized File Setups  
 Use Digitized Multiple Setups

**Run Selected Samples**

Delete Selected Samples

Delete Selected Positions

Import from ASCII (\*.POS File)

Export Selected Samples (to \*.POS)


Row	X	Y	Z	W	Grain #	Focus
1	27.85589	30.60919	11.28489	4	1	0
2	27.83589	30.60919	11.28489	4	1	0
3	27.81589	30.60919	11.28489	4	1	0
4	27.79589	30.60919	11.28489	4	1	0
5	27.77588	30.60919	11.28489	4	1	0
6	27.73489	30.62919	11.28489	4	1	0
7	27.75489	30.62919	11.28489	4	1	0
8	27.77489	30.62919	11.28489	4	1	0
9	27.79489	30.62919	11.28489	4	1	0
10	27.81489	30.62919	11.28489	4	1	0
11	27.83489	30.62919	11.28489	4	1	0
12	27.85489	30.62919	11.28489	4	1	0
13	27.87489	30.62919	11.28489	4	1	0
14	27.89489	30.62919	11.28489	4	1	0

KeV = 15 Curr = 30 Size = 10 Mag = 400 Mode = Analog Spot Sample Setup (row) Number = 0  
 MagAnal = 2000 MagImag = 40 ImgShift = -2.3


File Setup = C:\Probe Operators\Kremser\element-sample setup.MDB (12)

Multiple Setups = NONE Replicates = 1

Click the **Plot** button in the **Automate!** dialog box to open the **Position Plot** window and view the locations of all of the digitized points in this sample. In this example, the 10 micron spacing creates too many points to be individually visible on this view.



The user may click and drag the mouse to zoom on the plot to expand the scale.



Click the **Close** button of the **Position Plot** window to return to the **Automate!** dialog box.



The user should proceed with calibration and standardization of the elements in the probe run and checking the accuracy of the standardization.

Then, to run the just digitized polygon grid sample from the **Automate!** window, highlight it in the *Position List*. Under the *Automation Actions*, click the *Acquire Unknown Samples* check box. Finally, click the **Run Selected Samples** button.

The screenshot shows the Automate! software interface. The main window is titled "Automate!". On the left, there is a "Position List (multi-select) (double-click to see data)" panel. It contains radio buttons for "Standards", "Unknowns", "Wavescans", and "All Samples". The "Unknowns" option is selected. Below these are buttons for "Select Stds", "Select All", "Go", "Auto Focus", "Update", "Delete All", and "Re-Load". The "Position List" table is currently empty. Below the table are buttons for "Delete Selected Samples", "Delete Selected Positions", "Import from ASCII (\*.POS File)", and "Export Selected Samples (to \*.POS)".

In the center, there is a "Move" and "Stage" section with buttons for "Digitize", "Plot", "Fiducials", "Replicates", "Conditions", "Sample Setups", "File Setups", and "Multiple Setups".

On the right, there are two panels: "Automation Actions" and "Automation Options".

**Automation Actions:**

- Confirm Standard Positions
- Confirm Unknown Positions
- Confirm Wavescan Positions
- Peak Spectrometers **Peaking**
- Acquire Standard Samples
- Acquire Unknown Samples
- Acquire Wavescan Samples
- Acquire Standard Samples (again)

**Automation Options:**


- Peak on Assigned Standards
- Use "Quick" Standards
- Use Filament Standby Afterwards
- Use Confirm During Acquisition
- Use Beam Deflection For Position
- Suppress ROM Based Backlash
- Confirm All Positions In Sample
- Combine Multiple Sample Setups
- Use ROM Auto Focus
  - New Sample
  - Every Point
  - Digitized
  - Interval
- Standard Points To Acquire
- Automate Confirm Delay (sec)
- Standard X Increment (um)
- Re-Standard Y Increment (um)
- Re-Standard Interval (hrs)
- Use Last Unknown Sample
- Use Digitized Conditions
- Use Digitized Sample Setups
- Use Digitized File Setups
- Use Digitized Multiple Setups

At the bottom of the interface, there is a status bar with the following information:

- KeV = 15 Curr = 30 Size = 10 Mag = 400 Mode = Analog Spot
- MagAnal = 2000 MagImag = 40 ImgShift = -2, 3
- Sample Setup (row) Number = 0
- File Setup = C:\Probe Operators\Kremser\element-sample setup.MDB (12)
- Multiple Setups = NONE
- Replicates = 1


A large yellow button labeled "Run Selected Samples" is located at the bottom right of the interface.

The **AutomateConfirmSelected** window opens and the user clicks the **Yes** button to activate the acquisition.




Note, acquisition time is now calculated.

Upon completion of the data acquisition, open the **Plot!** window.




Select (highlight) all of the unknown digitized points (Un6-35 in this example). Click the *Minimum Total* check box to skip low points (analyses in holes, etc). Select the *3-D* check button under *Graph Type*.

Click the *Send Data to ASCII File* check button. This activates the other check boxes listed here. Click the *SURFER.BAS File* check box and select an *X-Axis*, *Y-Axis* and *Z-Axis (multi-select)* parameters to plot.




Click the **Output** button. Calculations happen for all samples.

The **Open File To Save ASCII Data To** window opens. Adjust the *Save in:* location if required. Enter a *File name:* in the text field provided.




Click the **Save** button.

The **PlotDoFile** window opens, click the **OK** button.




Another **PlotDoFile** window appears.



Click the **OK** button to create these files.


The SW123.BAS script file created above contains the OLE code for generating contour and surface plots of the digitized probe data.

To run the SW123.BAS script file, double click on the **GS Scriptor32** icon in the PROBE FOR EPMA Software folder on the desktop. Select the **File | Open** menu.




The **Open Document** window appears.

Edit the *Look in:* directory to identify the location of the SW123.BAS file.



Click the **Open** button.

GS Scripiter now details the open SW123.BAS file, of which a portion is illustrated below.



```
sw-123 (script) - Scripiter [design]
File Edit View Script Debug Sheet Help
Object: (General) Proc: (declarations)
1 Directory$ = "C:\UserData\Probe Projects FY11\"
  File$ = "SW-123"
  Sample$ = "Un      6  SW-123 Pyroxene"
  MaxCol% = 6
  Dim ZLabel$(MaxCol%) As String
  XLabel$ = "X Stage Coordinates"
  YLabel$ = "Y Stage Coordinates"
  ZLabel$( 0) = "SiO2 Oxide Percents"
  ZLabel$( 1) = "Al2O3 Oxide Percents"
  ZLabel$( 2) = "FeO Oxide Percents"
  ZLabel$( 3) = "MnO Oxide Percents"
  ZLabel$( 4) = "MgO Oxide Percents"
  ZLabel$( 5) = "CaO Oxide Percents"

  ' ATTENTION! ATTENTION! ATTENTION ATTENTION! ATTENTION! Surfer 6.0 version
  ' The last UNcommented "OutputType$" variable will determine the output
  OutputType$ = "TEST"
  'OutputType$ = "SURF"

  If OutputType$ = "" Then
  Print "No output type specified"
  End
  endif


  ' The last UNcommented "GridMethodType%" variable will determine the grid method
  ' GridMethodType% = 0   ' Inverse Distance
  GridMethodType% = 1   ' Kriging
  ' GridMethodType% = 2   ' Minimum Curvature
  ' GridMethodType% = 3   ' Polynomial Regression
  ' GridMethodType% = 4   ' Radial Bias
  ' GridMethodType% = 5   ' Shepard's
  ' GridMethodType% = 6   ' Triangulation

For Help, press F1
```

The default output mode of the script file is "TEST", which will only output the plots to the screen. To produce output to the default printer, comment out the line `OutputType$ = "TEST"` by placing a single quote in front of the line and uncomment the line `OutputType$ = "SURF"` by removing the single quote in front of it (highlighted in next screen capture).




Click **Run | Start** menu to begin the automated plotting.



Basic contour and surface maps will be output to the printer. Raw data concentration (\*.GRD) files will also be created; these may be opened in SURFER for further modification and output.


An example of a basic contour map for calcium is shown below. The perimeter of the pyroxene grain is visible. Regions of higher calcium concentrations appear dark in this view.

### SW-123 PYROXENE: CALCIUM CONTOUR MAP



The next screen capture illustrates a 3-D surface map for iron in the pyroxene. Here, the image of iron concentration (vertical scale) has been rotated and tilted slightly.

### SW-123 PYROXENE: IRON SURFACE MAP



(This page intentionally left blank)

## Stage Bit Maps and Picture Snap! Feature

Unknown or standard samples loaded into the electron microprobe can present some difficulty to the user in terms of rapid and precise positioning or the location of small phases or specific areas of interest to analyze upon a large sample. On the JEOL 733 microprobe the user has several options for searching for analysis or standard locations. An optical image (reflected or transmitted light) or a video feed of the same image is available but at only one relatively high magnification, about 400 times. Additionally, one can search for the area of interest utilizing the secondary or backscattered detectors on the microprobe at variable magnifications (ranging as low as 40 times), but this can be time consuming. Still the entire sample may not be in one field of view upon observation in the chamber.

Another device employed to aid in feature location and rapid positioning is a gridding device that holds a sample mounted in a standard holder under a moveable grid system. The rough coordinates of a region on the sample may be read off and used to effectively narrow the search for the analysis position.

Now, navigation around and exact positioning is easily accomplished using the stage bit map and Picture Snap! features in PROBE FOR EPMA. The Stage Bit Map feature will be discussed first. The **Stage** button is located in many locations in PROBE FOR EPMA programs; such as in the **Acquire!** window or from any **Move Motors and Change Crystals (Move button)** window.

**Move Motors and Change Crystals!**

**Stage Target Positions**

X: 26.8989    Y: 31.0604

Z: 11.0996    W: 4

Increment: .100

Use Stage Backlash

**Remove Faraday**    **Go All**    **Go Spectros**

**Z Axis Adjust**

Increment: .0010

**Update Positions**    **Free/Clear**

**Positions**    **Stage**

**Auto Focus**

**Exchange Sample**    **Filament Standby**    **Close**


**Spectrometer Target Positions (Load Element Setups From Acquire Elements/Cations Button)**

1	2	3			
PET	TAP	LIF			
228.140	89.5139	191.106			
Si ka	2	Ti ka	4	5	6

Use Spectrometer Backlash

**Jog Spectrometers**    **Park Spectrometers**

Clicking the **Stage** button opens the **Stage Map** window. Two different maps are displayed below.




To select another *Holder Selection* image, simply select the file from the drop-down list box. Image files (windows metafiles (\*.WMF)) and coordinate limits are specified in the Standards section of the PROBEWIN.INI file. The entire map maybe reduced or enlarged retaining scale using the **x0.8** or **x1.2** buttons or to re-size the Stage Bit Map window simply drag any corner of the window to the desired size and shape. The **minus** and **plus** button (upper left) minimizes the stage bitmap selection and cursor position display.

The current position is indicated as a small red-purple circle on the map. To move from one location to another, simply double-click on the spot you wish the stage to travel to. The current position (X and Y stage coordinates) is displayed above the **Remove Faraday/ Insert Faraday** button. Digitized positions of various samples can also be viewed by selecting the appropriate radio button.


To create stage drawing maps of your standard holders, for instance, use a vector based drawing program (Micrografx Designer or the shareware program Metafile Companion<sup>1</sup>) and the exact dimensions of your holders to build dimensionally correct drawings. These can be exported as windows metafiles and directly loaded into the graphical stage move feature in PROBE FOR EPMA.

<sup>1</sup>Mention of specific third party software products does not imply their endorsement.

Each stage map must be calibrated in coordinate space for accurate movement to features on the map. Typically two diagonally located points near the edge of the map are chosen for calibration. Initiate the calibration routine by clicking the @ button (upper right) in the **Stage Map** window. The **MoveStageMapCalibrate** window appears.



Click the **Yes** button to open the **Calibrate Stage Bit Map Min and Max** window for calibration.



Click the top **Select** button, opening the **Stage Select** window.



Click on the unique position on the stage map to identify the stage coordinates.

These values appear next to the **Select** button chosen.

**Calibrate Stage Bit Map Min and Max**

**Calibration Points**

Select  Update   1st

Select  Update   2nd

Click each "Select" button and click the StageMap image on a position for which the actual stage coordinates are known (coordinates displayed are calculated from corner positions)

Enter actual stage positions above (using the "Update" buttons or entering values manually) based on the optical stage coordinates.

Edit in the values in the [Standards] section of your PROBEWIN.INI file for the calculated corner values. Depending on the polarity of your stage axes, the min and max may need to be swapped.

**Calculate Corners**   **Min**

**Max**

**Re-Load INI Values**

**Close**

Next, activate the microprobe imaging and zoom up in magnification to locate the exact spot you just selected. Either click the **Update** button or manually enter the stage coordinate information for the 1<sup>st</sup> calibration point.

**Calibrate Stage Bit Map Min and Max**

**Calibration Points**

Select  Update   1st

Select  Update   2nd

Click each "Select" button and click the StageMap image on a position for which the actual stage coordinates are known (coordinates displayed are calculated from corner positions)

Enter actual stage positions above (using the "Update" buttons or entering values manually) based on the optical stage coordinates.

Edit in the values in the [Standards] section of your PROBEWIN.INI file for the calculated corner values. Depending on the polarity of your stage axes, the min and max may need to be swapped.

**Calculate Corners**   **Min**

**Max**

**Re-Load INI Values**

**Close**



Click the lower **Select** button and repeat the process. Click on the second position on the image. Activate the imaging and find this exact point and update the position. The **Calibrate Stage Bit Map Min and Max** window will appear as below.

**Calibrate Stage Bit Map Min and Max**

**Calibration Points**

	X	Y	
Select	3.46154, 18.2526		
Update	3.6000	19.2500	1st
Select	28.4615, 40.6140		
Update	28.5600	41.2499	2nd

Click each "Select" button and click the StageMap image on a position for which the actual stage coordinates are known (coordinates displayed are calculated from corner positions)

Enter actual stage positions above (using the "Update" buttons or entering values manually) based on the optical stage coordinates.

Edit in the values in the [Standards] section of your PROBEWIN.INI file for the calculated corner values. Depending on the polarity of your stage axes, the min and max may need to be swapped.

**Calculate Corners**

	X	Y	
Min			
Max			

**Re-Load INI Values**

**Close**

Click the **Calculate Corners** button to obtain the correct corner values to calibrate your Stage Map. These values min and max values are entered into the Standards section of the PROBEWIN.INI file.

**Calibrate Stage Bit Map Min and Max**

**Calibration Points**

	X	Y	
Select	3.46154, 18.2526		
Update	3.6000	19.2500	1st
Select	28.4615, 40.6140		
Update	28.5600	41.2499	2nd

Click each "Select" button and click the StageMap image on a position for which the actual stage coordinates are known (coordinates displayed are calculated from corner positions)

Enter actual stage positions above (using the "Update" buttons or entering values manually) based on the optical stage coordinates.

Edit in the values in the [Standards] section of your PROBEWIN.INI file for the calculated corner values. Depending on the polarity of your stage axes, the min and max may need to be swapped.

**Calculate Corners**

	X	Y	
Min	-3.8496	1.29248	
Max	36.0864	61.3063	


**Re-Load INI Values**

**Close**


Now this image is calibrated, so the user can easily “drive around” on the image.

The new **Picture Snap!** feature allows the user to incorporate images of your unknown thin section or polished mounts into PROBE FOR EPMA to aid in navigation and the digitizing of analysis locations. Images (BMP, JPEG, GRD) taken with a flatbed scanner or other camera system can be entered into **Picture Snap!**, then calibrated and used for analysis.


**Picture Snap!** dialog can be accessed from the **Digitize Sample Positions** window.




**Picture Snap!** can also be accessed from the STAGE program by accessing the **Window | Picture Snap!** menu.




The main **Picture Snap!** window appears.



Select the **F**ile menu and open the appropriate image file.




The **Open File To Input Bitmap Image From** window opens.




Select the appropriate directory and file to open and click the **Open** button.

The image is displayed in the **Picture Snap!** window.





Select the **Window | Calibrate** menu.





The **Picture Snap Calibration** window appears.

**Picture Snap Calibration**

**Point #1 Calibration**

X Pixel Coordinate

Y Pixel Coordinate

Pick Pixel Coordinate on Picture

X Stage Coordinate

Y Stage Coordinate

Read Current Stage Coordinate    Move To

**Point #2 Calibration**

X Pixel Coordinate

Y Pixel Coordinate

Pick Pixel Coordinate on Picture

X Stage Coordinate

Y Stage Coordinate

Read Current Stage Coordinate    Move To

Close

Two Points  
 Three Points

To calibrate the picture select two diagonal coordinates (for rectangular samples) or three points (for round samples subject to rotation) on the picture for which you can accurately locate the stage positions. Then click the Calibrate Picture button.

Calibrate Picture

Image Is NOT Calibrated

Display Calibration Points

**Light Mode**

Refl	Tran
On	Off

X and Y Pixel Coordinates are actually given in "Twip" units! (1440 twips per logical inch)

Image calibration is accomplished using a two point method for rectangular mounts or a three point calibration when importing round images subject to rotation. Click the Point #1 Calibration **Pick Pixel Coordinate on Picture** button. The **Picture Select Point** window appears, select the first unique point on the image.

**Picture Select Point**

Select a position on the image to calibrate the picture

Cancel

The X,Y Pixel Coordinates are entered.

**Picture Snap Calibration**

**Point #1 Calibration**

X Pixel Coordinate: 11040  
Y Pixel Coordinate: 612  
Pick Pixel Coordinate on Picture

X Stage Coordinate: 28.4616  
Y Stage Coordinate: 40.61312  
Read Current Stage Coordinate Move To

**Point #2 Calibration**

X Pixel Coordinate:   
Y Pixel Coordinate:   
Pick Pixel Coordinate on Picture

X Stage Coordinate: 28.4616  
Y Stage Coordinate: 40.61312  
Read Current Stage Coordinate Move To

**Close**

**Two Points**  
 **Three Points**

To calibrate the picture select two diagonal coordinates (for rectangular samples) or three points (for round samples subject to rotation) on the picture for which you can accurately locate the stage positions. Then click the Calibrate Picture button.

**Calibrate Picture**

Image Is NOT Calibrated

**Display Calibration Points**

**Light Mode**

Refl	Tran
On	Off

X and Y Pixel Coordinates are actually given in "Twip" units! (1440 twips per logical inch)

The values shown in the X, Y Stage Coordinates text boxes are the current stage location. Drive the stage to the same unique location and click the **Read Current Stage Coordinate** button.

The stage location for the first calibration point is entered.

**Picture Snap Calibration**

**Point #1 Calibration**

X Pixel Coordinate

Y Pixel Coordinate

X Stage Coordinate

Y Stage Coordinate

**Two Points**

**Three Points**

To calibrate the picture select two diagonal coordinates (for rectangular samples) or three points (for round samples subject to rotation) on the picture for which you can accurately locate the stage positions. Then click the Calibrate Picture button.

Image Is NOT Calibrated

**Point #2 Calibration**

X Pixel Coordinate

Y Pixel Coordinate

X Stage Coordinate

Y Stage Coordinate

**Light Mode**

Refl	Tran
On	Off

X and Y Pixel Coordinates are actually given in "Twip" units! (1440 twips per logical inch)

Repeat these steps for the second calibration point, resulting in the following window.

**Picture Snap Calibration**

**Point #1 Calibration**

X Pixel Coordinate: 11040  
Y Pixel Coordinate: 612  
**Pick Pixel Coordinate on Picture**

X Stage Coordinate: 9.700018  
Y Stage Coordinate: 32.19991  
**Read Current Stage Coordinate** **Move To**

**Point #2 Calibration**

X Pixel Coordinate: 2460  
Y Pixel Coordinate: 9852  
**Pick Pixel Coordinate on Picture**

X Stage Coordinate: 24.60079  
Y Stage Coordinate: 41.70004  
**Read Current Stage Coordinate** **Move To**

**Close**

**Two Points**  
 **Three Points**

To calibrate the picture select two diagonal coordinates (for rectangular samples) or three points (for round samples subject to rotation) on the picture for which you can accurately locate the stage positions. Then click the Calibrate Picture button.

**Calibrate Picture**

Image Is NOT Calibrated


**Display Calibration Points**

**Light Mode**

Refl	Tran
On	Off


X and Y Pixel Coordinates are actually given in "Twip" units! (1440 twips per logical inch)

Finally, click the **Calibrate Picture** button opening the **PictureSnapSaveCalibration** window.



Click the **OK** button to save the picture calibration.


The **Picture Snap Calibration** window now indicates that the *Image Is Calibrated*.



Close the **Picture Snap Calibration** window.

The operator can digitize analysis locations for later unattended work. Click the **Digitize** button.

This opens the familiar **Digitize Sample Positions** window.



Create a new unknown. Double click on the spot for the first analysis point on the just calibrated image to drive the stage to those coordinates. Click the **Random Point(s)** button to digitize that sample location. Additional points maybe saved. All analysis locations can be viewed from the **Display | Unknown Position Samples** menu. All analysis locations can be acquired via the **Automate!** Window.

# Modal Analysis

Modal analysis is a statement of the composition of a sample expressed in terms of the relative amounts of phases or minerals present. These volumetric proportions can be estimated from quantitative measurements made on the specimen by point counting analysis. This quantitative modal analysis on unknown compositions is based on a defined set of modal phases, selected from a standard database. Any database of standard compositions may be used to define the phases.

There are three basic steps involved in the modal analysis routine. This procedure involves initially the acquisition of a large set of compositional data acquired using either multiple traverses or large area gridding. It is assumed that this data set is statistically representative of the sample. In the example illustrated below, a large representative area of a fine-grained sandstone thin section was gridded and some 324 quantitative analysis points were collected.

The second step involves the creation of an input file to load into STANDARD for the actual modal analysis calculation. The simplest method of generating this input file is to use the **Plot!** window in PROBE FOR EPMA to output a \*.DAT file of the elemental or oxide weight percent compositions to disk.


After data collection has been completed, open the **Plot!** window.

Highlight the compositional dataset in the *Sample List*.

Select the first oxide for the *X-Axis* and the remainder in the *Y-Axis (multi-select)* range.

Activate the *Send Data to ASCII File (X, Y, (Z)...) check button and the ASCII File Column Labels check box. The labels are required so that the modal analysis routine can identify the elements in the input file.*


The input file can come from any source as long as the element or oxide symbols are in the first line, enclosed within double quotes, and the data is in weight percent. The weight percent data can be in any format. Do not include a totals column.



Click the **Output** button.




The **Open File To Save ASCII Data To** window appears. Locate the appropriate directory under *Save in:* and type in a *File name:* in the text field provided.



Click the **Save** button.

The **PlotDoFile** window appears, indicating that the data was saved.




Click the **OK** button.

The data saved to the \*.DAT file may be viewed using an editor such as Notepad. Here a portion of the CITYHALL.DAT file is displayed.


%CO2 Oxide Percents*	%Na2O Oxide Percents*	%FeO Oxide Percents*	%CaO Oxide Percents*	%Al2O3 Oxide Percents*	%SiO2 Oxide Percents*	%MgO Oxide Percents*	%TiO2 Oxide Percents*
3.212856	084278	057376	434429	3.858789	11.172727	001379	018428
029580	000000	005280	000000	080792	98.047482	000000	073651
010581	005889	050280	005228	000000	97.680529	018088	000000
015457	000000	028683	023569	000000	100.941032	000000	002045
188080	000000	218054	291782	3.940241	81.151772	024884	000000
005303	000000	025585	003580	000000	102.645638	000000	000000
000000	000000	084513	040705	1.15520	83.195908	000000	000000
007942	008243	023802	000000	000000	104.154282	000000	083880
088808	008581	221087	178842	581872	22.384785	023032	000018
083358	004089	075925	148523	500850	81.950417	044080	009942
009298	001881	002843	017831	000000	104.538147	000000	000000
086057	000000	078280	043218	1.223209	100.848114	001021	058022
012381	006702	000000	000000	041418	105.016945	000000	007092
1.388605	429426	10.343257	1.418728	15.828440	38.047485	1.831048	208419
15.718748	323534	015711	1.830254	16.165831	86.888898	024850	588850
2.184575	040035	2.433222	1.179088	18.731405	84.802748	547108	128404
088051	040871	044086	151371	154820	78.895401	029395	042584
825883	005728	375800	881520	8.238878	28.404783	142400	052258
005781	000086	000000	000394	000000	107.231512	057088	046888
883888	041141	318287	197038	2.921450	100.892586	184181	089282
824219	000000	623511	302535	378020	25.882721	000000	148085
202848	017280	324789	107448	1.514305	98.554527	088814	038251
154851	013822	351124	520921	1.338327	26.883922	043340	045254
000000	000000	000000	002255	000000	105.854887	020756	007077
027884	000000	041379	031258	000000	108.888490	001384	028788
064887	2.848057	1.18134	321882	4.812264	19.388606	000000	000000
008830	000000	080733	018485	002112	106.278352	048057	000000
3.74252	018848	8.948800	854708	2.285708	3.738187	225787	000000
004888	000000	018887	000804	000000	108.388805	015886	058281
000000	000000	000000	000000	000000	105.473846	000000	000000
024729	001342	213426	289289	648775	98.803633	014366	038330
328807	008876	282103	102549	778127	100.378180	058084	002847
038830	000000	248512	510183	287839	87.087351	048615	000000
711713	075025	388320	10.025884	27.285225	35.033842	288004	000000
5.031882	018844	884308	845424	1.058110	34.880588	003153	018432
1.88255	718878	19.230335	1.933186	5.358851	11.788387	288482	038944
062418	013558	438782	888028	770783	53.322881	017525	000000
018270	008575	011488	858180	522753	50.575278	000000	015835
233282	028784	444485	1.041038	2.588887	52.403080	034018	002833
003549	000000	017840	019700	000000	102.357887	000000	000000
15.388772	458128	088816	127881	17.728888	83.582888	050283	048788
011483	000000	080118	003187	073101	213154	002630	008184
148872	238388	584182	208818	4.384525	83.453882	081257	000000
033851	010748	218487	081888	000000	34.851734	000000	000000
783514	037701	328822	520783	2.183847	83.350385	070883	014188
8.148172	000000	7.837807	340207	18.815746	51.830448	3.062872	488256
035288	017275	085572	234828	311152	8.418805	012384	000000
000000	000000	000000	840770	000000	101.852583	000000	023814
000000	000180	028048	000578	000000	100.883883	000000	028342
028278	008007	023818	077138	098718	86.557487	022888	015818
018123	018542	003528	000000	000000	100.848188	000000	021288
188222	000000	048482	133885	488882	88.002083	008835	000000
288548	047021	083484	008788	888440	98.107885	000000	188438
054181	013174	088807	142538	204400	76.857407	000000	021238
000885	000000	012333	000000	000000	98.854172	000000	028358

The third and final step involves the setup and running of the modal analysis calculation. The modal analysis routine is located in STANDARD. Open the program **Standard** from the PROBE for EPMA Software folder on the desktop.



Select (highlight) a standard database that will be used to define the modal phases. Click the **Open** button to load this database.

Select **Options** | **Modal Analysis** from the menu.



The **Modal Analysis** window opens.

**Modal Analysis**

**Group Definitions**

Groups	Phases	Standards
<b>New</b> Delete	<b>New</b> Delete	<b>Add</b> Remove

**Close**

**Start**

**Group Options**

Minimum Total for Input

Do End-Member Calculations  
 Normalize Concentrations For Fit  
 Weight Concentrations For Fit

**Update Group**

**Phase Options**

Minimum Vector

None     Olivine     Pyroxene  
 Feldspar     Garnet

**Update Phase**

**Data Files**

Input Data File


 **Browse**

Output Data File

 **Browse**

Start by defining an overall Group, click the **New** button under *Groups*.


The **ModalGroupNew** window opens. Enter a descriptive name for the group of phases.



Click the **OK** button. Default *Group* and *Phase Options* are loaded; these will be discussed and modified shortly.


Click the **New** button under *Phases*.

The **ModalPhaseNew** window opens. Enter the first modal phase. In this example, the sandstone is composed of mostly quartz with two minor feldspars; an alkali (sodium-potassium) phase and a plagioclase phase along with iron oxides and other trace accessory minerals. The first modal phase is entered into the text field.



Click the **OK** button.

Select the **Add** button under *Standards*, opening the **Add Standards to Run** window.



Choose standards to define this modal phase. These are the phase compositions that the program will use to match against the unknown point analyses. Try to avoid over-determining the phase. For example, when defining a sodium-potassium feldspar, select the two end-members (albite and microcline).

The **Modal Analysis** window would now appear as below.

The screenshot shows the 'Modal Analysis' window with the following sections:

- Group Definitions:** Contains three columns: 'Groups' (with 'Sandstone'), 'Phases' (with 'Quartz'), and 'Standards' (with '319 SiO2 Quartz Tay'). Below each column are 'New' and 'Delete' buttons. A 'Start' button is located to the right of this section.
- Group Options:** Includes a 'Minimum Total for Input' field set to '95', and checkboxes for 'Do End-Member Calculations' (unchecked), 'Normalize Concentrations For Fit' (checked), and 'Weight Concentrations For Fit' (checked). An 'Update Group' button is at the bottom.
- Phase Options:** Includes a 'Minimum Vector' field set to '4', and radio buttons for 'None' (selected), 'Olivine', 'Pyroxene', 'Feldspar', and 'Garnet'. An 'Update Phase' button is at the bottom.
- Data Files:** Includes 'Input Data File' and 'Output Data File' fields, both containing 'C:\UserData\modal.out', with 'Browse' buttons to the right.

A green 'Close' button is located in the top right corner of the window.



Continue and enter all phases, defining the phase compositions (standards) to match. The Alkali Feldspar entry is illustrated below.

**Modal Analysis**

**Group Definitions**

Groups	Phases	Standards
Sandstone	Quartz Plagioclase Alkali Feldspar	81 Albite 116 Microcline

**Buttons:** New, Delete (under Groups); New, Delete (under Phases); Add, Remove (under Standards); Close, Start

**Group Options**

Minimum Total for Input: 95

Do End-Member Calculations  
 Normalize Concentrations For Fit  
 Weight Concentrations For Fit

**Update Group**

**Phase Options**

Minimum Vector: 4

None     Olivine     Pyroxene  
 Feldspar     Garnet

**Update Phase**

**Data Files**

Input Data File: C:\UserData\modal.out **Browse**

Output Data File: C:\UserData\modal.out **Browse**

Once all of the phases have been identified and standards defined for matching, adjust the *Group* and *Phase Options*.

**Modal Analysis**

**Group Definitions**

Groups	Phases	Standards
Sandstone	Quartz Plagioclase Alkali Feldspar Opagues	41 Ilmenite 115 Magnetite

**Group Options**

Minimum Total for Input:

Do End-Member Calculations  
 Normalize Concentrations For Fit  
 Weight Concentrations For Fit

**Phase Options**

Minimum Vector:

None
  Olivine
  Pyroxene  
 Feldspar
  Garnet

**Data Files**

Input Data File:

Output Data File:

The *Minimum Total for Input* is the rejection sum for the unknown compositions, sums below this value will not be used in the modal analysis. Typically 90-95% are good cutoffs.

Select the *Do End-Member Calculations* option and check the appropriate mineral name under *Phase Options* to perform end-member calculations as listed.

The *Normalize Concentrations For Fit* option is used to specify whether the standard and unknown concentrations (above the just defined minimum input total) should be normalized to 100% before the vector fit is calculated.

The *Weight Concentrations For Fit* option is used to specify if the element concentrations for the standards should be weighted, based on the composition of the element in that phase. Select this option if the major elements in a phase should have greater influence in determining the vector fit. Leave unselected, if all concentrations, regardless of their abundance should have equal weight in the vector fit.

The *Minimum Vector* number (default is 4.0) is basically the tolerance for the match to a defined phase. If a closer match is desired for one or more phases in the group, decrease the vector value for that phase. See the User's Guide and Reference documentation for specific details on the calculation of this vector.

Finally, under *Data Files*, select the appropriate *Input* and *Output Data File* locations.

**Modal Analysis**

**Group Definitions**

Groups	Phases	Standards
Sandstone	Quartz Plagioclase Alkali Feldspar Opagues	41 Ilmenite 115 Magnetite

New Delete      New Delete      Add Remove      Start

**Group Options**

Minimum Total for Input: 90

Do End-Member Calculations  
 Normalize Concentrations For Fit  
 Weight Concentrations For Fit

Update Group

**Phase Options**

Minimum Vector: 4

None       Olivine       Pyroxene  
 Feldspar       Garnet

Update Phase

**Data Files**


Input Data File: C:\UserData\Probe Projects FY11\cityhall.dat      Browse

Output Data File: C:\UserData\Probe Projects FY11\modal.out      Browse

Close

Click the **Start** button to initiate the modal analysis calculation on each data point.

After the calculation finishes the **ModalStartModal** window appears, stating that the output data has been saved to the specified \*.OUT file.



Click this **OK** button.

The modal analysis data may now be viewed in the log window in STANDARD or simply open the newly created \*.OUT file. The output file contains the vector fit, matched phase, end-member calculation (if checked), totals column and composition of each line in the input file. Lines that do not meet the minimum total requirement are excluded from the output, if those lines are desired either cut and paste the entire output from the main STANDARD log window or capture the entire output by EARLIER selecting the **Output | Save To Disk Log** menu.

The results of the modal analysis are also tabulated and summarized. The end summary lists the total number of analyses, the minimum total for a valid composition, number of valid points that sum above the minimum sum, the number of matched points and the percentage of points that were matched.

For each phase, the summary output then lists the phase name, the number of matches for that phase, the percent of matched points, valid points and total matched points for the matches in that phase. This is followed by the average end-member (if selected), weight percent sum and composition for that phase and the standard deviation for each element.

The last page of the just run output file is displayed below.

Line	Vector	Phase	Sum	K2O	Na2O	FeO	CaO	Al2O3	SiO2	MgO	TiO2
.....268	.05	Quartz	100.30	.03	.02	.10	.10	.08	99.93	.00	.04
269	.04	Quartz	99.50	.02	.01	.12	.00	.12	99.17	.00	.06
270	.00	Quartz	99.66	.00	.00	.02	.00	.00	99.60	.00	.04
271	.06	Quartz	99.97	.02	.01	.20	.00	.12	99.62	.00	.00
272	-----	-----	.21	.04	.00	.00	.02	.03	.12	.01	.00
273	.00	Quartz	99.10	.01	.00	.02	.01	.00	99.04	.00	.02
274	.02	Quartz	100.07	.00	.00	.01	.00	.19	99.87	.00	.00
275	.02	Alkali	99.50	.23	11.50	.32	.03	18.85	68.57	.00	.00
276	-----	-----	36.50	5.25	.13	1.22	1.35	6.90	21.61	.04	.00
277	.00	Quartz	99.68	.00	.00	.01	.00	.00	99.68	.00	.00
278	.04	Quartz	98.30	.07	.00	.00	.03	.15	97.96	.03	.06
279	-----	-----	.13	.00	.00	.05	.02	.01	.06	.00	.00
280	.02	Quartz	99.17	.00	.00	.13	.03	.00	99.00	.01	.00
281	.06	Alkali	98.43	15.22	.24	.07	.01	18.23	64.65	.01	.00
282	.00	Quartz	99.79	.01	.00	.03	.00	.00	99.75	.00	.00
283	.01	Plagioc	98.09	.13	.02	.03	18.91	35.54	43.35	.07	.03
284	.00	Quartz	99.75	.00	.00	.04	.00	.00	99.68	.02	.01
285	.00	Quartz	99.59	.01	.00	.04	.00	.00	99.53	.02	.00
286	.67	Quartz	92.61	.07	.01	.29	.54	.16	91.46	.02	.06

287	.00	Quartz	99.43	.00	.00	.00	.03	.03	99.32	.00	.05
288	.02	Opagues	92.35	.03	.04	88.00	.01	.12	.42	.38	3.36
289	-----	-----	31.13	.23	.01	.15	1.66	1.45	27.55	.08	.00
290	.02	Quartz	99.65	.00	.00	.13	.00	.00	99.52	.00	.00
291	.13	Plagioc	99.23	.87	.08	.03	18.34	36.21	43.59	.09	.01
292	-----	-----	13.55	.04	.27	.07	.30	.19	12.55	.11	.03
293	.03	Plagioc	99.23	.03	.30	.12	18.88	36.27	43.54	.02	.07
294	-----	-----	9.42	.45	.00	.47	.56	.07	7.72	.00	.15
295	-----	-----	35.75	.30	.01	.61	.72	11.33	22.68	.09	.00
296	-----	-----	.86	.00	.00	.01	.01	.00	.85	.00	.00
297	.00	Quartz	99.88	.01	.00	.04	.01	.03	99.78	.00	.01
298	.21	Quartz	93.18	.13	.02	.20	.07	.20	92.48	.01	.07
299	.10	Quartz	98.14	.05	.03	.14	.12	.12	97.63	.03	.02
300	.00	Quartz	100.08	.00	.00	.04	.00	.00	100.05	.00	.00
301	.01	Quartz	99.49	.01	.00	.04	.02	.00	99.34	.02	.06
302	.02	Quartz	100.11	.00	.00	.08	.04	.00	99.91	.02	.07
303	.00	Quartz	100.32	.00	.00	.02	.02	.00	100.28	.00	.01
304	.00	Quartz	99.95	.00	.00	.00	.00	.00	99.95	.00	.00
305	.00	Quartz	99.97	.00	.00	.00	.00	.00	99.92	.00	.05
306	.02	Quartz	100.39	.01	.00	.09	.01	.00	100.21	.00	.06
307	-----	-----	7.62	.79	.00	.12	.33	1.09	5.27	.01	.00
308	-----	-----	15.00	.14	.01	.22	.22	.72	13.64	.04	.01
309	-----	-----	2.08	.06	.01	.14	.23	.36	1.26	.02	.00
310	.01	Plagioc	98.90	.14	.15	.01	19.20	35.63	43.65	.06	.06
311	.00	Quartz	101.36	.02	.00	.00	.02	.00	101.32	.00	.01
312	.02	Plagioc	99.02	.27	.01	.00	19.01	35.78	43.84	.11	.00
313	-----	-----	19.98	.32	.02	1.38	1.25	4.13	12.56	.30	.01
314	-----	-----	44.17	.04	.00	.15	.26	17.47	26.15	.01	.08
315	-----	-----	6.51	.00	.01	.05	.00	.10	6.34	.00	.00
316	.01	Quartz	100.85	.01	.00	.07	.00	.00	100.76	.00	.00
317	-----	-----	.40	.00	.00	.02	.02	.02	.32	.01	.01
318	.01	Quartz	100.95	.03	.00	.05	.02	.00	100.81	.04	.01
319	.00	Quartz	101.53	.00	.00	.03	.01	.00	101.46	.00	.03
320	.56	Quartz	98.51	.06	.00	.01	.33	.71	97.31	.07	.01
321	.01	Quartz	101.04	.02	.00	.07	.00	.00	100.92	.03	.00
322	.01	Quartz	101.03	.02	.00	.05	.02	.01	100.88	.00	.06
323	.00	Quartz	100.46	.00	.00	.00	.00	.00	100.45	.00	.01
324	-----	-----	.13	.00	.01	.02	.00	.00	.08	.00	.02

Results of Modal Analysis

InputFile : C:\UserData\Probe Projects FY11\cityhall.dat  
OutputFile : C:\UserData\Probe Projects FY11\modal.out  
Date and Time: 1/24/2011 7:46:51 PM

Group Name : Sandstone  
Total Number of Points in File : 324  
Valid Number of Points in File : 240  
Match Number of Points in File : 237

Minimum Total for Valid Points : 90.00  
Percentage of Valid Points : 74.1  
Percentage of Match Points : 73.1

Phase	#Match	%Total	%Valid	%Match	AvgVec					
Quartz	190	58.6	79.2	80.2	.09					
Sum		K2O	Na2O	FeO	CaO	Al2O3	SiO2	MgO	TiO2	
Average:	99.27	.03	.01	.06	.04	.10	98.99	.01	.03	
Std Dev:	2.00	.08	.03	.08	.08	.21	2.12	.03	.03	
Minimum:	90.60	.00	.00	.00	.00	.00	89.53	.00	.00	
Maximum:	102.77	.89	.34	.50	.54	1.52	101.46	.32	.20	

Phase	#Match	%Total	%Valid	%Match	AvgVec					
Plagiocl	20	6.2	8.3	8.4	.05					
	Sum	K2O	Na2O	FeO	CaO	Al2O3	SiO2	MgO	TiO2	
Average:	99.15	.31	.05	.08	18.73	36.01	43.90	.05	.02	
Std Dev:	.78	.34	.07	.10	.47	.28	.40	.03	.02	
Minimum:	97.58	.02	.00	.00	17.95	35.54	43.32	.00	.00	
Maximum:	100.25	1.31	.30	.39	19.42	36.48	44.88	.11	.07	

Phase	#Match	%Total	%Valid	%Match	AvgVec					
Alkali F	21	6.5	8.8	8.9	.06					
	Sum	K2O	Na2O	FeO	CaO	Al2O3	SiO2	MgO	TiO2	
Average:	99.11	13.19	1.83	.07	.22	18.40	65.33	.02	.03	
Std Dev:	.61	5.39	4.01	.07	.19	.28	1.50	.03	.04	
Minimum:	97.60	.11	.02	.00	.01	18.02	63.91	.00	.00	
Maximum:	100.05	15.86	11.65	.32	.71	18.85	69.24	.08	.13	

Phase	#Match	%Total	%Valid	%Match	AvgVec					
Opagues	6	1.9	2.5	2.5	.16					
	Sum	K2O	Na2O	FeO	CaO	Al2O3	SiO2	MgO	TiO2	
Average:	91.69	.09	.05	88.86	.18	.39	.94	.43	.75	
Std Dev:	1.10	.07	.04	.64	.22	.46	.43	.33	1.31	
Minimum:	90.33	.01	.02	88.00	.00	.12	.42	.22	.03	
Maximum:	93.13	.18	.11	89.56	.59	1.29	1.38	1.09	3.36	

Click the **Close** button on the **Modal Analysis** window.

Finish by exiting STANDARD.

(This page intentionally left blank)



## Deadtime Calculations

This section describes how to calibrate the deadtime constants for each spectrometer and where to enter them so that PROBE FOR EPMA will utilize these factors.

Deadtime ( $\tau$ ) is defined as the time interval (after arrival of a pulse) when the counting system does not respond to additional incoming pulses (Reed, 1993). The equation normally used to correct for deadtime losses is given as:

$$n = \frac{n'}{(1 - \tau n')} \quad (1)$$

Where:  $n$  is the deadtime corrected count rate in counts per second  
 $n'$  is the measured count rate in counts per second  
 $\tau$  is the deadtime constant in seconds

The time interval when the counting system is dead to additional pulses is defined as  $\tau n'$ . The live time then, is  $(1 - \tau n')$ . The true count rate ( $n$ ) is proportional to the beam current ( $i$ ) by a constant factor, designated  $k$ . Thus, equation (1) may be rewritten as:

$$\frac{n'}{i} = k(1 - \tau n') \quad (2)$$

A plot of  $n'/i$  (cps/nA) versus  $n'$  (cps) will yield a straight line with slope of  $(-k\tau)$ . The intercept on the  $n'/i$  axis will be the constant,  $k$ , and thus the deadtime factor ( $\tau$ ) may be determined.

A second deadtime correction option is also available in PROBE FOR EPMA. This is a high precision expression for use with very high count rates (Willis, 1993). This expression differs from the normal equation **only** when very high count rates (>50K cps) are achieved. The precision deadtime expression is:

$$n = \frac{n'}{1 - \left( \tau n' + n'^2 \left( \frac{\tau^2}{2} \right) \right)} \quad (3)$$

The deadtime correction option and type is selected from the **Analysis Calculation Options** window. Click **Analytical | Analysis Options** menu from the main PROBE FOR EPMA log window. Click the **OK** button to confirm the selections.

STARTWIN can be used to obtain the x-ray intensities required for the deadtime calculation. The procedure involves collecting precise beam current and count rate data over a wide range of beam currents. This data set can then be loaded into the supplied Excel template to automatically calculate the deadtime factor for your spectrometers. Paul Carpenter has put together an excellent but slightly more elaborate Excel template, contact Probe Software, Inc. for further details on obtaining his spreadsheet and related documentation.

To calibrate the deadtime factors for your WDS system use high purity, homogeneous metal standards. Depending on the microprobe configuration one standard may be employed to collect data on all spectrometers. Here, a silicon metal standard will be used.

Open the **Count Times** window and disable both the *Use Beam Drift Correction* and the *Normalize To Counts Per Second* options to allow raw intensity data to be collected. Set an *On Peak Count Time* that will give a precise measurement of intensities.

Peak each spectrometer to the x-ray line that will be used (Si K $\alpha$  on the PET and TAP crystals in the JEOL 733). Upon completion of the peak center routine, park the spectrometers on the new peak positions.

Prior to collecting data, run PHA scans on each spectrometer for Si  $K\alpha$ , check the pulse height distribution at low and very high beam currents (ideally duplicating the range of beam currents for the deadtime measurements). At very high count rates (large beam currents), significant pulse pileup and gain shifts do occur. Fully open your pulse height windows, optimize your gain settings to see all the signal over the range of beam currents employed.

Data collection and analysis is straightforward. Select **Output | Open Link To Excel** menu from the main STARTWIN log window. Collect three replicate intensity measurements and beam current data. Each time count rate data is acquired, it will automatically be sent to an Excel spreadsheet along with column labels. Measure the replicate count intensities at ten different beam currents; ranging from a few nanoamps to several hundred nanoamps.

Create a count time column, prior to the beam current column in the Excel raw data spreadsheet and enter the relevant count times (in this example 10 seconds was used). The resulting spreadsheet may look similar to the one printed below except you may have data from more than three spectrometers.

Time	Beam	1	2	3
10	3.259	22083	79765	19470
10	3.263	22069	79995	19590
10	3.247	21755	80091	19679
10	7.08	42503	154130	37952
10	7.06	42642	154665	38087
10	7.072	42168	154182	38282
10	14.547	83163	293002	74572
10	14.539	83315	292602	74281
10	14.543	82998	293326	74636
10	29.911	162904	547744	147492
10	29.917	164053	549493	147209
10	29.935	163684	548386	147078
10	50.539	266841	838976	240665
10	50.562	266672	837625	240860
10	50.58	267751	837948	240463
10	80.844	409290	1169741	370608
10	80.856	410142	1169175	370821
10	80.933	409098	1168965	370368
10	103.225	507351	1352944	460976
10	103.229	507971	1354039	460353
10	103.235	508408	1352991	460165
10	153.811	709360	1640458	647802
10	153.828	711086	1640015	647158
10	153.871	710654	1641034	648664
10	199.604	873187	1790253	800511
10	199.545	871582	1788206	799268
10	199.402	873093	1788780	799117
10	248.192	1027320	1878180	945455
10	248.884	1027517	1878519	945061
10	248.947	1026681	1878716	945783

Open the DEADTIME\_CALC.XLS file from the floppy disk supplied. Copy and paste count times, beam current information and counts for the first spectrometer into the raw data template starting in cell **A26**.

By placing data into this template, the program will automatically calculate the following items: the average of three replicate time counts, the average of three replicate beam current measurements, the %RSD on the average beam current, the average of three replicate raw intensity measurements and the %RSD on the average raw intensity measurement.

Next, the counts per second (x-axis) and the counts per second per nanoamp (y-axis) are determined. A least squares method is then used to calculate a straight line that best fits your data. The slope and Y-intercept are reported for a straight line fit to all 10 data pairs and also for the last 6. The latter being a more precise determination of deadtime.

Below is the calculation portion of the Excel template.

	A	B	C	D	E	F	G	H	I	J	K
1	<b>Deadtime Calculations</b>										
2	<b>Title: PET</b>	Spectrometer 3, SiKa 3/11/99									
3											
4	<b>Time</b>	<b>Beam</b>	<b>%RSD</b>	<b>Counts</b>	<b>STDEV</b>	<b>%RSD</b>	<b>CPS (x)</b>	<b>CPS/nA (y)</b>	<b>DT (usec)</b>	<b>DT (usec)</b>	
5	10.00	3.26	0.26	19580	104.88	0.54	1958	601.28	-55.23		
6	10.00	7.07	0.14	38107	165.91	0.44	3811	538.94	1.77		
7	10.00	14.54	0.03	74496	189.21	0.25	7450	512.25	7.51		
8	10.00	29.92	0.04	147260	211.60	0.14	14726	492.16	6.31		
9	10.00	50.56	0.04	240663	198.51	0.08	24066	475.99	5.10	2.66	
10	10.00	80.88	0.06	370599	226.63	0.06	37060	458.22	4.20	2.67	
11	10.00	103.23	0.00	460498	424.50	0.09	46050	446.09	3.86	2.67	
12	10.00	153.84	0.02	647875	755.63	0.12	64787	421.14	3.46	2.65	
13	10.00	199.52	0.05	799632	764.97	0.10	79963	400.78	3.27	2.65	
14	10.00	248.67	0.17	945433	361.50	0.04	94543	380.19	3.17	2.67	
15											
16	Least Squares Regression for all data points										
17	Slope:	-0.0019					<b>Mean DT</b>		4.29	2.66	
18	Y-intercept:	542.606					<b>STDEV</b>		1.76	0.01	
19	Least Squares Regression for last 6 points (more precise)										
20	Slope:	-0.0014									
21	Y-intercept:	508.526					<b>Regression Deadtime</b>		<b>3.44</b>	<b>2.66</b>	
22											
23											

The following screen capture illustrates the raw data template.

The screenshot shows a Microsoft Excel window titled "Microsoft Excel - Deadtime\_Calc.xls". The menu bar includes File, Edit, View, Insert, Format, Tools, Data, Window, Help, and Adobe PDF. The toolbar contains various icons for file operations and editing. The active cell is A24, containing the formula =Title: PET, Spectrometer 3, 3/11/99, SiKa. The spreadsheet data is as follows:

	A	B	C	D	E	F	G
22							
23	<b>Raw Data Template</b>						
24	<b>Title: PET,</b> Spectrometer 3, 3/11/99, SiKa						
25	<b>Time(sec)</b>	<b>Beam(nA)</b>	<b>Counts(cts)</b>				
26	10	3.259	19470				
27	10	3.263	19590				
28	10	3.247	19679				
29	10	7.08	37952				
30	10	7.06	38087				
31	10	7.072	38282				
32	10	14.547	74572				
33	10	14.539	74281				
34	10	14.543	74636				
35	10	29.911	147492				
36	10	29.917	147209				
37	10	29.935	147078				
38	10	50.539	240665				
39	10	50.582	240860				
40	10	50.58	240463				
41	10	80.844	370608				
42	10	80.856	370821				
43	10	80.933	370368				
44	10	103.225	460976				
45	10	103.229	460353				
46	10	103.235	460165				
47	10	153.811	647802				
48	10	153.828	647158				
49	10	153.871	648664				
50	10	199.604	800511				
51	10	199.545	799268				
52	10	199.402	799117				
53	10	248.192	945455				
54	10	248.884	945061				
55	10	248.947	945783				
56							
57							
58							

Calculate the deadtime factor for each spectrometer in turn, by overwriting the last column of count data in the raw data template portion of the Excel spreadsheet. Simply highlight the data in the Excel linked spreadsheet (from STARTWIN), use the copy function and paste it into the appropriate column. Edit cells **A2** and **A24** to update the title of the spreadsheet, for documentation and printout purposes. Calculations on the new data set will be automatically updated and output.

The deadtime constants are placed into the SCALARS.DAT file (line 13). Enter a value for each spectrometer (units of microseconds, as output from the Excel spreadsheet).

Deadtime may not be a constant and probably varies with the line energy of the x-ray being measured. One way to get around this is to place a pulse stretching circuit before the counter timer board to ensure that a forced deadtime is used to mask the actual deadtime range of the spectrometer. A pulse width (from the pulse stretcher) greater than the worst case deadtime found for the spectrometer is produced. Using this value will lead to a more accurate deadtime correction at all energies.

(This page intentionally left blank)



# Calculation Options

Prior to analyzing collected x-ray data, the user may wish to specify various output calculation options. These choices may be found by clicking the **Calculation Options** button in the **Analyze!** window.

**Analyze!**

Sample List (multi-select) (double-click to see intensity data)

- Standards
- Unknowns
- Wavescans
- All Samples

Select All

Add To Setup

Save Setups

Un 5 Gore Mt. Garnet Standard

Un 6 93GH16B c1-1

Un 7 93GH16B c1-2

Un 8 93GH16B c1-3

Un 9 93GH16B c1-4

Un 10 93GH16B c1-5

Un 11 93GH16B c1-6

Un 12 93GH16B c1-7

Un 13 93GH16B c1-8

Un 14 93GH16B c1-1a

Analyze Data KRows

Combine Selected Samples >>Excel

List Report Calculation Options

Pause Between Samples

Use All Matrix Corrections

Report

Delete Selected Sample(s)

Undelete Selected Sample(s)

Match

Combine Analysis Lines From Selected Samples

Combine Data Lines From Selected Samples

Sort Stat and Data Grids In Geological or Atomic Number Order

Specified Concentrations Standard Assignments Name/Description Conditions Elements/Cations

Un 5 Gore Mt. Garnet Standard	.000	Total Oxygen	.000	Total Weight %
TO = 40, KeV = 15, Beam = 30, Size = 1	.000	Calculated Oxygen	.000	Z - Bar
(MaqAnal = 2000.), Mode = Analog Spot	.000	Excess Oxygen	.000	Atomic Weight
X-ray Counts (cps/30.43nA)				

Copy	mg ka Off	mn ka Off	ca ka Off	al ka Off	fe ka Off	ti ka Off	si ka Off	Beam
Average:	1807.9	20.8	1365.0	5792.0	999.3	14.5	2332.6	30.557
Std Dev:	9.4	.2	10.6	9.7	1.8	.3	11.2	.021
OneSigma:	6.8	.8	5.9	12.0	5.0	1.2	10.8	
Std Err:	6.6	.1	7.5	6.9	1.3	.2	7.9	
%Rel SD:	.52	.86	.78	.17	.18	2.36	.48	
Minimum:	1801.3	20.6	1357.5	5785.2	998.0	14.3	2324.7	30.542
Maximum:	1814.5	20.9	1372.5	5798.9	1000.6	14.8	2340.5	30.572

Delete Selected Line(s) Undelete Selected Line(s) Analyze Selected Line(s)

Copy	mg ka Off	mn ka Off	ca ka Off	al ka Off	fe ka Off	ti ka Off	si ka Off	Beam
22 G	1814.5	20.9	1372.5	5798.9	1000.6	14.3	2324.7	30.572
23 G	1801.3	20.6	1357.5	5785.2	998.0	14.8	2340.5	30.542

Cancel Next

The **Calculation Options** window opens.

**Calculation Options**

**Selected Samples**

Un 5 Gore Mt. Garnet Standard

**EDS Calculation Data**

Do Not Use EDS Element Data  
 Use EDS Spectrum Element Data

Assign EDS Spectral Elements

**Integrated Intensity Data Options**

Do Not Use Integrated Intensities  
 Use Integrated Intensities

**Sample Conductive Coating (need to explicitly turn on in Analytical | Analysis Options)**

Use Standard menu to specify standard coatings

Element	Density	Thickness (A)
c	2.1	200

Use Conductive Coating

**Calculations Options**

Display Results As Oxides  
 Calculate Atomic Percents  
 Calculate Detection Limits and Sensitivity  
 Calculate Projected Detection Limits  
 Calculate Homogeneity Ranges  
 Calculate Alternate Homogeneity Ranges  
 Calculate Pearson's Linear Correlation Coefficients

Calculate with Stoichiometric Oxygen  
 Calculate as Elemental

Use Particle/Film Calculations

Element By Difference: [ ]

Stoichiometry To Calculated Oxygen: [ ] Atoms Of [ ] To 1 Oxygen

Stoichiometry To Another Element: [ ] Atoms Of [ ] To [ ]

Hydrogen Stoichiometry To Excess Oxygen H:O Ratio .00 OH = 1, H2O = 2

**Formula and Mineral Calculations**

Calculate Formula Based On 24 Atoms Of 0 Add specified oxygen, etc. from the Elements/Cations button

No Mineral End-Member Calculation

Olivine  Feldspar  Pyroxene  Garnet (Ca,Mg,Fe,Mn)  Garnet (Al,Fe,Cr)

Amphibole (Ague, Auto Normalization)  Biotite (Brimhall and Ague, Halog Code)

Each of the selected options in the above window will be briefly discussed in conjunction with the data output for the selected sample *Un 5 Gore Mt. Garnet Standard*.

Un 5 Gore Mt. Garnet Standard  
 TakeOff = 40.0 KiloVolt = 15.0 Beam Current = 30.0 Beam Size = 1  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3

Std 40  
 Number of Data Lines: 2 Number of 'Good' Data Lines: 2  
 First/Last Date-Time: 11/21/1996 10:51:22 AM to 11/21/1996 10:55:35 AM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: 42.298 Average Total Weight%: 100.752  
 Average Calculated Oxygen: 42.298 Average Atomic Number: 13.652  
 Average Excess Oxygen: .000 Average Atomic Weight: 22.871  
 Average ZAF Iteration: 4.00 Average Quant Iterate: 2.00

Oxygen Calculated by Cation Stoichiometry and Included in the Matrix Correction

Un 5 Gore Mt. Garnet Standard, Results in Elemental Weight Percents

SPEC: O  
 TYPE: CALC

AVER: 42.298  
 SDEV: .027

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	LIN
TIME:	40.00	40.00	40.00	40.00	40.00	40.00	20.00
BEAM:	30.56	30.56	30.56	30.56	30.56	30.56	30.56

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si	SUM
22	5.229	.418	5.345	11.932	17.024	.065	18.474	100.765
23	5.189	.413	5.287	11.896	16.982	.067	18.588	100.740

AVER:	5.209	.415	5.316	11.914	17.003	.066	18.531	100.752
SDEV:	.028	.004	.041	.025	.029	.002	.081	.018
SERR:	.020	.003	.029	.018	.021	.001	.057	
%RSD:	.54	.85	.76	.21	.17	2.36	.44	
STDS:	40	26	40	40	40	41	40	

STKF:	.0330	.3086	.0499	.0841	.1468	.2898	.1370	
STCT:	1818.5	1807.9	1356.0	5863.6	993.9	7261.9	2319.7	

UNKF:	.0328	.0035	.0502	.0831	.1476	.0006	.1377	
UNCT:	1807.9	20.8	1365.0	5792.0	999.3	14.5	2332.6	
UNBG:	31.0	7.3	33.8	59.4	11.7	46.8	6.3	

ZCOR:	1.5892	1.1716	1.0580	1.4339	1.1517	1.1330	1.3457	
KRAW:	.9941	.0115	1.0066	.9878	1.0054	.0020	1.0055	
PKBG:	59.27	3.87	41.44	98.54	86.43	1.31	375.63	

Un 5 Gore Mt. Garnet Standard, Results in Oxide Weight Percents

SPEC: O  
 TYPE: CALC

AVER: .000  
 SDEV: .000

ELEM:	MgO	MnO	CaO	Al2O3	FeO	TiO2	SiO2	SUM
22	8.671	.539	7.478	22.545	21.901	.108	39.522	100.765
23	8.605	.533	7.398	22.478	21.848	.112	39.767	100.740
AVER:	8.638	.536	7.438	22.512	21.874	.110	39.645	100.752
SDEV:	.047	.005	.057	.048	.038	.003	.174	.018
SERR:	.033	.003	.040	.034	.027	.002	.123	
%RSD:	.54	.85	.76	.21	.17	2.36	.44	

Un 5 Gore Mt. Garnet Standard, Results in Atomic Percents

SPEC: O  
 TYPE: CALC

AVER: 60.010  
 SDEV: .026

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si	SUM
22	4.884	.173	3.027	10.040	6.921	.031	14.933	100.000
23	4.846	.170	2.994	10.007	6.902	.032	15.021	100.000
AVER:	4.865	.172	3.011	10.023	6.911	.031	14.977	100.000
SDEV:	.027	.001	.024	.023	.013	.001	.063	.000
SERR:	.019	.001	.017	.017	.009	.001	.044	
%RSD:	.56	.87	.78	.23	.19	2.34	.42	

Un 5 Gore Mt. Garnet Standard, Results Based on 24 Atoms of o

SPEC: O  
 TYPE: CALC

AVER: 24.000  
 SDEV: .000

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si	SUM
22	1.954	.069	1.211	4.016	2.769	.012	5.974	40.005
23	1.937	.068	1.197	4.001	2.759	.013	6.006	39.981
AVER:	1.946	.069	1.204	4.009	2.764	.012	5.990	39.993
SDEV:	.012	.001	.010	.011	.007	.000	.022	.017
SERR:	.008	.000	.007	.008	.005	.000	.016	
%RSD:	.61	.92	.83	.28	.24	2.30	.37	

Garnet Mineral End-Member Calculations (Ca, Mg, Fe, Mn)

	Gro	Pyr	Alm	Sp
22	20.2	32.6	46.1	1.2
23	20.1	32.5	46.3	1.1
AVER:	20.1	32.5	46.2	1.1
SDEV:	.1	.0	.1	.0

Detection limit at 99 % Confidence in Elemental Weight Percent (Single Line):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
22	.012	.039	.017	.012	.043	.023	.023
23	.012	.041	.017	.012	.044	.023	.020
AVER:	.012	.040	.017	.012	.044	.023	.021
SDEV:	.000	.001	.000	.000	.001	.000	.002
SERR:	.000	.001	.000	.000	.000	.000	.001

Percent Analytical Relative Error (One Sigma, Single Line):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
22	.4	5.1	.4	.2	.5	14.6	.5
23	.4	5.3	.4	.2	.5	14.4	.5
AVER:	.4	5.2	.4	.2	.5	14.5	.5
SDEV:	.0	.1	.0	.0	.0	.2	.0
SERR:	.0	.1	.0	.0	.0	.1	.0

Range of Homogeneity (t-test) in +/- Elemental Weight Percent (Average of Sample):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	.025	.004	.035	.023	.025	.001	.077
80ci	.055	.009	.078	.052	.056	.003	.172
90ci	.113	.019	.159	.107	.114	.007	.353
95ci	.228	.038	.321	.214	.230	.014	.710
99ci	1.142	.189	1.607	1.074	1.151	.069	3.556

Test of Homogeneity (t-test) at 1.0 % Precision (Average of Sample):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	yes	yes	yes	yes	yes	no	yes
80ci	no	no	no	yes	yes	no	yes
90ci	no	no	no	yes	yes	no	no
95ci	no	no	no	no	no	no	no
99ci	no	no	no	no	no	no	no

Level of Homogeneity (t-test) in +/- Percent (Average of Sample):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	.5	1.0	.7	.2	.1	2.3	.4
80ci	1.1	2.2	1.5	.4	.3	5.1	.9
90ci	2.2	4.5	3.0	.9	.7	10.4	1.9
95ci	4.4	9.1	6.0	1.8	1.4	21.0	3.8
99ci	21.9	45.5	30.2	9.0	6.8	105.2	19.2

Detection Limit (t-test) in Elemental Weight Percent (Average of Sample):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	---	.008	---	---	---	.008	---
80ci	---	.017	---	---	---	.018	---
90ci	---	.035	---	---	---	.037	---
95ci	---	.070	---	---	---	.075	---
99ci	---	.349	---	---	---	.374	---

Analytical Sensitivity (t-test) in Elemental Weight Percent (Average of Sample):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	.035	.006	.049	.033	.035	.002	.109
80ci	.078	.013	.110	.073	.079	.005	.243
90ci	.160	.027	.225	.151	.161	.010	.499
95ci	.322	.053	.454	.303	.325	.020	1.004
99ci	1.615	.267	2.272	1.519	1.628	.098	5.029

Range of Ideal Homogeneity (t-test) in +/- Elemental Weight Percent (Average of Sample)  
(Meisenkothen and Donovan):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	.026	.017	.031	.034	.116	.002	.118
80ci	.059	.038	.069	.075	.260	.004	.263
90ci	.121	.078	.141	.155	.532	.008	.539
95ci	.243	.157	.285	.311	1.072	.017	1.085
99ci	1.219	.788	1.427	1.560	5.368	.084	5.436

Range of Actual Homogeneity (t-test) in +/- Elemental Weight Percent (Average of Sample)  
(Meisenkothen and Donovan):

ELEM:	Mg	Mn	Ca	Al	Fe	Ti	Si
60ci	.035	.006	.049	.033	.035	.002	.109
80ci	.078	.013	.110	.073	.079	.005	.243
90ci	.160	.027	.225	.151	.161	.010	.499
95ci	.322	.053	.454	.303	.325	.020	1.004
99ci	1.615	.267	2.272	1.519	1.628	.098	5.029

**Calculation Options**

**Selected Samples**

Un 5 Gore Mt. Garnet Standard

**OK** **Cancel**

**EDS Calculation Data**

Do Not Use EDS Element Data  
 Use EDS Spectrum Element Data

**Assign EDS Spectral Elements**

**Integrated Intensity Data Options**

Do Not Use Integrated Intensities  
 Use Integrated Intensities

**Sample Conductive Coating (need to explicitly turn on in Analytical | Analysis Options)**

Use Standard menu to specify standard coatings

Element	Density	Thickness (A)
c	2.1	200

Use Conductive Coating

**Calculations Options**

Display Results As Oxides  
 Calculate Atomic Percents  
 Calculate Detection Limits and Sensitivity  
 Calculate Projected Detection Limits  
 Calculate Homogeneity Ranges  
 Calculate Alternate Homogeneity Ranges  
 Calculate Pearson's Linear Correlation Coefficients

Calculate with Stoichiometric Oxygen  
 Calculate as Elemental

**Use Particle/Film Calculations**

Element By Difference:

Stoichiometry To Calculated Oxygen:  Atoms Of  To 1 Oxygen  
 Stoichiometry To Another Element:  Atoms Of  To   
 Hydrogen Stoichiometry To Excess Oxygen H:O Ratio .00 OH = 1, H2O = 2

**Formula and Mineral Calculations**

Calculate Formula Based On 24 Atoms Of 0

Add specified oxygen, etc. from the Elements/Cations button

No Mineral End-Member Calculation  
 Olivine  Feldspar  Pyroxene  Garnet (Ca,Mg,Fe,Mn)  Garnet (Al,Fe,Cr)  
 Amphibole (Ague, Auto Normalization)  Biotite (Brimhall and Ague, Halog Code)

Selecting *Display Results As Oxides* check box permits the user to display the results of an analysis in oxide weight percents based on the cation ratios defined for each element in the **Element/Cations** dialog window. Results are also, always reported in elemental weight percents.

The *Calculate with Stoichiometric Oxygen* button allows the user to calculate oxygen by stoichiometry if oxygen is not an analyzed element in the routine. If oxygen is either measured or calculated by stoichiometry and the *Display Results As Oxides* check box is selected, then the program will automatically calculate and report the actual excess or deficit oxygen in the analysis. This information can be very useful in determining if the selected cation ratios are correct (iron bearing oxides, for example).

All elements to be calculated by stoichiometry, difference or formula basis must be listed in the sample setup. Add these elements using the **Elements/Cations** button. Each must be added as a “not analyzed” element; click any empty row in the element list, type in the element symbol and leave the x-ray line blank. Analyses can also be output as atomic percents if the *Calculate Atomic Percents* check box is marked. This calculation is based on the fraction of the atomic weight of each element and is normalized to a 100% total.

The *Formula and Mineral Calculation* fields at the base of the **Calculation Options** window allow the user to compute formulas based on any number of oxygens for oxide runs or any analyzed or specified element in elemental runs. Further, olivine, feldspar, pyroxene, and two garnet end-member calculations are written into the software. These formula calculations are based only on atomic weight and do not consider charge balance and site occupancy. See the appendix sections in *An introduction to The Rock-Forming Minerals* by Deer, Howie, and Zussman (1992) for details on calculating formulas for hydrous phases.

The user may also select the *Calculate Detection Limits and Homogeneity* check box. The calculation of the sample detection limits is based on the standard counts, the unknown background counts, and includes the magnitude of the ZAF correction factor. The calculation is adapted from Scott et al., (1995). This detection limit calculation is useful in that it can be used even on inhomogenous samples and can be quoted as the detection limit in weight percent for a single analysis line with a confidence of 99% (assuming 3 standard deviations).

$$C_{\text{MDL}} = (\text{ZAF}) \frac{3\sqrt{I_B/t}}{I_S} \cdot 100$$

Where:

ZAF	is the ZAF correction factor for the sample matrix
$I_S$	is the count rate on the analytical (pure element) standard
$I_B$	is the background count rate on the unknown sample
t	is the counting time on the unknown sample

After this, a rigorous calculation of the analytical error also for single analysis lines, is performed based on the peak and background count rates (Scott et al., 1995). The results of the calculation are displayed after multiplication by a factor of 100 to give a percent analytical error of the net



count rate. This analytical error result can be compared to the percent relative standard deviation (%RSD) displayed in the analytical calculation. The analytical error calculation is as follows:

$$\mathcal{E}_{P-B} = \sqrt{\frac{N_p}{t_p^2} + \frac{N_B}{t_B^2}} \left/ \left( \frac{N_p}{t_p} - \frac{N_B}{t_B} \right) \right.$$

Where:  $N_p$  is the total peak counts  
 $N_B$  is the total background counts  
 $t_p$  is the peak count time  
 $t_B$  is the background count time

A more comprehensive set of calculations for analytical statistics will also be performed. These statistics are based on equations adapted from *Scanning Electron Microscopy and X-Ray Microanalysis, Second Edition* by Goldstein, et al., (1992). All calculations are expressed for various confidence intervals from 60 to 99% confidence.

The calculations are based on the number of data points acquired in the sample and the measured standard deviation for each element. This is important because although x-ray counts theoretically have a standard deviation equal to square root of the mean, the actual standard deviation is usually larger due to variability of instrument drift, x-ray focusing errors, and x-ray production. The statistical calculations include:

The range of homogeneity in plus or minus weight percent.

$$W_{1-\alpha} = \pm C \left( \frac{t_{n-1}^{1-\alpha}}{n^{1/2}} \right) \frac{S_C}{\bar{N}}$$

The level of homogeneity in plus or minus percent of the concentration.

$$\pm \frac{W_{1-\alpha}}{C} = \pm \frac{(t_{n-1}^{1-\alpha}) S_C (100)}{n^{1/2} \bar{N}}$$

The trace element detection limit in weight percent.

$$C_{DL} = \frac{C_S}{\bar{N}_S - \bar{N}_{SB}} \frac{2^{1/2} (t_{n-1}^{1-\alpha}) S_C}{n^{1/2}}$$

The analytical sensitivity in weight percent.

$$\Delta C = C - C' \geq \frac{2^{1/2} C (t_{n-1}^{1-\alpha}) S_C}{n^{1/2} (\bar{N} - \bar{N}_B)}$$

Where:  $C'$  is the concentration to be compared with  
 $C$  is the actual concentration in weight percent of the sample

$C_S$  is the actual concentration in weight percent of the standard  
 $t_{n-1}^{1-\alpha}$  is the Student t for a  $1-\alpha$  confidence and  $n-1$  degrees of freedom  
 $n$  is the number of data points acquired  
 $S_C$  is the standard deviation of the measured values  
 $\bar{N}$  is the average number of counts on the unknown  
 $\bar{N}_B$  is the continuum background counts on the unknown  
 $\bar{N}_S$  is the average number of counts on the standard  
 $\bar{N}_{SB}$  is the continuum background counts on the standard

The homogeneity test compares the 99% confidence range of homogeneity value with 1% of the sample concentration for each element. If the range of homogeneity is less than 1% of the sample concentration then the sample may be considered to be homogenous within 1%. The detection limit calculation here is intended only for use with homogenous samples since the calculation includes the actual standard deviation of the measured counts. This detection limit can, however, be quoted for the sample average and of course will improve as the number of data points acquired increases. Note that the homogenous sample detection limit calculation are ignored for those elements which occur as minor or major concentrations (>1%).

Conversely, the analytical sensitivity calculation is ignored for elements whose concentrations are present at less than 1%.

Three other calculation options are available to the user; *Element By Difference*, *Stoichiometry To Calculated Oxygen*, and *Stoichiometry To Another Element*.

When the *Element By Difference* check box is selected, the user can include an element in the analysis to account for their affect on the other x-ray intensities. This element must be specified in the sample setup. Note this method causes the calculation to result in a 100% total.

The *Stoichiometry To Calculated Oxygen* option is often used in the analysis of carbonate or borate samples in an oxide run. This feature permits the user to analyze just the cations in the sample and have oxygen calculated by stoichiometry and another specified element (usually C in carbonates and B in borates) calculated relative to oxygen. In the carbonate scenario ( $\text{CaCO}_3$ ), carbon is always in the ratio of 1 to 3 to oxygen. If the user specifies carbon by stoichiometry relative to the stoichiometric element oxygen at 0.333 (1 divided by 3) the correct amount of both carbon and oxygen will be incorporated into the ZAF matrix correction and totals without analyzing for either element. This method should only be used with phases where the ratio to oxygen is both fixed and known.

**Calculation Options**

**Selected Samples**

Un 4 Siderite Standard

**EDS Calculation Data**

Do Not Use EDS Element Data  
 Use EDS Spectrum Element Data

Assign EDS Spectral Elements

**Integrated Intensity Data Options**

Do Not Use Integrated Intensities  
 Use Integrated Intensities

**Sample Conductive Coating (need to explicitly turn on in Analytical | Analysis Options)**

Use Standard menu to specify standard coatings

Element	Density	Thickness (A)
c	2.1	200

Use Conductive Coating

**Calculations Options**

Display Results As Oxides  
 Calculate Atomic Percents

Calculate Detection Limits and Sensitivity  
 Calculate Projected Detection Limits  
 Calculate Homogeneity Ranges  
 Calculate Alternate Homogeneity Ranges  
 Calculate Pearson's Linear Correlation Coefficients

Element By Difference:

**Stoichiometry To Calculated Oxygen:** 0.333 Atoms Of c To 1 Oxygen  
 **Stoichiometry To Another Element:** Atoms Of To  
 Hydrogen Stoichiometry To Excess Oxygen H:O Ratio .00 OH = 1, H2O = 2

Calculate Formula Based On Atoms Of Sum Add specified oxygen, etc. from the Elements/Cations button

No Mineral End-Member Calculation

Olivine  Feldspar  Pyroxene  Garnet (Ca,Mg,Fe,Mn)  Garnet (Al,Fe,Cr)

Amphibole (Ague, Auto Normalization)  Biotite (Brimhall and Ague, Halog Code)

OK Cancel

The following iron carbonate mineral (siderite) output illustrates oxygen calculated by cation stoichiometry with the element carbon is calculated at 0.333 atoms relative to 1.0 atom of oxygen.

```
Un      4 Siderite Standard
TakeOff = 40.0 KiloVolt = 15.0 Beam Current = 15.0 Beam Size = 10
(Magnification (analytical) = 2000), Beam Mode = Analog Spot
(Magnification (default) = 0, Magnification (imaging) = 40)
Image Shift (X,Y): -2, 3
Number of Data Lines: 3 Number of 'Good' Data Lines: 3
First/Last Date-Time: 07/26/1999 02:45:52 PM to 07/26/1999 02:51:12 PM
WARNING- Forcing negative k-ratios to zero
```

```
Average Total Oxygen:      41.449 Average Total Weight%: 100.042
Average Calculated Oxygen: 41.449 Average Atomic Number: 16.438
Average Excess Oxygen:     .000 Average Atomic Weight: 23.164
Average ZAF Iteration:     8.00 Average Quant Iterate: 2.00
```

Oxygen Calculated by Cation Stoichiometry and Included in the Matrix Correction  
Element C is Calculated .333 Atoms Relative To 1.0 Atom of Oxygen

Un 4 Siderite Standard, Results in Elemental Weight Percents

```
SPEC:      C      O
TYPE:      STOI   CALC

AVER:    10.356  41.449
SDEV:     .017   .013

ELEM:     Ca      Mg      Mn      Fe
BGDS:     LIN     LIN     LIN     LIN
TIME:    20.00   30.00   40.00   40.00
BEAM:    10.34   10.34   10.34   10.34

ELEM:     Ca      Mg      Mn      Fe      SUM
19      .000   .067   2.235  45.825  99.929
20      .014   .072   2.375  46.005 100.264
21      .013   .077   2.323  45.708  99.933

AVER:     .009   .072   2.311  45.846 100.042
SDEV:     .008   .005   .071   .150   .192
SERR:     .004   .003   .041   .086
%RSD:    86.77   6.95   3.06   .33
STDS:     130   131   132   132

STKF:     .3826  .0853  .0202  .4131
STCT:    3040.8 1502.9  43.1   955.3

UNKF:     .0001  .0003  .0205  .4124
UNCT:     .2     6.1   43.6   953.6
UNBG:     11.0   11.0   3.0    5.6

ZCOR:     .9873  2.0734 1.1296  1.1117
KRAW:     .0001  .0041  1.0111  .9982
PKBG:     1.03   1.55   16.50  172.56
```

Un 4 Siderite Standard, Results in Oxide Weight Percents

```
SPEC:     CO2      O
TYPE:     STOI   CALC

AVER:    37.946  .000
SDEV:     .062   .000

ELEM:     CaO      MgO      MnO      FeO      SUM
19      .000   .111   2.886  58.954  99.929
20      .019   .119   3.066  59.185 100.264
21      .018   .128   3.000  58.803  99.933
```

AVER:	.012	.119	2.984	58.981	100.042
SDEV:	.011	.008	.091	.192	.192
SERR:	.006	.005	.053	.111	
%RSD:	86.77	6.95	3.06	.33	

Another interesting example demonstrating this feature is nicely documented in the User's Guide and Reference documentation (see Stoichiometry to Oxygen section). There, several trace metals are analyzed for in a stoichiometric Al<sub>2</sub>O<sub>3</sub> matrix without measuring aluminum or oxygen, BUT the correct amount of Al<sub>2</sub>O<sub>3</sub> is added to the matrix correction!

The *Stoichiometry To Another Element* option gives the user another recalculation method similar to the *Stoichiometry To Calculated Oxygen* option just discussed. Here, the user may select any other analyzed or specified element as the stoichiometric basis element.

The example below calculates CO<sub>2</sub> on the basis of moles of CaO, rather than by stoichiometry to oxygen.

The setup is shown in the **Calculation Options** window below.

**Calculation Options**

**Selected Samples**

Un 3 Dolomite Standard

**EDS Calculation Data**

Do Not Use EDS Element Data  
 Use EDS Spectrum Element Data

Assign EDS Spectral Elements

**Integrated Intensity Data Options**

Do Not Use Integrated Intensities  
 Use Integrated Intensities

**Sample Conductive Coating (need to explicitly turn on in Analytical | Analysis Options)**

Use Standard menu to specify standard coatings

Element	Density	Thickness (A)
c	2.1	200

Use Conductive Coating

**Calculations Options**

Display Results As Oxides  
 Calculate Atomic Percents

Calculate Detection Limits and Sensitivity  
 Calculate Projected Detection Limits  
 Calculate Homogeneity Ranges  
 Calculate Alternate Homogeneity Ranges  
 Calculate Pearson's Linear Correlation Coefficients

Element By Difference:

Stoichiometry To Calculated Oxygen:  Atoms Of  To 1 Oxygen

Stoichiometry To Another Element:  Atoms Of  To

Hydrogen Stoichiometry To Excess Oxygen H:O Ratio  OH = 1, H2O = 2

Use Particle/Film Calculations

Calculate Formula Based On  Atoms Of  Add specified oxygen, etc. from the Elements/Cations button

No Mineral End-Member Calculation

Olivine  Feldspar  Pyroxene  Garnet (Ca,Mg,Fe,Mn)  Garnet (Al,Fe,Cr)

Amphibole (Ague, Auto Normalization)  Biotite (Brimhall and Ague, Halog Code)

The resulting carbonate output is seen next.

Un 3 Dolomite Standard  
 TakeOff = 40.0 KiloVolt = 15.0 Beam Current = 15.0 Beam Size = 10  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 3 Number of 'Good' Data Lines: 3  
 First/Last Date-Time: 07/26/1999 02:36:45 PM to 07/26/1999 02:42:40 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: 52.237 Average Total Weight%: 100.376  
 Average Calculated Oxygen: 52.237 Average Atomic Number: 10.882  
 Average Excess Oxygen: .000 Average Atomic Weight: 18.446  
 Average ZAF Iteration: 3.00 Average Quant Iterate: 2.00

Oxygen Calculated by Cation Stoichiometry and Included in the Matrix Correction  
 Element C is Calculated 2 Atoms Relative To 1.0 Atom of Ca

Un 3 Dolomite Standard, Results in Elemental Weight Percents

SPEC:	C	O			
TYPE:	RELA	CALC			
AVER:	13.070	52.237			
SDEV:	.042	.111			
ELEM:	Ca	Mg	Mn	Fe	
BGDS:	LIN	LIN	LIN	LIN	
TIME:	20.00	30.00	40.00	40.00	
BEAM:	10.34	10.34	10.34	10.34	
ELEM:	Ca	Mg	Mn	Fe	SUM
16	21.744	13.226	.044	.033	100.210
17	21.881	13.159	.023	.030	100.560
18	21.793	13.246	.000	.028	100.356
AVER:	21.806	13.210	.022	.030	100.376
SDEV:	.069	.045	.022	.003	.176
SERR:	.040	.026	.013	.002	
%RSD:	.32	.34	98.33	8.66	
STDS:	130	131	132	132	
STKF:	.3826	.0853	.0202	.4131	
STCT:	3040.8	1502.9	43.1	955.3	
UNKF:	.2045	.0847	.0002	.0003	
UNCT:	1625.7	1492.4	.3	.6	
UNBG:	8.6	9.0	2.1	3.4	
ZCOR:	1.0661	1.5597	1.2224	1.2016	
KRAW:	.5346	.9930	.0070	.0006	
PKBG:	189.50	167.57	1.19	1.18	

Un 3 Dolomite Standard, Results in Oxide Weight Percents

SPEC:	CO2	O			
TYPE:	RELA	CALC			
AVER:	47.890	.000			
SDEV:	.152	.000			
ELEM:	CaO	MgO	MnO	FeO	SUM
16	30.425	21.932	.056	.043	100.210
17	30.616	21.822	.030	.038	100.560
18	30.493	21.966	.000	.036	100.356
AVER:	30.511	21.907	.029	.039	100.376
SDEV:	.097	.075	.028	.003	.176
SERR:	.056	.043	.016	.002	
%RSD:	.32	.34	98.33	8.66	

(This page intentionally left blank)



# Linear Calibration - Curve Method


PROBE FOR EPMA offers a sophisticated calibration curve (multi-standard) method for correcting x-ray data. It is based on a second order polynomial fit to multiple standard intensity data. This option has been utilized in special situations such as the analysis of trace carbon in steels and when a suitable set of well characterized standards are available.

The example outlined below will document this calibration curve method for the specific analysis of Al, Ti and O in titanium aluminides doped with oxygen. This data and mdb file was generously supplied from research conducted by Jim Smith at the NASA Glenn Research Center. These low density, high strength alloys are part of an ongoing study of the transport kinetics of oxygen in these metals in conjunction with the development of superior alloys for aircraft engine gas turbine turn blades.


A series of nine titanium-aluminides (varying Ti/Al ratio) were carefully prepared, each doped with a specific concentration of oxygen, ranging from 0 to 3.21%, thereby bracketing the expected unknowns range of oxygen concentration. Each standard alloy was analyzed by other techniques to verify the nominal compositions. The nine standard compositions were then entered into the STANDARD.MDB database and the positions of each standard digitized in PROBE FOR EPMA.

Aluminum, titanium and oxygen were peaked on the appropriate standards and count rate data (five spots each) were acquired on each of the nine standards. The count rate data was then examined in the **Analyze!** window to ascertain the precision of the five data points on each standard, deleting any selected lines as deemed appropriate.

After standard collection, the user must select the calibration curve approach as the matrix correction method. From the main PROBE FOR EPMA log window, select **Analytical | ZAF Selections** from the menu.



The **Matrix Correction Methods** window opens. Select the *Calibration Curve (multi-standard)* button for the *Correction Method*.




Click the **OK** button, returning to the main PROBE FOR EPMA log window.



All of the selected standards are analyzed and reported in the main PROBE FOR EPMA log window.

The **Plot Graph Data** window appears.




The *Show Grid Lines* check box has been marked to facilitate reading the percent and count values. This calibration curve may be printed out by clicking the **Print** button next to the **OK** button. If the user selects the *Force Black and White Print* check box in the **Plot!** window, then the corresponding output will be a black and white print, if not, then the printer will output the above gray background.

The user may evaluate the data using the **Zoom Full** capabilities (click and drag mouse over region of interest on graph) to expand the scaling. Here, in the center group, two data points clearly overlap. Placing the mouse cursor over any selected point on the graph returns the x and y values of that position (read above the **Zoom Full** button).


When finished, click the **OK** button to return to the **Plot!** window to next review the other calibration curves.

The output for aluminum is plotted similarly. Select the *Al Elemental Percents* versus *Al On Peak Counts (P+B)*, (the counts per second determined on peak) and click the **Output** button in the **Plot!** window. This curve is viewed below.




The graph may be modified by selecting any of the buttons across the top to enter the **Graph Control** module.

**Graph Control** (see tabs in display below) can be used to edit numerous graph parameters. Under the **Trends** tab, for instance, are the curve fitting options. Select *Curve Fit* under *Statistical Lines* and edit the *Curve Fit Type* and *Order* as required.



Click the **Apply Now** or **OK** button to see the changes in the **Graph Data** window.

The modified **Plot Graph Data** window returns. Here, a title, different, larger symbols, standard numbers and the second order polynomial curve have been added to the graph below.



Upon creating the previous plots, the selected standards were analyzed, each of the nine standards were reported in the main PROBE FOR EPMA log window along with the second order polynomial fit parameters.



The following output displays the log window output for one of these: St 617 Set 5 TiAl-7 standard.

```
St 617 Set 5 TiAl-7
TakeOff = 52.5 KiloVolt = 15.0 Beam Current = 30.0 Beam Size = 1
(Magnification (analytical) = 2000), Beam Mode = Analog Spot
(Magnification (default) = 0, Magnification (imaging) = 40)
Image Shift (X,Y): -2, 3
```

```
Titanium aluminide 3.03 Oxygen
Number of Data Lines: 5 Number of 'Good' Data Lines: 5
First/Last Date-Time: 12/11/2000 04:03:01 PM to 12/11/2000 04:06:28 PM
WARNING- Using Calibration Curve Matrix Corrections
WARNING- Forcing negative k-ratios to zero
```

```
Average Total Oxygen: .000 Average Total Weight%: 100.651
Average Calculated Oxygen: .000 Average Atomic Number: 20.177
Average Excess Oxygen: .000 Average Atomic Weight: 40.571
```

St 617 Set 5 TiAl-7, Results in Elemental Weight Percents

```
ELEM: Al Ti O
BGDS: LIN LIN LIN
TIME: 10.00 10.00 10.00
BEAM: 300.00 300.00 300.00
```

```
ELEM: Al Ti O SUM
978 15.452 81.618 2.986 100.056
979 15.680 81.966 3.114 100.761
980 15.864 81.880 3.053 100.796
981 15.788 81.957 3.000 100.744
982 15.747 82.245 2.908 100.900
```

```
AVER: 15.706 81.933 3.012 100.651
SDEV: .157 .224 .077 .338
SERR: .070 .100 .035
%RSD: 1.00 .27 2.56
```

```
PUBL: 15.730 81.240 3.030 100.000
%VAR: -.15 .85 -.59
DIFF: -.024 .693 -.018
```

```
UNCT: 13020.8 11807.2 329.4
UNBG: 215.5 29.2 116.0
KRAW: 1.0000 1.0000 1.0000
PKBG: 61.44 406.68 3.84
```

```
FIT1: .4671-363.4302 -.5152
FIT2: .0013 .0712 .0098
FIT3: .000000-.000003 .000003
DEV: 2.2 .9 7.1
```

The coefficients for the second order polynomial are listed last (Fit 1, Fit 2, Fit 3 and Dev). The three Fit terms represent the intercept, the slope and the second order curvature factor, respectively. The DEV term represents the total deviation (sum of the residuals) between the calculated curve and the original data. The smaller the number, the better here. The software prints a warning line if this correction method is active.

The analysis of unknown samples is straightforward. Create a new sample and collect the x-ray intensity data on the unknown. Select the **Analyze!** window and click **Analyze**.

An example is printed next.

```
Un 25 Sample 3-9 LowOx 0-30u 2u increments
TakeOff = 52.5 KiloVolt = 15.0 Beam Current = 30.0 Beam Size = 1
(Magnification (analytical) = 2000), Beam Mode = Analog Spot
(Magnification (default) = 0, Magnification (imaging) = 40)
Image Shift (X,Y): -2, 3
Number of Data Lines: 16 Number of 'Good' Data Lines: 16
First/Last Date-Time: 12/11/2000 01:27:41 PM to 12/11/2000 01:40:41 PM
WARNING- Using Calibration Curve Matrix Corrections
WARNING- Forcing negative k-ratios to zero
```

```
Average Total Oxygen: .000 Average Total Weight%: 99.869
Average Calculated Oxygen: .000 Average Atomic Number: 19.838
Average Excess Oxygen: .000 Average Atomic Weight: 39.903
```

Un 25 Sample 3-9 LowOx 0-30u 2u increments, Results in Elemental Weight Percents

```
ELEM:      Al      Ti      O
BGDS:      LIN      LIN      LIN
TIME:      10.00    10.00    10.00
BEAM:      300.24   300.24   300.24

ELEM:      Al      Ti      O      SUM
851  21.260  76.350  1.956  99.566
852  21.236  76.469  2.003  99.708
853  21.193  76.689  1.957  99.839
854  21.128  76.174  1.941  99.244
855  21.222  76.599  1.842  99.663
856  21.120  76.425  1.755  99.300
857  21.126  76.494  1.838  99.459
858  21.125  77.135  1.767  100.027
859  21.142  76.686  1.791  99.619
860  21.153  76.835  1.796  99.784
861  21.139  77.248  1.703  100.089
862  21.256  77.439  1.685  100.379
863  21.190  77.360  1.833  100.383
864  21.216  77.503  1.779  100.499
865  21.286  77.127  1.643  100.055
866  21.264  77.541  1.482  100.287

AVER:      21.191  76.880  1.798  99.869
SDEV:      .058    .454    .133    .393
SERR:      .015    .113    .033
%RSD:      .28     .59     7.42

UNCT:      18339.0  11021.6  222.4
UNBG:      218.2     26.6     113.8
KRAW:      1.0000    1.0000    1.0000
PKBG:      85.05     417.09    2.96

FIT1:      .4718-363.4019  -.5152
FIT2:      .0013    .0712    .0098
FIT3:      .0000000-.0000003 .0000003
DEV:      2.2     .9     7.1
```

The second order polynomial coefficients are always listed last (just above). Elements calculated by difference or stoichiometry can be calculated along with calibration curve corrected elements. And if both off-peak and MAN acquired data are present, PROBE FOR EPMA will construct separate sets of calibration curves and fit a second order polynomial expression that is used in the iteration procedure to determine the concentration of the element.

The user may elect to run standards after completing unknown sample acquisition and then correct for any standard intensity drift. From the main PROBE FOR EPMA log window select **Analytical** and then choose **Analysis Options** from the drop-down menu.

The **Analysis Calculation Options** window appears, remember to check that the *Use Automatic Drift Correction on Standard Intensities* check box is marked.

**Analysis Calculation Options**

**Quantitative Acquisition Options**

- Use Deadtime Correction
  - Use Normal Deadtime Correction (single term factorial)
  - Use Precision Deadtime Correction (two term factorial for > 50K cps)
- Use Beam Drift Correction
- Use Automatic Drift Correction on Standard Intensities

**Quantitative Analysis Options**

- Use Assigned Interference Corrections on Standards and Unknowns
  - Do Not Use Matrix Correction Term For Full Quant Interference Correction
  - Do Not Use Full Quant Interference Corrections (use traditional Gilfrich, et. al.)
- Use Assigned or Self Time Dependent Intensity (TDI) Corrections on Unknowns
  - Use Linear Fit (slope coefficient only) for TDI Extrapolation
  - Use Quadratic Fit (two coefficient) for TDI Extrapolation
- Use Time Weighted Data for TDI Fit (weight intensities based on elapsed time)
  - Time Weighted Data Weight Factor
- Use Absorption Corrected MAN Continuum Intensities
- Use Particle or Thin Film Correction Parameters
- Check For Same Peak Positions in Unknown and Standard
- Check For Same PHA Settings in Unknown and Standard
- Use Zero Point For Calibration Curve (off-peak elements only)
- Use Conductive Coating Correction For Beam Energy Loss
- Use Conductive Coating Correction For X-ray Absorption

**MAC (mass absorption coefficient) and APF (area peak factor) Options**

- Use Empirical MAC Values
- Use Empirical APF Values
  - Use Empirical APF Factors (calculated from elemental composition)
  - Use Specified APF Factors (based on a fixed composition)

**Calculation Options**

- Use Aggregate Intensities for Duplicate Quantitative Elements
- Use Blank Calibration Sample Trace Element Accuracy Corrections
- Force Negative K-Ratios To Zero in ZAF Calculations
- Calculate Electron and Xray Ranges for Sample Compositions
- Use Oxygen From Halogens (F, Cl, Br and I) Correction
- Use Nth Point Calculation For Off-Peak Intensities (for testing only)
- Use Count Overwrite Intensity Table for Data Calculations
- Force Negative Interference Intensities To Zero in Corrections
- Use Chemical Age Calculation (U, Th, Pb)

**Formatting Options**

- Use Automatic Format For Quantitative Results
  - Display the Maximum Number of Numerical Digits
  - Display Only Statistically Significant Number of Numerical Digits
- Use Detailed Printout For Data and Analytical Results
- Print Analyzed And Specified On Same Line
- Display Count Intensities Unnormalized To Time (in Analyze!)

**Output Options**

- Display Charge Balance Calculation

Elemental Output Sort Order For JJD-2 and HW Custom Output:

- Use Traditional Geological Sort Order (SiO<sub>2</sub>, TiO<sub>2</sub>, etc.)
- Use Low To High Atomic Number Sort Order
- Use High To Low Atomic Number Sort Order

OK Cancel

In some instances, it may be useful to add to the acquired data set a zero point (off-peak elements only) to improve the polynomial fit. To include a zero point, check the box (prior to analyzing the standards) labeled *Use Zero Point For Calibration Curve (Off-Peak Elements Only)*. This choice is also found in the **Analysis Calculation Options** window.

(This page intentionally left blank)

# Time Dependent Intensity (TDI) Corrections


PROBE FOR EPMA gives the user two powerful methods for the correction of time dependent intensity element effects unknown samples. Each correction option provides a means to automatically correct data for the time dependent loss of x-ray intensity due to continuous electron bombardment (heating damage/charge buildup) or carbon contamination that occurs during normal electron microprobe work. The time dependent intensity element extrapolation may be applied to any degradation (or enhancement) of the x-ray intensity over time. Elements most susceptible include sodium, potassium, perhaps sulfur and bound water. This correction is most useful for samples that are too small to utilize a defocused beam and also allows the operator to run higher than normal sample currents to improve analytical sensitivity.

Two different methods are available for volatile element corrections; the *Self Calibration Time Dependent Intensity (TDI) Acquisition* and the *Assigned Calibration Time Dependent Intensity (TDI) Acquisition*. Each will be documented.

In the *Self Calibration Time Dependent Intensity (TDI) Acquisition* method, the program acquires the on-peak count data for the volatile element **during** the normal sample acquisition process for the unknown sample. This method works well when samples to be analyzed have widely differing compositions as the calibration is determined with every sample acquisition.

Open PROBE FOR EPMA, and proceed through the normal calibration and standardization routine. Check suitable standards for accuracy, these should not be volatile or beam sensitive.

Move to your first unknown sample. Open the **New Sample** window from the **Acquire!** window and create a new unknown sample.



Click the **Special Options** button in the **Acquire!** window.

The **Volatile Calibration and Quick Wavescan Samples** window opens. Select the *Self Calibration Time Dependent Intensity (TDI) Acquisition* option. The time dependent intensity element correction is only applied to the first element analyzed for on each spectrometer. Enter a number into the *Volatile Count Time Intervals* text field (up to 50 intervals may be used). The program will use this interval number and the on-peak count time to create a calibration curve. In this example, Na is counted for 40 seconds on peak and with an interval of 10 entered, the program will automatically count ten 4 second intervals. Each element listed first on each spectrometer is treated in this way. The off peak counts are not affected.

**Volatile Calibration and Quick Wavescan Samples**

**Special Sample Acquisition Options**

**Normal Acquisition**  
Acquire a normal standard, unknown or wavescan sample.

**Self Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample that uses a calibration curve based on itself. The TDI assignment will be made automatically. This method will append a new analysis line for each acquisition on an existing unknown sample.

Note: time dependent intensity (TDI) elements must be the first element acquired on each spectrometer. The acquisition order may be modified in the Acquisition Options dialog.

**Time Dependent Intensity (TDI) Count Time Intervals**

**Acquire TDI Data on Standard Samples**

**Consider using "synchronous" spectrometer acquisition for TDI to minimize beam exposure before starting count integration.**

**Assigned Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample for use as a TDI calibration curve for another sample(s). This assignment is made subsequently in the Standard Assignments dialog in the Analyze! window. This method will create a new time dependent intensity calibration sample for each acquisition.

**Time Dependent Intensity (TDI) Sample Name**

**Time Dependent Intensity (TDI) Count Time Intervals**

**Stage X Increment (um)**

**Stage Y Increment (um)**

**Quick Wavescan Acquisition**  
Quick wavescan samples are created based on the first element on each spectrometer used in the last unknown sample (or wavescan sample if present) or a selected sample setup. This method will create a new quick wavescan sample for each acquisition.

**Quick Wavescan Sample Name**

**Quick Scan Speed % (.001 - 100)**

**Use ROM Based Spectrometer Scan**

**OK**  
**Cancel**  
**Setups**

Click the **OK** button to return to the **Acquire!** window.

Click the **Start Standard or Unknown Acquisition** button to initiate count data collection on the unknown sample.

The screenshot shows the 'Acquire!' software interface. At the top, there are columns labeled 1, 2, 3, X, Y, Z, and W. Below these are several rows of numerical data. A 'Current Sample' field contains 'Un 3 \* Glass V-1-1'. A prominent yellow button labeled 'Start Standard or Unknown Acquisition' is highlighted. Below this are fields for 'Data Rows: 0' and 'Good Data Rows: 0', and a 'Start Wavescan' button. A grid of buttons includes 'New Sample', 'PHA', 'Move', 'Stage', 'Elements/Cations', 'Peak/Scan Options', 'Acquisition Options', 'Imaging', 'Analytical Conditions', 'Count Times', 'Special Options', 'Peaking Options', 'Combined Conditions', 'Locate', 'Rate Meter', and 'Start Peaking'. On the right side, there is a 'Progress/Beam Deflection' section with a red crosshair and numerical values for 'um' and 'px'. Below that is a table of acquisition parameters: Magnification (400), Beam Mode (Analog Spot), Kilovolts (15), Beam Current (20), and Beam Size (30).

1	2	3	X	Y	Z	W
228.065	129.262	191.143	8.00024	29.9993	11.0996	4.00000
<b>1-PET</b>	<b>2-TAP</b>	<b>3-LIF</b>	<b>Faraday</b>			
4.00	4.00	8.00	1.00			
60.	59.	209.	1999.98			

Current Sample: Un 3 \* Glass V-1-1

Self Time Dependent Intensity (TDI) Correction

Data Rows: 0    Good Data Rows: 0

Start Wavescan

**Start Standard or Unknown Acquisition**

New Sample	PHA	Move	Stage
Elements/Cations	Peak/Scan Options	Acquisition Options	Imaging
Analytical Conditions	Count Times	Special Options	Peaking Options
Combined Conditions	Locate	Rate Meter	Start Peaking

Progress/Beam Deflection


3.87630    27.9220

um .000000    .000000

px 0    0

Magnification	400
Beam Mode	Analog Spot
Kilovolts	15
Beam Current	20
Beam Size	30
Acquire String	

After collecting an unknown sample the user may display the volatile correction from the **Analyze!** window.



Click the **Standard Assignments** button.



The **Standard and Interference Assignments** window opens.

**Standard and Interference Assignments**

**Selected Samples**

Un 3 Glass V-1-1

OK Cancel

Save Element Setup

Save Sample Setup

Add/Remove Standards

Reload Standard Assignments

Remove TDI Correction

1 2 3 4 5 6

**Click Element Row to Edit Standard/Interference/Time Dependent Intensity (TDI) Assignments**

Channel	Element	X-Ray	Analyzed	Standard	Interf-Ele	Interf-Std
1	si	ka	Yes	128	----	0,0,0,0,0
2	ti	ka	Yes	212	----	0,0,0,0,0
3	al	ka	Yes	81	----	0,0,0,0,0
4	fe	ka	Yes	128	----	0,0,0,0,0
5	mn	ka	Yes	220	----	0,0,0,0,0
6	mg	ka	Yes	128	----	0,0,0,0,0
7	ca	ka	Yes	128	----	0,0,0,0,0
8	na	ka	Yes	81	----	0,0,0,0,0
9	k	ka	Yes	116	----	0,0,0,0,0

Select the element row (sodium in this example) to view the TDI calibration data.

The **Assignment Properties** window opens and lists all possible TDI element calibration samples.

**Assignment Properties**

Enter Standard Assignments for: na ka

Element	X-Ray	Assigned (Primary) Standard
na	ka	81 Albite

Use Virtual Standard For Standard Intensity Calculation (see Analytical menu)

OK

Cancel

---

Interference Standard Assignments for Interfered Element: na ka

Intf Elem	Intf Order	Interference Standard	
1st			Remove
2nd			Remove
3rd			Remove
4th			Remove
5th			Remove

Check All Interfering Elements

Calculate Interferences

Minimum Overlap Intensity: .1

The standard used for the interference correction must contain a known concentration of the interfering element and none of the interfered element, nor any other interfering elements.

---

Time Dependent Intensity (TDI) Calibration Assignment (select unknown sample for calibration)

**TDI Correction Type (Self or Assigned)**

No TDI Calibration Correction  
 Use TDI "Self" Calibration Correction  
 Use TDI "Assigned" Calibration Correction

Both "assigned" and "self" calibration Time Dependent Intensity (TDI) element samples can be acquired. See the Special Options dialog in the Acquire window.

Both "assigned" and "self" Time Dependent Intensity (TDI) element calibrations can be assigned or unassigned here. Assigned Time Dependent Intensity (TDI) corrections are assigned to samples acquired with the "assigned" flag in Special Options. "Self" Time Dependent Intensity (TDI) corrections are assigned to themselves. Time Dependent Intensity (TDI) element "self" calibrations are automatically assigned to themselves at the time of acquisition.

Display TDI Fit

Use Log-Linear (exponential) Fit  
 Use Log-Quadratic (hyper-exponential) Fit

---

Blank Correction Sample Assignment

Un	2	Test Calibration—Rhyolite Glass
Un	3	Glass V-1-1


Assign a sample to be used for a "blank" trace element correction. The blank sample should be a similar matrix to the unknown sample and should have a zero or known trace of the element present.

Blank Level (wt. %): .000000

Remove Blank Assignment

Highlight the appropriate *Time Dependent Intensity (TDI) Calibration Assignment* sample and click the **Display TDI Fit** button. Two types of fits are possible; Log-Linear or Log -Quadratic.

The **Time Dependent Intensity (Self) Fit Plot** window opens.



The drop in sodium count intensity with time appears to fit an exponential function (Nielson and Sigurdsson, 1981). A plot of the natural log of the intensity data versus time should yield a straight line relationship as seen above.


The extrapolation correction is quantitative in that the correction is based on a self calibration curve acquired during the run. The correction uses the actual elapsed time for all calculations. It is applied during the ZAF or Phi-Rho-Z iteration phase of the analysis to correct for changes in the matrix correction due to the extrapolation correction.

The next example will illustrate the *Assigned Calibration Time Dependent Intensity Acquisition* method. Here, small rhyolitic glass inclusions of similar composition will be analyzed.

Start PROBE FOR EPMA in the normal manner. Go through the calibration and standardization process, then check standards. Save the analysis routine as a sample setup. Next, create a separate sample setup with a subset of elements to which the TDI correction will be assigned, in this example silicon and sodium. Note TDI element calibration corrections can only be applied to elements that are the first element collected by each spectrometer.

Locate the sample to obtain the TDI element correction acquisition upon. This should be either the unknown sample or a material similar to the unknowns.

Click the **Special Options** button in the **Acquire!** window.



The **Volatile Calibration and Quick Wavescan Samples** window opens. The default acquisition choice is *Normal Acquisition*.

**Volatile Calibration and Quick Wavescan Samples**

**Special Sample Acquisition Options**

**Normal Acquisition**  
Acquire a normal standard, unknown or wavescan sample.

**Self Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample that uses a calibration curve based on itself. The TDI assignment will be made automatically. This method will append a new analysis line for each acquisition on an existing unknown sample.

Note: time dependent intensity (TDI) elements must be the first element acquired on each spectrometer. The acquisition order may be modified in the Acquisition Options dialog.

**Time Dependent Intensity (TDI) Count Time Intervals**   
 Acquire TDI Data on Standard Samples

**Assigned Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample for use as a TDI calibration curve for another sample(s). This assignment is made subsequently in the Standard Assignments dialog in the Analyze! window. This method will create a new time dependent intensity calibration sample for each acquisition.

**Time Dependent Intensity (TDI) Sample Name**

**Time Dependent Intensity (TDI) Count Time Intervals**   
**Stage X Increment (um)**   
**Stage Y Increment (um)**

**Quick Wavescan Acquisition**  
Quick wavescan samples are created based on the first element on each spectrometer used in the last unknown sample (or wavescan sample if present) or a selected sample setup. This method will create a new quick wavescan sample for each acquisition.


**Quick Wavescan Sample Name**

**Quick Scan Speed % (.001 - 100)**   
 Use ROM Based Spectrometer Scan

**OK**  
**Cancel**  
**Setups**

*Consider using "synchronous" spectrometer acquisition for TDI to minimize beam exposure before starting count integration.*

Click the **Setups** button in the **Volatile Calibration and Quick Wavescan Samples** window, opening the **Sample Setup Selection** dialog box. Check that the appropriate volatile setup is active.



Click the **OK** button.

Select the *Assigned Calibration Time Dependent Intensity Acquisition* method. Enter text into the *Time Dependent Intensity Sample Name* field. Edit the *TDI Count Time Intervals* (number of steps in the calibration line) and adjust the *Stage X and Y Increment (um)* values if required, resulting in the following window.

**Volatile Calibration and Quick Wavescan Samples**

**Special Sample Acquisition Options**

**Normal Acquisition**  
Acquire a normal standard, unknown or wavescan sample.

**Self Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample that uses a calibration curve based on itself. The TDI assignment will be made automatically. This method will append a new analysis line for each acquisition on an existing unknown sample.

Note: time dependent intensity (TDI) elements must be the first element acquired on each spectrometer. The acquisition order may be modified in the Acquisition Options dialog.

**Time Dependent Intensity (TDI) Count Time Intervals**   
 Acquire TDI Data on Standard Samples

**Assigned Calibration Time Dependent Intensity (TDI) Acquisition**  
Acquire a time dependent intensity (TDI) element sample for use as a TDI calibration curve for another sample(s). This assignment is made subsequently in the Standard Assignments dialog in the Analyzel window. This method will create a new time dependent intensity calibration sample for each acquisition.

**Time Dependent Intensity (TDI) Sample Name**

**Time Dependent Intensity (TDI) Count Time Intervals**   
**Stage X Increment (um)**   
**Stage Y Increment (um)**

**Quick Wavescan Acquisition**  
Quick wavescan samples are created based on the first element on each spectrometer used in the last unknown sample (or wavescan sample if present) or a selected sample setup. This method will create a new quick wavescan sample for each acquisition.

**Quick Wavescan Sample Name**

**Quick Scan Speed % (.001 - 100)**   
 Use ROM Based Spectrometer Scan

**Consider using "synchronous" spectrometer acquisition for TDI to minimize beam exposure before starting count integration.**

**OK**  
**Cancel**  
**Setups**

The TDI element effect will be calibrated precisely on the first unknown sample. One important assumption is that the TDI element behavior on the calibration sample is similar to all the unknowns to be analyzed. Therefore, because a consistent TDI correction is used, variation in composition represents real differences in composition (or volatilization) not precision of the analyses. Each element in the method (sample setup) will be acquired one element at a time. In this example, 10 time intervals are specified and the default on-peak count times for silicon and sodium are 20 and 40 respectively. Thus the calibration curve for each will be composed of 10 spots of 2 and 4 seconds each. The program will acquire each element in the sample (at a fresh spot) at a new stage position based on the X and Y increments specified to allow the volatile element effect to be calibrated precisely.


Click the **OK** button to return to the **Acquire!** window.

The TDI sample acquisition will start automatically upon clicking the **Start Standard or Unknown Acquisition** button, using the sample name entered in the *TDI Sample Name* field.

The user may now collect unknown data points. **REMEMBER** to load in the appropriate sample setup **AND** under the **Special Options** button of the **Acquire!** window, click the *Normal Acquisition* button.



After collecting an unknown sample the user may assign the volatile correction from the **Analyze!** window.



The volatile element calibration can be assigned to any element in a sample provided that it was acquired as the **FIRST** element on that spectrometer.

Click the **Standard Assignments** button. The **Standard and Interference Assignments** window opens.

**Standard and Interference Assignments**

**Selected Samples**

Un 5 Unk Inclusion M-1-1

OK Cancel

Save Element Setup

Save Sample Setup

Add/Remove Standards

Reload Standard Assignments

Remove TDI Correction

1 2 3 4 5 6

**Click Element Row to Edit Standard/Interference/Time Dependent Intensity (TDI) Assignments**

Channel	Element	X-Ray	Analyzed	Standard	Interf-Ele	Interf-Std
1	si	ka	Yes	128	---	0,0,0,0,0
2	na	ka	Yes	81	---	0,0,0,0,0

To evaluate the correction effect on sodium, click that element's row to view and edit its' TDI assignment.

The **Assignment Properties** window opens and lists all possible volatile element calibration samples.

**Assignment Properties**

Enter Standard Assignments for: na ka

Element	X-Ray	Assigned (Primary) Standard
na	ka	81 Albite

Use Virtual Standard For Standard Intensity Calculation (see Analytical menu)

---

Interference Standard Assignments for Interfered Element: na ka

Intf Elem	Intf Order	Interference Standard	
1st			Remove
2nd			Remove
3rd			Remove
4th			Remove
5th			Remove

Check All Interfering Elements

Calculate Interferences

Minimum Overlap Intensity: .1

The standard used for the interference correction must contain a known concentration of the interfering element and none of the interfered element, nor any other interfering elements.

---

Time Dependent Intensity (TDI) Calibration Assignment (select unknown sample for calibration)

**TDI Correction Type (Self or Assigned)**

No TDI Calibration Correction  
 Use TDI "Self" Calibration Correction  
 Use TDI "Assigned" Calibration Correction

Un 4 Rhyolite Unknown for Calibration Curve

Both "assigned" and "self" calibration Time Dependent Intensity (TDI) element samples can be acquired. See the Special Options dialog in the Acquire window.

Both "assigned" and "self" Time Dependent Intensity (TDI) element calibrations can be assigned or unassigned here. Assigned Time Dependent Intensity (TDI) corrections are assigned to samples acquired with the "assigned" flag in Special Options. "Self" Time Dependent Intensity (TDI) corrections are assigned to themselves. Time Dependent Intensity (TDI) element "self" calibrations are automatically assigned to themselves at the time of acquisition.

Display TDI Fit

Slope of Log Counts= -3.514931E-03, Ln Intercept= 5.80298, Intercept= 331.285

Use Log-Linear (exponential) Fit  
 Use Log-Quadratic (hyper-exponential) Fit

---

Blank Correction Sample Assignment

Un 2 Glass std check  
 Un 3 Volatile Element setup  
 Un 4 Rhyolite Unknown for Calibration Curve  
 Un 5 Unk Inclusion M-1-1  
 Un 11 sample 1


Assign a sample to be used for a "blank" trace element correction. The blank sample should be a similar matrix to the unknown sample and should have a zero or known trace of the element present.

Blank Level (wt. %): .000000

Remove Blank Assignment

Highlight the appropriate *Time Dependent Intensity Calibration Assignment* sample and click the **Display TDI Fit** button.

The **Time Dependent Intensity (Assigned) Fit Plot** window opens.



The drop in sodium count intensity with time appears to fit an exponential function (Nielson and Sigurdsson, 1981). A plot of the natural log of the intensity data versus time should yield a straight line relationship as seen above.

The extrapolation correction is quantitative in that the correction is based on a calibration curve acquired during the run. The correction uses the actual elapsed time for all calculations. It is applied during the ZAF or Phi-Rho-Z iteration phase of the analysis to correct for changes in the matrix correction due to the extrapolation correction.

Along with alkali loss, the operator may notice an increase in count intensity from the non-mobile elements (such as silicon and aluminum) in the sample. The possible enhancement with time may be corrected for as well.

**Close the Time Dependent Intensity (Assigned) Fit Plot window, returning to the Assignment Properties window. Click the OK button.**

To evaluate the possible enhancement of intensity with time on silicon in the rhyolitic glass, click on the element row for silicon in the **Standard and Interference Assignments** window.

**Standard and Interference Assignments**

**Selected Samples**

Un 5 Unk Inclusion M-1-1

**OK** **Cancel**

**Save Element Setup**

**Save Sample Setup**

**Add/Remove Standards**

**Reload Standard Assignments**

**Remove TDI Correction**


1 2 3 4 5 6

**Click Element Row to Edit Standard/Interference/Time Dependent Intensity (TDI) Assignments**

Channel	Element	X-Ray	Analyzed	Standard	Interf-Ele	Interf-Std
1	si	ka	Yes	128	----	0,0,0,0,0
2	na	ka	Yes	81	----	0,0,0,0,0

The **Assignment Properties** window opens. Again, highlight the appropriate *TDI Calibration Sample Assignment* sample and click the **Display TDI Fit** button.

The **Time Dependent Intensity (Assigned) Fit Plot** window opens displaying the natural log of silicon x-ray intensity versus time. The enhancement is evident and maybe corrected for quantitatively.



Return to the **Analyze!** window by clicking the **Close** button above.

Click the **OK** button of the **Assignment Properties** and the **Standard and Interference Assignments** windows, respectively.

## Advanced Interference Corrections

PROBE FOR EPMA permits the user to select a fully quantitative correction for spectral interferences (Donovan et al., 1993). The user can correct for up to four interfering elements per channel. The program requires that both the interfered and interfering elements be analyzed for. Further, an interference calibration standard must be acquired that contains a major concentration of the interfering element and none of the interfered element nor any other elements that interfere with the interfered element.

Most interferences between a pair of elements work in one direction. Consider a phase with high concentrations of manganese in the presence of the iron. Here the Mn  $K\beta$  line interferes with the Fe  $K\alpha$  analytical line. The reverse situation does not cause a problem, iron does not interfere with manganese. However, there are some cases where both elements interfere with each other! These dual interfering elements are extremely difficult to correctly quantify. Fortunately, PROBE FOR EPMA can handle this situation because its quantitative interference correction is an iterated solution (see Donovan et al., 1993 for details).

The following example (analyzed at U.C. Berkeley by John Donovan) involves the dual interference of zinc and rhenium in a natural organo-metallic phase. Both lines interfere with each other (Zn  $K\alpha$  and Re  $L\alpha$ ) and both lines are used for quantitative analysis. Other elements analyzed for are cesium, iron, sulfur, and selenium. Oxygen, nitrogen, carbon and hydrogen are also in the samples. Each is listed in the **Element/Cations** window for use in the matrix correction routine but are not analyzed.

Solving these interferences requires the analyses of both rhenium and zinc and two interference standards. To correct for the interference on zinc, a standard that contains rhenium but no zinc is required. Likewise, to handle the interference on rhenium, a standard will be needed that contains zinc but no rhenium.

The procedure to specify interferences and the correction thereof was documented in the silicate chapter in the Users Guide to Getting Started manual, cogent details will briefly mentioned here.

To view the various interference assignments, click the **Standard Assignment** button in the **Analyze!** window.

The **Standard and Interference Assignments** window opens.

**Standard and Interference Assignments**

**Selected Samples**

Un 10 Zn-ReSCN gr2

OK Cancel

Save Element Setup

Save Sample Setup

Add/Remove Standards

Reload Standard Assignments

Remove TDI Correction

1 2 3 4 5 6

**Click Element Row to Edit Standard/Interference/Time Dependent Intensity (TDI) Assignments**

Channel	Element	X-Ray	Analyzed	Standard	Interf-Ele	Interf-Std
1	cs	la	Yes	834	....	0,0,0,0,0
2	fe	ka	Yes	730	....	0,0,0,0,0
3	zn	ka	Yes	660	re,....	575,0,0,0,0
4	re	la	Yes	575	zn,....	660,0,0,0,0
5	s	ka	Yes	730	cs,....	834,0,0,0,0
6	se	ka	Yes	660	....	0,0,0,0,0
7	o		No	0	....	0,0,0,0,0
8	n		No	0	....	0,0,0,0,0
9	c		No	0	....	0,0,0,0,0
10	h		No	0	....	0,0,0,0,0

In John's routine, Cs  $L\alpha$  is measured on the standard CsBr, likewise Fe  $K\alpha$  is done on a pyrite ( $FeS_2$ ) standard and the Se  $K\alpha$  line is measured on the ZnSe standard.

Zinc, rhenium and sulfur require additional discussion.

Highlighting the zinc element row opens the **Assignment Properties** window for that element. The user has the option to display all possible interferences based on the current set of analyzed and interfered elements. Clicking the **Calculate Interference** button displays these. The program calculates the interference based on a gaussian peak shape assuming a worst case



scenario of 0.1% of the analyzed element and 100% of each of the other analyzed and possibly interfering elements.

As mentioned earlier and seen below, Re  $L\alpha$  interferes with the Zn  $K\alpha$  on the LIF crystal. To correct for this overlap, a standard with no zinc present (rhenium metal) is employed for the calculation.

**Assignment Properties**

Enter Standard Assignments for: zn ka

Element	X-Ray	Assigned (Primary) Standard
zn	ka	660 ZnSe (synthetic)

Use Virtual Standard For Standard Intensity Calculation (see Analytical menu)

OK  
Cancel

---

Interference Standard Assignments for Interfered Element: zn ka

	Intf Elem	Intf Order	Interference Standard	
1st	re	1	575 rhenium metal	Remove
2nd				Remove
3rd				Remove
4th				Remove
5th				Remove

Check All Interfering Elements

Calculate Interferences

Minimum Overlap Intensity: .1

The standard used for the interference correction must contain a known concentration of the interfering element and none of the interfered element, nor any other interfering elements.

Interference by Re SLA^Y	at 1.43050 =	539.7%
Interference by Re LA1	at 1.43310 =	62354.2%
Interference by Re SLA3^Z	at 1.43480 =	656.0%
Interference by Re SLAS	at 1.43800 =	657.1%
Interference by Re LA2	at 1.44420 =	5265.9%

---

Time Dependent Intensity (TDI) Calibration Assignment (select unknown sample for calibration)

TDI Correction Type (Self or Assigned)

No TDI Calibration Correction

Use TDI \*Self\* Calibration Correction

Use TDI \*Assigned\* Calibration Correction

Both "assigned" and "self" calibration Time Dependent Intensity (TDI) element samples can be acquired. See the Special Options dialog in the Acquire window.

Both "assigned" and "self" Time Dependent Intensity (TDI) element calibrations can be assigned or unassigned here. Assigned Time Dependent Intensity (TDI) corrections are assigned to samples acquired with the "assigned" flag in Special Options. "Self" Time Dependent Intensity (TDI) corrections are assigned to themselves. Time Dependent Intensity (TDI) element "self" calibrations are automatically assigned to themselves at the time of acquisition.

Display TDI Fit

Use Log-Linear (exponential) Fit

Use Log-Quadratic (hyper-exponential) Fit

---

Blank Correction Sample Assignment

Un 2 Fe-ReSCN
Un 3 Fe-ReSCN gr2
Un 4 Fe-ReSCN gr3
Un 5 Fe-ReSCN gr4
Un 6 Fe-ReSeCN gr1
Un 7 Fe-ReSeCN gr2

Assign a sample to be used for a "blank" trace element correction. The blank sample should be a similar matrix to the unknown sample and should have a zero or known trace of the element present.

Blank Level (wt. %): .000000

Remove Blank Assignment

The **Assignment Properties** window for Re  $L\alpha$  is shown next. Here, Zn  $K\alpha$  interferes with the Re  $L\alpha$  x-ray position. To make the quantitative correction the standard ZnSe (with no rhenium) is used.

**Assignment Properties**

Enter Standard Assignments for: re la

Element	X-Ray	Assigned (Primary) Standard
re	la	575 rhenium metal

Use Virtual Standard For Standard Intensity Calculation (see Analytical menu)

OK Cancel

---

Interference Standard Assignments for Interfered Element: re la

	Intf Elem	Intf Order	Interference Standard	
1st	zn	1	660 ZnSe (synthetic)	Remove
2nd				Remove
3rd				Remove
4th				Remove
5th				Remove

Check All Interfering Elements

Minimum Overlap Intensity: .1

The standard used for the interference correction must contain a known concentration of the interfering element and none of the interfered element, nor any other interfering elements.

Calculate Interferences

Interference by Zn SKA3`	at 1.42920 =	921.2%
Interference by Zn SKA3	at 1.43000 =	950.8%
Interference by Zn SKA`	at 1.43060 =	968.7%
Interference by Zn KA1	at 1.43550 =	96051.6%
Interference by Zn KA2	at 1.43930 =	40269.6%

---

Time Dependent Intensity (TDI) Calibration Assignment (select unknown sample for calibration)

TDI Correction Type (Self or Assigned)

No TDI Calibration Correction  
 Use TDI "Self" Calibration Correction  
 Use TDI "Assigned" Calibration Correction

Both "assigned" and "self" calibration Time Dependent Intensity (TDI) element samples can be acquired. See the Special Options dialog in the Acquire window.

Both "assigned" and "self" Time Dependent Intensity (TDI) element calibrations can be assigned or unassigned here. Assigned Time Dependent Intensity (TDI) corrections are assigned to samples acquired with the "assigned" flag in Special Options. "Self" Time Dependent Intensity (TDI) corrections are assigned to themselves. Time Dependent Intensity (TDI) element "self" calibrations are automatically assigned to themselves at the time of acquisition.

Display TDI Fit

Use Log-Linear (exponential) Fit  
 Use Log-Quadratic (hyper-exponential) Fit

---

Blank Correction Sample Assignment

Un 2 Fe-ReSCN
Un 3 Fe-ReSCN gr2
Un 4 Fe-ReSCN gr3
Un 5 Fe-ReSCN gr4
Un 6 Fe-ReSeCN gr1
Un 7 Fe-ReSeCN gr2

Assign a sample to be used for a "blank" trace element correction. The blank sample should be a similar matrix to the unknown sample and should have a zero or known trace of the element present.

Blank Level (wt. %): .000000

Remove Blank Assignment

Finally, the **Assignment Properties** window for S K $\alpha$  is displayed. In this unique mineral, Cs L $\beta$  overlaps with S K $\alpha$ . Therefore, a cesium standard without sulfur (CsBr) is also required to complete the spectral interference deconvolution!

**Assignment Properties**

Enter Standard Assignments for: s ka

Element	X-Ray	Assigned (Primary) Standard
s	ka	730 pyrite UC # 21334

Use Virtual Standard For Standard Intensity Calculation (see Analytical menu)

**Interference Standard Assignments for Interfered Element: s ka**

	Intf Elem	Intf Order	Interference Standard	
1st	cs	I	834 CsBr (synthetic)	Remove
2nd				Remove
3rd				Remove
4th				Remove
5th				Remove

Check All Interfering Elements

**Calculate Interferences**

Minimum Overlap Intensity: .1

The standard used for the interference correction must contain a known concentration of the interfering element and none of the interfered element, nor any other interfering elements.

For S ka PET at 5.37310 angstroms, concentration of .100 wt. %

Interference by Re LG1	V	at 5.30620 =	11.5%
Interference by Cs LB4	II	at 5.33370 =	878.2%
Interference by Cs LB1	II	at 5.36860 =	22426.5%

**Time Dependent Intensity (TDI) Calibration Assignment (select unknown sample for calibration)**

**TDI Correction Type (Self or Assigned)**

No TDI Calibration Correction  
 Use TDI "Self" Calibration Correction  
 Use TDI "Assigned" Calibration Correction

Both "assigned" and "self" calibration Time Dependent Intensity (TDI) element samples can be acquired. See the Special Options dialog in the Acquire window.

Both "assigned" and "self" Time Dependent Intensity (TDI) element calibrations can be assigned or unassigned here. Assigned Time Dependent Intensity (TDI) corrections are assigned to samples acquired with the "assigned" flag in Special Options. "Self" Time Dependent Intensity (TDI) corrections are assigned to themselves. Time Dependent Intensity (TDI) element "self" calibrations are automatically assigned to themselves at the time of acquisition.

**Display TDI Fit**

Use Log-Linear (exponential) Fit  
 Use Log-Quadratic (hyper-exponential) Fit

**Blank Correction Sample Assignment**

Un 2 Fe-ReSCN
Un 3 Fe-ReSCN gr2
Un 4 Fe-ReSCN gr3
Un 5 Fe-ReSCN gr4
Un 6 Fe-ReSeCN gr1
Un 7 Fe-ReSeCN gr2

Assign a sample to be used for a "blank" trace element correction. The blank sample should be a similar matrix to the unknown sample and should have a zero or known trace of the element present.

Blank Level (wt. %): .000000

**Remove Blank Assignment**

After setting up the parameters for the analysis session, calibration and standardization was accomplished, then several tiny and complex organometallic phases were probed for their chemistry.

The **Analysis Calculation Options** window allows the user to enable or disable the interference correction routine to view its effect on quantitative analysis results. This window is activated from the **Analytical | Analysis Options** menu in the main PROBE FOR EPMA log window.

Analysis Calculation Options

**Quantitative Acquisition Options**

- Use Deadtime Correction
  - Use Normal Deadtime Correction (single term factorial)
  - Use Precision Deadtime Correction (two term factorial for > 50K cps)
- Use Beam Drift Correction
- Use Automatic Drift Correction on Standard Intensities

**Quantitative Analysis Options**

- Use Assigned Interference Corrections on Standards and Unknowns
  - Do Not Use Matrix Correction Term For Full Quant Interference Correction
  - Do Not Use Full Quant Interference Corrections (use traditional Gilfrich, et. al.)
- Use Assigned or Self Time Dependent Intensity (TDI) Corrections on Unknowns
  - Use Linear Fit (slope coefficient only) for TDI Extrapolation
  - Use Quadratic Fit (two coefficient) for TDI Extrapolation
  - Use Time Weighted Data for TDI Fit (weight intensities based on elapsed time)
    - Time Weighted Data Weight Factor:
- Use Absorption Corrected MAN Continuum Intensities
- Use Particle or Thin Film Correction Parameters
- Check For Same Peak Positions in Unknown and Standard
- Check For Same PHA Settings in Unknown and Standard
- Use Zero Point For Calibration Curve (off-peak elements only)
- Use Conductive Coating Correction For Beam Energy Loss
- Use Conductive Coating Correction For X-ray Absorption

**MAC (mass absorption coefficient) and APF (area peak factor) Options**

- Use Empirical MAC Values
- Use Empirical APF Values
- Use Empirical APF Factors (calculated from elemental composition)
- Use Specified APF Factors (based on a fixed composition)

**Calculation Options**

- Use Aggregate Intensities for Duplicate Quantitative Elements
- Use Blank Calibration Sample Trace Element Accuracy Corrections
- Force Negative K-Ratios To Zero in ZAF Calculations
- Calculate Electron and Xray Ranges for Sample Compositions
- Use Oxygen From Halogens (F, Cl, Br and I) Correction
- Use Nth Point Calculation For Off-Peak Intensities (for testing only)
- Use Count Overwrite Intensity Table for Data Calculations
- Force Negative Interference Intensities To Zero in Corrections
- Use Chemical Age Calculation (U, Th, Pb)

**Formatting Options**

- Use Automatic Format For Quantitative Results
  - Display the Maximum Number of Numerical Digits
  - Display Only Statistically Significant Number of Numerical Digits
- Use Detailed Printout For Data and Analytical Results
- Print Analyzed And Specified On Same Line
- Display Count Intensities Unnormalized To Time (in Analyze!)

**Output Options**

- Display Charge Balance Calculation

Elemental Output Sort Order For JJD-2 and HW Custom Output:

- Use Traditional Geological Sort Order (SiO2, TiO2, etc.)
- Use Low To High Atomic Number Sort Order
- Use High To Low Atomic Number Sort Order

OK

Cancel

Data illustrating the power of the PROBE FOR EPMA interference correction routine is shown below. The results are for a real unknown sample that has some variability and was rather small to analyze. The following unknown analysis is without any interference corrections.

```
Un 10 Zn-ReSCN gr2
TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0
(Magnification (analytical) = 2000), Beam Mode = Analog Spot
(Magnification (default) = 0, Magnification (imaging) = 40)
Image Shift (X,Y): -2, 3
Number of Data Lines: 5 Number of 'Good' Data Lines: 3
First/Last Date-Time: 06/15/1998 04:39:29 PM to 06/15/1998 04:55:23 PM
WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 122.912
Average Calculated Oxygen: .000 Average Atomic Number: 53.080
Average Excess Oxygen: .000 Average Atomic Weight: 54.153
Average ZAF Iteration: 4.00 Average Quant Iterate: 2.00
```

Un 10 Zn-ReSCN gr2, Results in Elemental Weight Percents

SPEC:	O	N	C	H
TYPE:	SPEC	SPEC	SPEC	SPEC
AVER:	1.900	5.000	4.200	.200
SDEV:	.000	.000	.000	.000

ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
53	.000	.000	19.553	73.669	17.465	.000	121.988
55	.000	.007	20.556	74.517	16.507	.000	122.886
56	.000	.019	19.668	74.715	18.159	.000	123.862
AVER:	.000	.009	19.926	74.301	17.377	.000	122.912
SDEV:	.000	.010	.549	.556	.830	.000	.937
SERR:	.000	.006	.317	.321	.479	.000	
%RSD:	.10	113.30	2.75	.75	4.77	.13	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	627.07	3670.93	3712.94	3672.80	4423.12	1315.94	
UNKF:	.0000	.0001	.2100	.6516	.0943	.0000	
UNCT:	-2.69	.33	1559.22	2393.02	884.28	-3.16	
UNBG:	11.62	28.03	67.30	61.54	6.57	64.00	
ZCOR:	1.2016	.9276	.9487	1.1404	1.8421	1.0254	
KRAW:	-.0043	.0001	.4199	.6516	.1999	-.0024	
PKBG:	.77	1.01	24.18	40.00	135.89	.95	

Un 10 Zn-ReSCN gr2, Results Based on 6 Atoms of re

SPEC:	O	N	C	H			
TYPE:	SPEC	SPEC	SPEC	SPEC			
AVER:	1.786	5.368	5.258	2.983			
SDEV:	.013	.040	.039	.022			
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
53	.000	.000	4.536	6.000	8.260	.000	34.323
55	.000	.002	4.714	6.000	7.718	.000	33.784
56	.000	.005	4.499	6.000	8.468	.000	34.281
AVER:	.000	.002	4.583	6.000	8.149	.000	34.129
SDEV:	.000	.003	.115	.000	.387	.000	.300
SERR:	.000	.002	.067	.000	.224	.000	
%RSD:	.77	113.20	2.51	.00	4.75	.79	
%RSD:	.8	113.2	2.5	.0	4.8	.8	

The user immediately realizes that there is a problem since the unknown sums to a total of 123%.

The following is the same unknown but with the iterated interference software applied.

Un 10 Zn-ReSCN gr2  
 TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 5 Number of 'Good' Data Lines: 3  
 First/Last Date-Time: 06/15/1998 04:39:29 PM to 06/15/1998 04:55:23 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 101.367  
 Average Calculated Oxygen: .000 Average Atomic Number: 54.332  
 Average Excess Oxygen: .000 Average Atomic Weight: 50.124  
 Average ZAF Iteration: 4.00 Average Quant Iterate: 13.00

Un 10 Zn-ReSCN gr2, Results in Elemental Weight Percents

SPEC:	O	N	C	H			
TYPE:	SPEC	SPEC	SPEC	SPEC			
AVER:	1.900	5.000	4.200	.200			
SDEV:	.000	.000	.000	.000			
ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
53	.000	.000	6.237	65.850	17.508	.000	100.895
55	.000	.007	7.365	65.342	16.513	.000	100.527
56	.000	.020	6.100	67.048	18.210	.000	102.678
AVER:	.000	.009	6.568	66.080	17.410	.000	101.367
SDEV:	.000	.010	.694	.876	.853	.000	1.150
SERR:	.000	.006	.401	.506	.492	.000	
%RSD:	.13	113.29	10.57	1.33	4.90	.03	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	627.07	3670.93	3712.94	3672.80	4423.12	1315.94	
UNKF:	.0000	.0001	.0691	.5781	.0943	.0000	
UNCT:	-2.69	.33	512.67	2123.23	884.28	-3.16	
UNBG:	11.62	28.03	67.30	61.54	6.57	64.00	
ZCOR:	1.2162	.9395	.9511	1.1431	1.8455	1.0235	
KRAW:	-.0043	.0001	.1381	.5781	.1999	-.0024	
PKBG:	.77	1.01	8.61	35.59	135.89	.95	
INT%:	----	----	-67.17	-11.27	----	----	

Un 10 Zn-ReSCN gr2, Results Based on 6 Atoms of re

SPEC:	O	N	C	H			
TYPE:	SPEC	SPEC	SPEC	SPEC			
AVER:	2.008	6.036	5.913	3.355			
SDEV:	.026	.080	.078	.044			
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
53	.000	.000	1.619	6.000	9.264	.000	34.252
55	.000	.002	1.926	6.000	8.805	.000	34.239
56	.000	.006	1.555	6.000	9.463	.000	34.083
AVER:	.000	.003	1.700	6.000	9.177	.000	34.191
SDEV:	.000	.003	.199	.000	.337	.000	.094
SERR:	.000	.002	.115	.000	.195	.000	
%RSD:	1.21	112.34	11.69	.00	3.68	1.30	

The total now is acceptable, slightly over 100%.

Finally to demonstrate with standard samples (of known composition), both rhenium metal and the ZnSe standard will be reproduced without the benefit of the interference correction.

St 575 Set 1 rhenium metal  
 TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 5 Number of 'Good' Data Lines: 5  
 First/Last Date-Time: 06/15/1998 01:40:14 PM to 06/15/1998 01:48:18 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 124.489  
 Average Calculated Oxygen: .000 Average Atomic Number: 67.420  
 Average Excess Oxygen: .000 Average Atomic Weight: 142.028  
 Average ZAF Iteration: 3.00 Average Quant Iterate: 2.00

St 575 Set 1 rhenium metal, Results in Elemental Weight Percents

ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
1	.000	.000	21.083	103.240	.022	.000	124.344
2	.161	.005	20.061	102.988	.008	.000	123.223
3	.000	.000	21.644	103.274	.000	.176	125.093
4	.000	.023	21.565	103.834	.000	.156	125.578
5	.000	.000	19.873	104.122	.000	.213	124.207
AVER:	.032	.006	20.845	103.492	.006	.109	124.489
SDEV:	.072	.010	.833	.469	.009	.102	.902
SERR:	.032	.004	.372	.210	.004	.045	
%RSD:	223.61	177.20	3.99	.45	159.05	93.18	
PUBL:	n.a.	n.a.	n.a.	100.000	n.a.	n.a.	100.000
%VAR:	.00	.00	.00	3.49	.00	.00	
DIFF:	.000	.000	.000	3.492	.000	.000	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	628.01	3665.44	3726.73	3679.31	4367.49	1316.18	
UNKF:	.0003	.0001	.2405	1.0000	.0000	.0012	
UNCT:	-.02	-1.11	1792.16	3679.35	-.42	1.76	
UNBG:	11.57	41.92	95.49	93.54	8.43	93.41	
ZCOR:	1.2008	.8713	.8667	1.0349	2.1623	.9212	
KRAW:	.0000	-.0003	.4809	1.0000	-.0001	.0013	
PKBG:	1.00	.97	19.77	40.37	.95	1.02	

St 660 Set 1 ZnSe (synthetic)  
 TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 5 Number of 'Good' Data Lines: 5  
 First/Last Date-Time: 06/15/1998 01:50:41 PM to 06/15/1998 01:58:46 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 151.911  
 Average Calculated Oxygen: .000 Average Atomic Number: 47.924  
 Average Excess Oxygen: .000 Average Atomic Weight: 92.977  
 Average ZAF Iteration: 3.00 Average Quant Iterate: 2.00

St 660 Set 1 ZnSe (synthetic), Results in Elemental Weight Percents

ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	LIN
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
6	.000	.000	44.441	55.414	.019	51.417	151.290
7	.000	.014	44.022	56.029	.036	51.668	151.768
8	.000	.000	44.551	56.249	.035	51.099	151.934
9	.000	.015	44.701	56.069	.037	50.975	151.796
10	.000	.000	45.281	56.025	.052	51.409	152.767
AVER:	.000	.006	44.599	55.957	.036	51.314	151.911
SDEV:	.000	.008	.457	.317	.012	.277	.537
SERR:	.000	.003	.204	.142	.005	.124	
%RSD:	.02	137.06	1.02	.57	33.50	.54	
PUBL:	n.a.	n.a.	45.290	n.a.	n.a.	54.710	100.000
%VAR:	.00	.00	-1.53	.00	.00	-6.21	
DIFF:	.000	.000	-.691	.000	.000	-3.396	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	628.01	3665.44	3726.64	3678.98	4367.49	1316.18	
UNKF:	.0000	.0001	.5002	.5193	.0002	.5162	
UNCT:	-23.07	-1.38	3726.73	1910.63	1.72	1316.18	
UNBG:	29.32	22.96	48.34	47.66	5.61	47.46	
ZCOR:	1.1332	.8826	.8917	1.0775	1.9149	.9940	
KRAW:	-.0367	-.0004	1.0000	.5193	.0004	1.0000	
PKBG:	.21	.95	78.18	41.17	1.31	28.85	

The rhenium standard displays an apparent rhenium concentration of 103% and a zinc total of an additional 20%. The ZnSe is even more interesting in that the total is approaching 152% with 44% zinc, 56% rhenium and 51% selenium. The true composition is 45% zinc, and 55% selenium with no rhenium! The normal matrix correction comes close with respect to both zinc and selenium abundances but also reports a whopping 56% rhenium concentration.



Below, both standards are rerun with the interference corrections applied.

St 575 Set 1 rhenium metal  
 TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 5 Number of 'Good' Data Lines: 5  
 First/Last Date-Time: 06/15/1998 01:40:14 PM to 06/15/1998 01:48:18 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 100.138  
 Average Calculated Oxygen: .000 Average Atomic Number: 74.736  
 Average Excess Oxygen: .000 Average Atomic Weight: 184.273  
 Average ZAF Iteration: 2.00 Average Quant Iterate: 13.80

St 575 Set 1 rhenium metal, Results in Elemental Weight Percents

ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
1	.000	.000	.375	99.332	.023	.000	99.730
2	.164	.005	.000	99.625	.007	.000	99.801
3	.000	.000	1.096	98.532	.000	.168	99.796
4	.000	.023	.847	99.364	.000	.149	100.383
5	.000	.000	.000	100.776	.000	.203	100.979
AVER:	.033	.006	.463	99.526	.006	.104	100.138
SDEV:	.074	.010	.496	.810	.010	.097	.539
SERR:	.033	.004	.222	.362	.004	.043	
%RSD:	223.61	177.20	107.02	.81	164.63	93.19	
PUBL:	n.a.	n.a.	n.a.	100.000	n.a.	n.a.	100.000
%VAR:	.00	.00	.00	-.47	.00	.00	
DIFF:	.000	.000	.000	-.474	.000	.000	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	628.01	3665.44	3726.73	3679.31	4367.49	1316.18	
UNKF:	.0003	.0001	.0055	.9942	.0000	.0012	
UNCT:	-.02	-1.11	10.26	3657.79	-.43	1.76	
UNBG:	11.57	41.92	95.49	93.54	8.43	93.41	
ZCOR:	1.2302	.8692	.8409	1.0011	2.3003	.8800	
KRAW:	.0000	-.0003	.0028	.9942	-.0001	.0013	
PKBG:	1.00	.97	1.11	40.14	.95	1.02	
INT%:	----	----	-99.57	-.59	-2.86	----	

St 660 Set 1 ZnSe (synthetic)  
 TakeOff = 40.0 KiloVolt = 20.0 Beam Current = 20.0 Beam Size = 0  
 (Magnification (analytical) = 2000), Beam Mode = Analog Spot  
 (Magnification (default) = 0, Magnification (imaging) = 40)  
 Image Shift (X,Y): -2, 3  
 Number of Data Lines: 5 Number of 'Good' Data Lines: 5  
 First/Last Date-Time: 06/15/1998 01:50:41 PM to 06/15/1998 01:58:46 PM  
 WARNING- Forcing negative k-ratios to zero

Average Total Oxygen: .000 Average Total Weight%: 100.261  
 Average Calculated Oxygen: .000 Average Atomic Number: 32.316  
 Average Excess Oxygen: .000 Average Atomic Weight: 72.276  
 Average ZAF Iteration: 3.00 Average Quant Iterate: 15.00

St 660 Set 1 ZnSe (synthetic), Results in Elemental Weight Percents

ELEM:	Cs	Fe	Zn	Re	S	Se	
BGDS:	LIN	LIN	LIN	LIN	LIN	LIN	LIN
TIME:	10.00	10.00	10.00	10.00	10.00	10.00	
BEAM:	.20	.20	.20	.20	.20	.20	
ELEM:	Cs	Fe	Zn	Re	S	Se	SUM
6	.000	.000	45.121	.000	.017	54.799	99.937
7	.000	.014	44.451	1.080	.032	54.981	100.558
8	.000	.000	45.153	.477	.032	54.459	100.120
9	.000	.015	45.415	.000	.033	54.368	99.831
10	.000	.000	46.003	.000	.047	54.807	100.857
AVER:	.000	.006	45.229	.311	.032	54.683	100.261
SDEV:	.000	.008	.560	.477	.011	.258	.434
SERR:	.000	.003	.251	.213	.005	.116	
%RSD:	.05	137.06	1.24	153.10	33.46	.47	
PUBL:	n.a.	n.a.	45.290	n.a.	n.a.	54.710	100.000
%VAR:	.00	.00	-.14	.00	.00	-.05	
DIFF:	.000	.000	-.061	.000	.000	-.027	
STDS:	834	730	660	575	730	660	
STKF:	.5978	.4297	.5001	1.0000	.4719	.5162	
STCT:	628.01	3665.44	3726.64	3678.98	4367.49	1316.18	
UNKF:	.0000	.0001	.4995	.0028	.0002	.5162	
UNCT:	-23.07	-1.38	3721.55	2.34	1.72	1316.18	
UNBG:	29.32	22.96	48.34	47.66	5.61	47.46	
ZCOR:	1.0893	.8842	.9055	1.1054	1.7229	1.0593	
KRAW:	-.0367	-.0004	.9986	.0006	.0004	1.0000	
PKBG:	.21	.95	78.08	1.05	1.31	28.85	
INT%:	----	----	-.14	-99.88	----	----	

Now, the apparent zinc in the rhenium metal standard is gone and the average total sums correctly near 100%. The ZnSe standard is perfect, matching the published standard composition for both zinc and selenium, and effectively removing all of the apparent 56% rhenium.

## Light Element Analysis - Empirical APFs

Quantitative analysis of light elements (beryllium to fluorine) is difficult to do correctly with the electron microprobe. Numerous issues impede the analysis of light elements (see Appendix B in the User's Guide and Reference documentation as well as Goldstein et al., 1992 for further discussion). A few brief comments will be made here, as an introduction to this section.

Typically, for x-ray analysis in this energy range (0.1 to 0.7 keV), a large absorption correction is necessary. This large correction in conjunction with the fact that the mass absorption coefficients for the low energy x-rays are very large and not very well known (see Appendix C of the User's Guide and Reference documentation) place a severe demand on the established ZAF and Phi-Rho-Z models for light element matrix corrections. Some reduction in this large x-ray absorption factor is possible by analyzing at higher take-off angles and at lower acceleration voltages. The former aids by shortening the path length for absorption in the sample while the latter involves a decrease in the electron beam penetration so that x-rays are generated closer to the surface and can escape to be detected.

Low count rates for these light elements are often found. This is due to the low fluorescent yields from the  $K\alpha$  x-ray lines and the inefficient nature of WDS counting systems. Count rates can be increased by increasing the beam current substantially but this may then lead to deadtime problems for metal lines that interfere with the light element lines of interest. The use of new layered dispersive element (LDE) synthetic multi-layer crystals with large d-spacings can also improve light element peak count rates.

There is also a strong possibility of interferences from higher order metal lines such as titanium, chromium, manganese, iron, nickel, zirconium, niobium, and molybdenum with the light element lines. These interferences are often severe for minor and trace level measurements. It is critical to eliminate peak overlaps and understand the background intensity around the light element peak position in both the sample and standard. The new LDE multi-layer crystals also help here by strongly suppressing these higher order reflections.

Finally, chemical bonding effects can result in wavelength shifts, increases or decreases in the relative intensities of various lines and alteration of the shape of the analytical line. These effects are most significant for the lightest (lowest energy) elements. Polarization phenomena and crystallographic orientation may also cause variations in peak shape and intensity especially for boron. Therefore, to measure the intensity of the light elements, one measures the integral intensity (area) under the characteristic peak rather than just the peak intensity.

Bastin and Heijliers (1984, 1991) pioneered the area-peak factor (APF) analysis method. The APF is defined as the ratio between the integral intensity (area) k-ratio from the sample and the standard and the peak intensity k-ratio from the same sample and the standard. This factor is only valid for a given compound with respect to a given standard and a specific spectrometer setup. These integral measurements can be made by scanning the spectrometer in small steps across the characteristic peak with the wavelscan feature in PROBE FOR EPMA. After acquiring the peak shape profiles for a primary standard and a number of secondary standards and correcting for the background and removal of interfering peaks, the APF can be calculated as:

$$\text{APF} = \frac{I_U^I \cdot I_S^P}{I_U^P \cdot I_S^I}$$

Where :

$I_U^I$	is the integrated intensity of the secondary standard
$I_S^P$	is the peak intensity of the primary standard
$I_U^P$	is the peak intensity of the secondary standard
$I_S^I$	is the integrated intensity of the primary standard

After an APF has been determined, future measurements of that compound can be accomplished by measuring only the peak intensity in the sample and standard. Then, multiplication of the peak k-ratio with the appropriate APF will yield the correct integral k-ratio.

PROBE FOR EPMA allows the user to select an APF correction for use in correcting x-ray intensities for peak shift and shape changes between the standard and the multi-element unknown. This is critical when the user is analyzing the  $K\alpha$  lines of the light elements (boron, carbon, nitrogen and oxygen). This correction may also be of use for other elements such as S  $K\alpha$  that also exhibits peak shift and shape changes when comparing sulfate and sulfide peaks.

The power of this correction will be documented in the analysis of oxygen using several simple oxide standards. Open a new PROBE FOR EPMA run. Here oxygen  $K\alpha$  will be analyzed using MgO as the calibration standard. The spectrometer is equipped with a standard P-10 gas flow counter and a LDE (59.8Å) W/Si reflecting crystal. The other two elements to be determined are magnesium on MgO and silicon on  $\text{SiO}_2$ . Peak the three elements and acquire standard samples for each.

Create a new unknown sample and collect data on the SiO<sub>2</sub> standard. Analyze the sample from the **Analyze!** window.

The screenshot shows the 'Analyze!' software window. The 'Sample List' on the left contains two entries: 'Un 1 \* setup' and 'Un 2 SiO2 std check'. The 'Standard Assignments' tab is active, showing a table of results for 'Un 2 SiO2 std check'. The table includes columns for 'Copy', 'Si', 'Mg', 'O', and 'Total'. The 'Average' row shows Si: 46.997, Mg: .005, O: 50.987, and Total: 97.989. Below this, there are rows for 'Std Dev', 'ZAF Corr', 'Std Err', '%Rel SD', 'Minimum', and 'Maximum'. At the bottom, there are three buttons: 'Delete Selected Line(s)', 'Undelete Selected Line(s)', and 'Analyze Selected Line(s)'. The 'Analyze Selected Line(s)' button is highlighted, and the table below it shows three rows of data for '21 G', '22 G', and '23 G'.

Copy	Si	Mg	O	Total
Average:	46.997	.005	50.987	97.989
Std Dev:	.311	.009	.302	.533
ZAF Corr:	1.1318	1.3928	2.0346	
Std Err:	.179	.005	.174	.308
%Rel SD:	.66	173.20	.59	.54
Minimum:	46.715	.000	50.663	97.378
Maximum:	47.330	.016	51.261	98.365

Copy	Si	Mg	O	Total
21 G	46.715	.000	50.663	97.378
22 G	46.946	.016	51.261	98.223
23 G	47.330	.000	51.035	98.365

A low total for the analysis is found. The nominal composition for the SiO<sub>2</sub> standard is silicon: 46.74 and oxygen: 53.26. Here, the discrepancy in the total rests primarily with the oxygen concentration. The low oxygen concentration is independent of the matrix correction (and mass absorption coefficient) and can only be corrected for by using the appropriate APFs.

Select **Analytical | Analysis Options** from the main PROBE FOR EPMA log window to open the **Analysis Calculation Options** window. Click the *Use Empirical APF Values* check box to activate this option.

**Analysis Calculation Options**

**Quantitative Acquisition Options**

- Use Deadtime Correction
  - Use Normal Deadtime Correction (single term factorial)
  - Use Precision Deadtime Correction (two term factorial for > 50K cps)
- Use Beam Drift Correction
- Use Automatic Drift Correction on Standard Intensities

**Quantitative Analysis Options**

- Use Assigned Interference Corrections on Standards and Unknowns
  - Do Not Use Matrix Correction Term For Full Quant Interference Correction
  - Do Not Use Full Quant Interference Corrections (use traditional Gilfrich, et. al.)
- Use Assigned or Self Time Dependent Intensity (TDI) Corrections on Unknowns
  - Use Linear Fit (slope coefficient only) for TDI Extrapolation
  - Use Quadratic Fit (two coefficient) for TDI Extrapolation
- Use Time Weighted Data for TDI Fit (weight intensities based on elapsed time)
  - Time Weighted Data Weight Factor:
- Use Absorption Corrected MAN Continuum Intensities
- Use Particle or Thin Film Correction Parameters
- Check For Same Peak Positions in Unknown and Standard
- Check For Same PHA Settings in Unknown and Standard
- Use Zero Point For Calibration Curve (off-peak elements only)
- Use Conductive Coating Correction For Beam Energy Loss
- Use Conductive Coating Correction For X-ray Absorption

**MAC (mass absorption coefficient) and APF (area peak factor) Options**

- Use Empirical MAC Values
- Use Empirical APF Values
  - Use Empirical APF Factors (calculated from elemental composition)
  - Use Specified APF Factors (based on a fixed composition)

**Calculation Options**

- Use Aggregate Intensities for Duplicate Quantitative Elements
- Use Blank Calibration Sample Trace Element Accuracy Corrections
- Force Negative K-Ratios To Zero in ZAF Calculations
- Calculate Electron and Xray Ranges for Sample Compositions
- Use Oxygen From Halogens (F, Cl, Br and I) Correction
- Use Nth Point Calculation For Off-Peak Intensities (for testing only)
- Use Count Overwrite Intensity Table for Data Calculations
- Force Negative Interference Intensities To Zero in Corrections
- Use Chemical Age Calculation (U, Th, Pb)

**Formatting Options**

- Use Automatic Format For Quantitative Results
  - Display the Maximum Number of Numerical Digits
  - Display Only Statistically Significant Number of Numerical Digits
- Use Detailed Printout For Data and Analytical Results
- Print Analyzed And Specified On Same Line
- Display Count Intensities Unnormalized To Time (in Analyze!)

**Output Options**

- Display Charge Balance Calculation

Elemental Output Sort Order For JJD-2 and HW Custom Output:

- Use Traditional Geological Sort Order (SiO<sub>2</sub>, TiO<sub>2</sub>, etc.)
- Use Low To High Atomic Number Sort Order
- Use High To Low Atomic Number Sort Order

**OK** **Cancel**

Click the **OK** button.

Next, select **Analytical | Empirical APFs** from the main PROBE FOR EPMA log window to open the **Add Empirical APFs (area peak factors) to Run** window.

Find the appropriate correction. Two choices are given for oxygen in the presence of SiO<sub>2</sub>, Bastin's value of 1.04440 using Fe<sub>2</sub>O<sub>3</sub> as the calibration standard and Donovan's correction factor of 1.070 when using MgO as the calibration standard. Although the values seen in the text field below are relative to Fe<sub>2</sub>O<sub>3</sub>, the APF for MgO relative to Fe<sub>2</sub>O<sub>3</sub> is also 1.000, which means that these correction factors apply equally well relative to MgO. Hence, the use of MgO as the calibration standard.

The APF correction values are defined in the EMPAPF.DAT file in the PROBE FOR EPMA directory. The file contains some 80 values that may or may not be applicable, depending on the analyzing crystals and standards available. It is strongly recommended that the user measure the integral intensities and peak intensities and calculate the APFs for your particular spectrometer setup. And if necessary, edit the EMPAPF.DAT file using any ASCII text editor such as NotePad to insert their own measurements. See the User's Guide and Reference documentation for editing format and details.

Click the **Add To Run >>** button to place the selected empirical APF into the run.

**Add Empirical APFs (area peak factors) to Run**

Available Empirical APFs from C:\Program Files\Probe Software\Probe for EPMA\EMPAPF.DAT

Element	APF	Standard
o ka in b	1.06280	B60/Fe203/Wsi/59.8
o ka in na	1.00000	----/Fe203/Wsi/59.8
o ka in mg	1.00000	Mg0/Fe203/Wsi/59.8
o ka in al	1.02130	Al203/Fe203/Wsi/59.8, Bastin
o ka in al	1.02850	Al203/Mg0/Wsi/59.8, Donovan
o ka in si	1.04440	Si02/Fe203/Wsi/59.8, Bastin
o ka in si	1.07000	Si02/Mg0/Wsi/59.8, Donovan
o ka in p	1.05000	----/Fe203/Wsi/59.8
o ka in s	1.20000	----/Fe203/Wsi/59.8
o ka in k	1.00000	----/Fe203/Wsi/59.8
o ka in ca	.970000	----/Fe203/Wsi/59.8
o ka in ti	.979600	Ti02/Fe203/Wsi/59.8

Current Empirical APFs in Run

o ka in si	1.07000	Si02/Mg0/Wsi/59.8, Donovan
------------	---------	----------------------------

The Area Peak Factor (APF) is the ratio of the integrated to peak intensity of the unknown sample relative to the primary standard for that element. This dialog allows the user to select binary (two components, e.g., O Ka in Si) APFs for the analyzed elements to correct for peak shape changes in unknowns of arbitrary compositions containing those elements. For example, if one desired to correct for peak shape changes for measuring oxygen in a variety of samples containing different amounts of Si, Al and Mg, one would load individual APFs for O Ka in SiO<sub>2</sub>, in Al<sub>2</sub>O<sub>3</sub> and in MgO. For situations where a limited range of composition is being measured and binary component APFs are not available, use the "Specified APFs" in the Elements/Cations dialog (Analyze! window) where a single APF can be specified.

**Add To Run >>**

**Delete From Run**

$$APF = \frac{I_U^I \cdot I_S^P}{I_U^P \cdot I_S^I}$$

**Re-Normalization Factor**

**Re-Normalization Standard**

If using a different primary standard for calculating empirical APF factors than the reference standard indicated above (usually B, Fe3C, AlN or MgO), enter the APF for the primary standard here before clicking the "Add To Run" button. The Re-Normalization standard is a string indicating the actual standard used for the primary standard and will be appended to the comment string for documentation purposes.

**OK** **Cancel**

Click the **OK** button.

Re-analyze the unknown SiO<sub>2</sub> sample. The total is now closer to 100% and a more reasonable oxygen concentration is calculated.

**Analyze!**

Sample List (multi-select) (double-click to see intensity data)

- Standards
- Unknowns
- Wavescans
- All Samples

Select All

Add To Setup

Save Setups

Un 1 \* setup

Un 2 SiO2 std check

Analyze

Data

KRaws

Combine Selected Samples

>>Excel

List Report

Calculation Options

Combine Analysis Lines From Selected Samples

Pause Between Samples

Use All Matrix Corrections

Report

Delete Selected Sample(s)

Undelete Selected Sample(s)

Match

Combine Data Lines From Selected Samples

Combined Conditions

Count Times

Sort Stat and Data Grids In Geological or Atomic Number Order

Specified Concentrations

Standard Assignments

Name/Description

Conditions

Elements/Cations

Un 2 SiO2 std check

TO = 40, KeV = 15, Beam = 30, Size = 2

(MagAnal = 2000.), Mode = Analq Spot

Results in Elemental Weight Percent

	.000	Total Oxygen	99.669	Total Weight %
	.000	Calculated Oxygen	10.834	Z - Bar
	.000	Excess Oxygen	20.082	Atomic Weight

Copy	Si	Mg	O	Total
Average:	47.078	.005	52.586	99.669
Std Dev:	.311	.009	.308	.543
ZAF Corr:	1.1337	1.3980	2.0170	
Std Err:	.180	.005	.178	.313
%Rel SD:	.66	173.20	.59	.54
Minimum:	46.796	.000	52.253	99.048
Maximum:	47.412	.016	52.859	100.057

Delete Selected Line(s)

Undelete Selected Line(s)

Analyze Selected Line(s)

Copy	Si	Mg	O	Total
21 G	46.796	.000	52.253	99.048
22 G	47.026	.016	52.859	99.901
23 G	47.412	.000	52.645	100.057

Cancel

Next



The APF values selected or entered are always measured relative to some standard sample. For example, if measuring oxygen  $K\alpha$  and using either MgO or  $Fe_2O_3$  as the primary standard for oxygen, then any oxygen  $K\alpha$  APF values used must be those measured relative to either MgO or  $Fe_2O_3$ . For the same reason, if using APF values for a particular light element and the user decides to re-assign the standard for that element, the APF values for that element must also be re-normalized to reflect the standard re-assignment.

Thus, it is most efficient to always use the same standard for each light element analyzed. Typically (in order to utilize the APF values in the supplied EMPAPF.DAT file) these will be:

- Oxygen : MgO or  $Fe_2O_3$
- Nitrogen : AlN
- Carbon :  $Fe_3C$
- Boron : B metal

The APF correction in PROBE FOR EPMA will allow the user to enter empirical APF values in each run. The user may enter one or more APF factors for each emitting element although they are generally applied to soft x-ray lines. The APF for each absorber will be summed according to its weight fraction in the composition and applied to the emitting element counts during the ZAF or Phi-Rho-Z iteration.

(This page intentionally left blank)

## References

- Bastin G.F., and Heijliers H.J.M., (1984) Quantitative Electron-Probe Microanalysis of Carbon in Binary Carbides., *Microbeam Analysis* (Editors: A.D. Romig, Jr. and J.I. Goldstein) pp 291-294.
- Bastin G.F., and Heijliers H.J.M., (1991) Quantitative Electron-Probe Microanalysis of Ultra-Light Elements (Boron-Oxygen)., *Electron Probe Quantitation* (Editors: K.F.J. Heinrich and D.E. Newbury), New York: Plenum Press pp 145-161.
- Deer W.A., Howie R.A., and Zussman J., (1992) *An introduction to The Rock-Forming Minerals, Second Edition*, New York: Wiley.
- Donovan J.J., Snyder D.A., and Rivers M.L., (1993) An Improved Interference Correction for Trace Element Analysis. *Microbeam Analysis* **2** 23-28.
- Goldstein J.I., Newbury D.E., Echlin P., Joy D.C., Romig A.D. Jr, Lyman C.E., Fiori C., and Lifshin E., (1992) *Scanning Electron Microscopy and X-Ray Microanalysis, Second Edition*, New York: Plenum.
- Nielsen C.H., and Sigurdsson H., (1981) Quantitative methods for electron microprobe analysis of sodium in natural and synthetic glasses. *Am. Mineral.*, **66** 547-552.
- Reed S.J.B., (1993) *Electron-Microprobe Analysis, Second Edition*, New York: Wiley.
- Scott V.D., Love G., and Reed S.J.B. (1995) *Quantitative Electron-Probe Microanalysis, Second Edition*, London: Ellis Horwood.
- Willis J.P., (1993) *Course on Theory and Practice of XRF Spectrometry*, University of Western Ontario, Department of Geology, Course Notes.