



Release Notes

DIFFRAC.SUITE TOPAS V.5.0

This file contains important supplementary information that may not appear in the main product documentation. Bruker AXS recommends reading this file in its entirety.

Contents

PROGRAM VERSION

- Date and Version Numbers

- Changes to Previous Version(s)

- Installation

SYSTEM REQUIREMENTS

- Operating system

- Minimum hardware configuration

- Recommended hardware configuration

CONTACT AND PROBLEM REPORTS

PROGRAM VERSION

Date and Version Numbers

TOPAS V.5.0, dated 20.11.2014

The installation file has the version number 5.0.0.13.

Changes to Previous Version(s)

See "TOPAS 5 Whats New.pdf" in the Documents folder.

Installation

1. Log on to Microsoft Windows as an administrator.
2. Start DIFFRACTOPAS.exe.
3. Click "Next".
4. Change the installation folder if required and click "Next".
5. Select the "COD database" feature if the structure files from the Crystallography Open Database should be installed and click "Next".
6. Click "Install".
7. Click "Finish".

SYSTEM REQUIREMENTS

Operating system

- Windows XP SP3, Windows 7 32bit and 64bit, Windows 8 or 8.1 32 bit and 64bit

Minimum hardware configuration

TOPAS can be installed with the minimum configuration, but performance will be lower than expected.

- 2.0 GHz or higher CPU
- 1 GB of RAM
- 0.5 GB of free hard disk space
- 800 x 600 pixel screen size

Recommended hardware configuration

- 3.0 GHz or higher multi-core CPU
- 4 GB RAM
- 1 GB free hard disk space
- 1920 x 1080 pixel screen size

CONTACT AND PROBLEM REPORTS

Email: software@bruker-axs.de

Phone: +49 (721) 50997 5200

Fax: +49 (721) 50097 5222

Please include the version number of the product and the operating system version including its service pack level in your bug report. Attach the files which caused problems, if possible.